

High Frequency Electrosurgical Unit  
Model YR02144  
**Service Manual**





OUR SERVICES

## Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



### Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



### Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



### #Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



### Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



### Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



### Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.





## COPYRIGHT

1. Without the written permission of our company, any content of the manual should by no means be duplicated or disclosed to any other individual or organization.
2. Individual or organization should in no way use the manual unless authorized by our company.
3. Our company has the sole right to change or modify this manual.  
Any other individual or organization thereby should in no case change or modify the manual.
4. Without our written permission, any organization or individual should not, by the name of sale or donation, provide the public the manual or the duplication of the manual.
5. Anyone who violates our copyright must be pursued according to the copyright law.



## Table of Contents

<b>General</b>	2
Performance Description	3
Operating Power Supply	3
Operating Environment	3
Transportation	3
Configuration and Weight	4
Inspection of Accessories	4
Application of Accessories	5
<b>Internal Structure</b>	7
Integral Structure	8
Panel Installation	9
Base Plate Installation	10
Diagram for Integral Machine	11
Service parts List	12
<b>Instruction to Panel</b>	13
Introduction to Control and Display Area	14
Introduction to Accessory Junction Area	19
Introduction to Back Panel	20
<b>Technical Parameters</b>	21
Maximum Output Parameter	22
Load Curve	22
Output Power and Open Circuit Voltage Crest Peak Value Curve	26
Safety Indexes	29
<b>Operating Principle</b>	30
Block Diagram for Principal Machine	31
Wiring Diagram for Integral Machine	32
Theoretical Description	33
Theoretical Description of Subsidiary Panel	34
<b>Debugging of Integral Machine</b>	45
Power Setting	46
Regulation of Short Circuit Current	49
Regulation of High Frequency Creepage Current	49
Regulation of Affixing Area	50
<b>Fault List</b>	50
Precise Fault List	51
Fault Code Table	54
<b>Appendix</b>	69
Diagram for Control Panel	60
Diagram for Display Panel	61
Diagram for Power Supply Panel	62
Diagram for Power Panel	63
Diagram for Output Panel 1	64
Diagram for Output Panel 2	65
Diagram for Back Panel	66

## General

This part includes the following seven sections:

- Performance Description
- Operating Power Supply
- Operational Environment
- Transportation
- Configuration and Weight
- Inspection of Accessories
- Application of Accessories

## Performance Description

YR02144 high frequency surgical applications are applicable to normal surgical, gynecological, facial, orthopedic, urological, cerebral, dermatological and thoracic surgeries. It is available for such work modes as simple cut, mixed cut, double pole cut, point coagulation, soft coagulation, surface coagulation and double pole coagulation.

The detail performance description is shown as follows:

1. It belongs to Class I CF equipment.
2. Operating Mode: intermittent loaded continuous operating
3. The double pole mode is available for both cut and coagulation, of which, the double pole coagulation can be divided into typical coagulation and precise coagulation.
4. The single pole mode is available for both cut and coagulation, of which the single pole cutting can be divided into simple cut, mixed cut 1, mixed cut 2 and mixed cut 3, and the single coagulation can be divided into soft coagulation, point coagulation and surface coagulation.
5. Testing system for affixing area of the pole plate
6. Single pole coagulation is available for concurrent output.
7. The volume is available for automatic regulation (except for alarm).
8. It can be equipped with argon device.

## Operating Power Supply

220-240Vac, 50Hz power supply properly grounded

## Operational Environment

Temperature: 5-40°C

Humidity: ≤80%RH

Range of Atmospheric Pressure: 700hPa~1060hPa

Clean Indoor environment which is free of any corrosion, flammable substance and has sound ventilation.


## Transportation

Temperature: -40 ~ +55°C

Humidity: 10~100%RH



Handle with care



Upwards



Waterproofing

Range of Atmospheric Pressure: 500hpa ~ 1060hpa

It is necessary to guard against any violent vibration, collision and rains during the progress of transportation of utilization.

## Configuration and Weight

Packing Size (mm): 550\*520\*400


Gross Weight (kg): 23

Overall Dimension (mm): 516\*410\*160

Net Weight (kg): 20

## Inspection of Accessories


### Manual Control Scalpel



- ① Cutting Button  
Cutting Terminal
- ② Coagulation Button  
b Coagulation Terminal
- ③ Scalpel Head  
c Common Terminal


- Manual control scalpel head ③ is switched into conduction with common terminal c.
- Start-up cutting button ① is switched into conduction with the cutting terminal of the manual control scalpel plug a as well as the common terminal c.
- Start-up coagulation button ② is switched into conduction with the coagulation terminal b of the manual control scalpel plug a as well as the common terminal c.

### Foot Control Scalpel




- ① Scalpel Head
- ② Plug
- Foot Scalpel Head ① is switched into conduction with plug ②.

### Nipper and Connecting Line



- Nipper ① is switched into conduction with plug a.
- Nipper ② is switched into conduction with plug b.

### Foot Switch



- ① Cutting Pedal  
a Cutting Terminal
- ② Coagulation Pedal  
b Coagulation Terminal
- ③ Plug  
c Common Terminal

- Start-up cutting pedal ① and cutting terminal a of plug ③ are switched into conduction with common terminal c.
- Start-up coagulation pedal ② and coagulation terminal b of plug ③ are switched into conduction with common terminal c.

## Application of Accessories

### Manual Control Scalpel

#### Directions for Use

1. Manual control scalpel is used in combination with principal highfrequency electric scalpel machine.
2. It should subject to the formaldehyde fumigation prior to operation.
3. It is a must to check whether the electrode is firmly fixed, any loose one should be replaced.
4. Insert the plug of the manual control scalpel head into the corresponding output plug jack.
5. Remove the protective cover of the scalpel head, and press the“CUT” or “COAG” to select either the electric cutting mode or electric coagulation mode.
6. It is proposed that the scalpel should be cleaned with medical tampon dipped with 70% alcohol for medical purpose during the operation.
7. Operational Environment  
Temperature: 5-40°C  
Relative Humidity: ≤80%

#### Note:

1. Please do not immerse the product in the detergent or sterilizing agent solution when performing the cleaning and disinfection.
2. As the product belongs to single pole electric scalpel head, one of itselectrodes should be connected with the body of the patient prior to operation.

### Foot Control Scalpel

#### Directions for Use

1. Foot control scalpel is used in combination with principal highfrequency electric scalpel machine
2. It should subject to the formaldehyde fumigation prior to operation.
3. It is a must to check whether the electrode is firmly fixed, any looseone should be replaced
4. Insert the plug of the foot control scalpel head into the correspondingoutput plug jack
5. Remove the protective cover of the scalpel head, and press the “CUT” or “COAG” to select either the electric cutting mode or electric coagulation mode
6. It is proposed that the scalpel should be cleaned with medical tampondipped with 70% alcohol for medical purpose during the operation
7. Operational Environment  
Temperature: 5-40°C Relative  
Humidity: ≤80%

#### Note:

1. Please do not immerse the product in the detergent or sterilizing agent solution when performing the cleaning and disinfection.
2. As the product belongs to single pole electric scalpel head, one of itselectrodes should be connected with the body of the patient prior to operation.



### Soft Pole Plate Connecting Line

#### Directions for Use

1. This product can be used to connect the principal high frequency electric scalpel machine with soft pole plate.
2. It should subject to formaldehyde fumigation prior to operation.
3. It is a must to check whether the plug or the clip is firmly fixed, and replace the loose connecting line of soft pole plate with new one.
4. Insert the plug of the connecting line of the soft pole plate into the corresponding output plug jack, and clamp the soft pole plate with the other end of the clip.
5. Operational Environment  
Temperature: 5-40°C  
Relative Humidity: ≤80%
6. It is proposed that the product should be cleaned with medical tampon dipped with 70% alcohol for medical purpose during the operation

#### Note:

Please do not immerse the product in the detergent or sterilizing agent solution when performing the cleaning and disinfection

### Nipper Connecting Line

#### Directions for Use

1. This product can be used to connect the principal high frequency electric scalpel machine and nipper.
2. It should subject to formaldehyde fumigation prior to operation
3. It is a must to check whether the plug or the plug jack is firmly fixed, and replace the loose connecting line of soft pole plate with new one
4. Insert the plug of nipper connecting line and the other end into the corresponding plug jack and nipper respectively.
5. Operational Environment  
Temperature: 5-40°C  
Relative Humidity: ≤80%
6. It is proposed that the product should be cleaned with medical tampon dipped with 70% alcohol for medical purpose during the operation

#### Note:


Please do not immerse the product in the detergent or sterilizing agent solution when performing the cleaning and disinfection

## Internal Structure



This part includes the following four sections:

- Integral Structure
- Panel Installation
- Base Plate Installation
- Diagram for Integral Machine

## Integral Structure



## Panel Installation







Diagram for Integral Machine




## Introduction to Panel

This part includes the following three sections:

- Introduction to Control and Display Area
- Introduction to Accessory Junction Area
- Introduction to Back Panel

### Introduction to Control and Display Area



**User's Application Status**

1	System Initiation Mode Programming Key	20	Equipment Model
2	System Initiation Mode Setting Key	21	Single Pole Plate Mode Indicator Lamp
3	System Initiation Mode Number DisplayWindow	22	Double Pole Plate Mode Indicator Lamp
4	Double Pole Cutting Power Set ValueDisplay Window		Alarm Indicator Lamp
5	Double Pole Cutting Start-up IndicatorLamp	24	Pole Plate Affixing Area Display Window(Light Arrangement Window)
6	Double Pole Power Setting Increment Key	25,26	Double Pole Precise Coagulation Mode Indicator Lamp and Selector Key
7	Double Pole Cutting Power SettingDecrement Key	27,28	Double Pole Typical Coagulation Mode Indicator Lamp and Selector Key
8	Double Pole coagulation Set Value Display Window	29	Foot Switch Double Pole Operation Indicator Lamp
9	Double Pole Coagulation Start-up IndicatorLamp	30	Foot Switch Working Mode Selector Key
10	Double Pole Coagulation Power Setting Increment Key	31	Foot Switch Single Pole Operation IndicatorLamp
11	Double Pole Coagulation Power Setting Decrement Key	32,33	Simple cutting Mode Selector Key and Indicator Lamp
12	Single Pole Cutting Power Set ValueDisplay Window	34,35	Mixed Cutting 1 Mode Selector Key and Indicator Lamp
13	Single Pole Cutting Start-up Indicator Lamp	36,37	Mixed Cutting 2 Mode Selector Key and Indicator Lamp
14	Single Pole Cutting Power Setting Increment Key	38,39	Mixed Cutting 3 Mode Selector Key andIndicator Lamp
15	Single Pole Cutting Power Setting Decrement Key	40,41	Soft Coagulation Mode Indicator Lampand Selector Key
16	Single Pole Coagulation Power Set Value Display Window	42,43	Point Coagulation Mode Indicator Lampand Selector Key
17	Single Pole Coagulation Start-up IndicatorLamp	44,46	Surface Coagulation Mode IndicatorLamp and Selector Key
18	Single Pole Coagulation Power Setting Increment Key	45	Trademark
19	Single Pole Coagulation Power Setting Decrement Key		



### Cross Reference List of Debugging and Application Statuses

No.	Application Status	Debugging Status
2	System Initiation Mode Setting Key	Feedback signal indicating switching between main and feedback CPU.
3	System Initiation Mode Number Display Window	Feedback Signal Display Window (U, I, P, u, r)
4 8	Double Pole Power Set Value Display Window Double Pole Coagulation Set Value Display Window	Feedback Value Display Window
6 7	Double Pole Cutting Power Setting Increment Key Double Pole Coagulation Power Setting Decrement Key	Feedback Signal Switching Key (switching between U, I, P, u and r)
26 28	Double Pole Precise Coagulation Mode Selector Key Double Pole Typical Coagulation Mode Selector Key	Setting Mode Switching Key
14 15	Single Pole Cutting Power Setting Increment Key Single Pole Cutting Power Setting Decrement Key	Set Point Selector Key (0-9)
12	Single Pole Cutting Power Set Value Display Window	Set Point Position Display Window
16	Single Pole Coagulation Power Set Value Display Window	PWM Display Window
32	Simple Cutting Mode Selector Key	Setting Mode Selector Key (current, power)
34	Mixed Cutting 1 Mode Selector Key	Switching between Feedback AD Value and Measured Value
41	Soft Coagulation Mode Selector Key	Power Setting Key
18 19	Single Pole coagulation Power Setting Increment Key Single Pole Coagulation Power Setting Decrement Key	PWM AD Value Increment Key PWM AD Value Decrement Key

### Description of Debugging Status

#### Debugging Status Access Mode:

Switch on the machine and press and hold the double pole precisecoagulation mode selection key and soft coagulation mode selection key for 5 seconds.



**Debugging Status Function:**

- a. It will produce a continuous PWM which is used to test the control panel, power supply panel and power panel.  
The single coagulation display window will display a continuous PWM value which can be changed by pressing the single coagulation power setting increment key or decrement key. At the time of initiation, the main CPU will send out a set PWM value.
- b. The feedback value indicating the voltage, current, power, mode and affixing area is used to check whether the machine is in normal operation.
- c. It will set the accurate power for the user's application status (for detail, please refer to Part 6).

**Specific Application of Feedback Value**

- a. At debugging status, the system initiation mode number display window indicates types of various feedback signals. By pressing the double pole power setting increment or decrement key to make switch between U, I, P, u, and r, the double pole cutting power set value display window and double pole coagulation power set value display window will display relevant feedback values.

U--- Output feedback voltage. I---

Output feedback current P---


Output feedback power









u--- Mode signal

r--- Affixing area of pole plate


- b. Press the system initiation mode setting key to make switch between main CPU and feedback CPU. At the time of initiation, the double pole cutting power set value key display window and double coagulation power set value display window will display the feedback value of main CPU. Once the system initiation mode setting key is pressed, the feedback signal displayed in the system initiation mode number display window will flash. At this point, the double pole power set value display window and double pole coagulation power set value display window will display the feedback value of feedback CPU.
- c. Value of feedback voltage, current or power varies in line with the variation of PWM and load. However, the variation of mode signal as well as the pole plate affixing area is so limited.
- d. It is applicable to make switch between feedback AD value and measured value by pressing and holding the mixed cutting 1 mode selector key for 2 seconds. At the time of initiation, the double pole cutting power set value display window and double coagulation power set value display window will display feedback values. Once the mixed cutting 1 mode selector key is pressed, a decimal point will be added to the feedback signal displayed in the system initiation mode number display window. At this point, measured values will be displayed in the double pole cutting power set value display window and double coagulation power set value display window.

## Introduction to Accessory Junction Area



- |   |   |   |   |   |  |
|---|---|---|---|---|--|
| 1 |    | Power Switch  |  | Single Pole Output Channel 2 (socket of manual control scalpel 2)                     |  |
| 2 |    | Pole Plate Socket   |  | Full Suspension Tuning Balance Output   |  |
| 3 |  | Nipper socket   | 7   |  | CF De-quiver Equipment   |
| 4 |  | Single Pole Output Channel 1 (manual control scalpel 1 and foot control scalpel socket) | 8   |  | It means that all output sockets possess hazardous high voltage. |

## Introduction to Back Panel



- 1 Nameplate
- 2 Foot Switch
- 3 , 4 4A Fuse Protector
- 5 Volume Tuning Knob
- 6 Power Socket

## Technical Parameters

This part includes the following four sections:

- Maximum Output Parameter
- Load Curve
- Output Power and Peak Value Curve of OpenCircuit Voltage
- Safety Indexes


### Maximum Output Parameter

WORKING MODE	MAXIMUM POWER (W)	RATED LOAD ( $\Omega$ )	MAXIMUM OPEN CIRCUIT VOLTAGE PEAK VALUE (V)	
Single Pole Cutting	Simple Cutting	400	2300	
	Mixed Cutting 1	300	2700	
	Mixed Cutting 2	200	3200	
	Mixed Cutting 3	150	3700	
	Single Pole Coagulation	Soft Coagulation Point	120	3400
		Coagulation Surface	100	5600
		Double Pole Coagulation	120	7250
Double Pole Cutting	Double Pole Cutting	70	1450	
	Precise Coagulation	Typical coagulation	70	240
			200	360

### Load Curve

The following chart shows the curve line of output power P at full power and half power settings in various modes, which changes with load  $R_L$  (the full power and half power is respectively indicated with solid line and dash line)

a. Simple Cutting (not applicable when the load is below 100 $\Omega$ )




400W

200W


2  
1




b. Mixed Cutting 1 (not applicable when the load is below 100Ω)




c. Mixed Cutting 2 (not applicable when the load is below 100Ω)




d. Mixed Cutting 3 (not applicable when the load is below 100Ω)




e. Soft Coagulation(not applicable when the load is below 100  $\Omega$ )




f. Point Coagulation(not applicable when the load is below 100 $\Omega$ )




g. Surface Coagulation (not applicable when the load is below 100 $\Omega$ )




h. Double Pole Cutting (not applicable when the load is below 100Ω)



g. Double Pole Precise Coagulation (not applicable when the load is below 100Ω)




i. Double Pole Typical Coagulation (not applicable when the load is below 100Ω)




## Output Power and Open Circuit Voltage Peak Value curve


a. Simple Cutting




b. Mixed Cutting 1




c. Mixed Cutting 2




d. Mixed Cutting 3




e. Soft Coagulation




f. Point Coagulation




g. Surface Coagulation




h. Double Pole Cutting



i. Double Pole Precise Coagulation



j. Double Pole Typical Coagulation



### Safety Indexes

- 1. Pressure Endurance** (It shall be deemed as being qualified if it subject to 5mA excessive creepage for 1 minute without any punch-through or flash)

No.	Item Tested
1	Grounded Network Power Supply: AC1500V
2	Grounded Shell Insulating Parts: AC4000V
3	Grounded Single Pole Surgery Electrode: AC4700V
4	Single Pole Surgery Electrode to Network Power Supply: AC5600V
5	Grounded Neutral Electrode: AC4700V
6	Neutral Electrode to Network Power Supply: AC5600V
7	Grounded Bipolar Electrode: AC3000V
8	Bipolar Electrode to Network Power Supply: AC4200V

### 2. Low Frequency Creepage Current (mA)

		Normal Status	Single Failure Status
Grounded Creepage Current		<0.5	<1
Shell Creepage Current		<0.1	<0.5
Patient Creepage Current		<0.01	<0.05
Patient Auxiliary Creepage Current	DC	<0.01	<0.05


## Operating Principle

This part includes the following three sections:

- Block Diagram for Principal Machine
- Wiring Diagram for Integral Machine
- Theoretical Description
- Introduction to Principle of Subsidiary Panel



## Block Diagram for Principal Machine



# Wiring Diagram



J6006


J6007 J6008

J6009

J6010

J  
6  
0  
1  
2

31



## Theoretical Description

Press the POWER (power supply control) button to switch on the power. At this point, the electric scalpel is at stand-by status, and the control circuit is under the control of operators (mode, power and so on). During the progress of stand-by and operation, the main CPU will continuously detect the pole plate signal. In case of broken circuit or other abnormal condition occurs to the pole plate, the main CPU will give out alarm to prohibit the start-up. If the pole plate signal is normal, the control panel will send a mode signal, a pulse width modulation signal and a sound indicating signal (not provided in the diagram) to the power panel. After receiving the mode signal from the control panel, the power panel will send a modulation signal of certain cycle and width through the mode selection and by producing a circuit. This signal will be transmitted to the impulse circuit for amplification prior to being converted to power amplifying grid signal through high frequency transformer. At the same time, the control panel will send out pulsewidth modulation power signal which is corresponding to the set power to the power panel. Through the amplification by the switching circuit and high frequency transformer, a stable and isolated DC current can be provided. Furthermore, the output terminal is also equipped with a current sampling circuit, through which, the sampled signal will be transmitted to the control switch circuit of comparator to limit the short circuit current.

During the power output of electric scalpel, the main CPU and the monitoring CPU will perform real-time inspection of the high frequency voltage, current and power signal from the output terminal so as to realize the safety redundancy. At the same time, with the software computation and compensation, the output power of the electric scalpel will become stable, which is helpful in ensuring the safe and reliable operation of the electric scalpel.

## Power Supply Panel

### Precise Operating Principle

The major function of the power supply panel is to make use of the principle of switching power supply to provide power amplification power. It receives a PWM signal from the control panel, and converts into DC signal which is to be sent to the U1004 (UC3825 High Speed PWM Controller) so as to produce two parallel alternating square waves. Driven by the U1002 and U1003 (7667), the two square waves shall be transmitted to the coil T1001 and T1002 which are used to drive 4 VMOS tubes (IRF840). Furthermore, the 220V AC voltage passing through the bridge rectifying circuit (BR1001) and filter circuit can obtain a DC voltage, which shall be switched into conduction with 4 VMOS tubes (840) by turns prior to the output from the transformer T1001. It will be converted into DC voltage through the bridge rectification by diode (U860) and output amplification by L and C filter.

### Interface Description

J1001: ~230V Input

J1002: Power amplifying DC voltage output, connected with power panel J2005

J1003: Connected with power panel J2002 J1005

, J1006: Short circuit plug


R1012: Load resistor installed on the (470Ω) base plate. M1001~M1004: VMOS tube (IRF840) on the large-size radiator.

For diagram for printing panel, please see appendix.

### Functional Detection

1. The J1001 (~230Vplug) should be pulled out prior to the detection of waveform.
2. Detect the two alternating square waves produced by the UC3825 high-speed PWM controller.


Make use of the dual trace oscilloscope to observe the testing points--- TP1018 and TP1019 as indicated in the following diagram.




Waveform: Alternating Square Wave, Overlapped high levels unacceptable.

3. Driving Waveform of VMOS Tube.

Observe the g-s waveform of 4 VMOS tubes as indicated in the following figure with oscilloscopes.





4. Connect the J1001 (~230V plug), pull out the short circuit of J1005 and make use of the multi-meter to measure pin 1 and 2 (1 refers to grounding pin, and the span of the multi-meter shall be higher than DC 200V) of J1005.
5. Switch on the machine to enter the power setting debugging status, initiate simple cutting, and the voltage value will increase with the increment of AD value of the PWM until it reaches 200V.

**Function of J1005 short circuit plug:** Pin 2 and 3 of J1005 in short circuit are used to send the power from the switch power source to the power amplifying circuit. When the short circuit plug is pulled out, it is applicable to check whether the switch power source is in normal operation.

## Control Panel

### Precise Operating Principle

The major functions of the control panel are performed by two CPU. UD1 is the main control chip, of which, the major function is to receive start-up signal and pole plate mode signal, send out mode signal and control the power signal of PWM, control various lights through manipulating the 7219 on the display panel and scan the keyboard. Furthermore, the main control CPU is also used to receive power, voltage and current feedback signals and control the final power output. UD2 is the feedback chip, of which, the major function is to receive the feedback signals from the relay as well as power, voltage and current feedback signals. It will give alarm in case of emergency. The control panel aims to perform the AC/DC conversion of voltage and current feedback from the output panel through AD536, and produce the power value through AD633. The control panel is equipped with pole plate detecting circuit. Any signal transmitted from the control panel to the power panel is isolated via the photo-coupler 521.

### Interface:

J1: To be connected with power panel J2001. J5: To be connected with display panel J20.  
For diagram for printing panel, please see appendix.

### Functional Detection

1. Each key of the display panel and display are normal.
2. Pole plate affixing area indication is normal (for detail, please see part 6).
3. Current, voltage, power and mode signal feedback is normal (for detail, please see part 6).
4. PWM signal is normal.
  - Pull out the J1005 short circuit plug from the power supply panel, and makes of the multi-meter to measure pin 1 and 2 of J1005 (1 refers to the grounding pin, and the span of the multi-meter should be more than the DC 200V).
  - Switch on the machine to enter the power setting debugging status, initiate simple cutting, and the voltage value will increase with the increment of AD value of the PWM until it reaches 200V.

### Precise Operating Principle

The major function of the display panel is to receive signals from main control CPU of the control panel, and control the display of the nixie light through two MAX7219. The display panel is also equipped with 4 reverse driver of 2803 model which is used to drive various lights and keys.

## Display Panel

### Interface

J20: Connected with J5 of control panel

For diagram for printing panel, please see appendix.

### Functional Detection

1. The display of each nixie light is normal.
2. The indicator lamp is normal.
3. Each key is in normal operation.
4. The start-up indicator lamp flashes when the single coagulation and double pole coagulation are initiated.

## Power Panel

### Precise Operating Principle

The major function of the power panel is to output relevant power tube driving signals in line with the mode provided by the control panel in addition to the power amplifying power and output power from the power supply panel. The power panel has two power output circuits, namely bridge power amplifying circuit and single tube power amplifying circuit. The single-tube power amplifying circuit is applicable to point and surface coagulation modes; whereas, the bridge power amplifying circuit is applicable to other modes. The power panel is also equipped with special circuit which is used to limit the short circuit current and the high frequency creepage current. The power panel will send a signal to the power supply panel for the purpose of reducing the voltage of the power amplifying circuit in case of excessive short circuit current and high frequency creepage current (except for surface coagulation mode). Similarly, when the high frequency creepage current exceeds the rated value in the surface coagulation mode, it will send a signal to the mode circuit to change the driving signal.

### Interface

J2001: Connected with control panel J1

J2002: Connected with power supply panel J1003 J2003:

Connected with 5V and  $\pm 12V$  power modules J2004: Connected with back panel J7003

J2005: Connected with power supply panel J1002

J2006: Connected with large-size coil and J6015 of output panel 2 J2007:

Connected with 12V power module


J2008: Connected with J5002 of output panel 1 J2009:

Connected with J5001 of output panel 1 M2001~M2008:


Connected with large-size radiator For diagram for printing panel, please see appendix.

### Functional Detection

1. Pull out the short circuit plug J1005 from the power supply panel.
2. Detect the driving signal of full-bridge power amplifying circuit.
  - a. Adjust the full bridge.
    - Observe the testing points TP2002 with oscilloscope and adjust the PR2001 so as to make sure that the pulse width is 100ns (as indicated in the following figure).



○ Observe the testing pints TP2006 and TP2008 (as indicated in the following figure) with dual trace oscilloscope.



Waveform: Alternating Square Wave, Overlapped high level is unacceptable.

b. Observe the testing pints TP2006 and TP2007 (as indicated in the following figure) with dual trace oscilloscope.





c. Observe the testing pints TP2008 and TP2009 (as indicated in the following figure) with dual trace oscilloscope.




d. Detect the driving waveform of the VMOS tube, initiate the simple cutting and observe the g-s waveform of power tube (P11NM60) with oscilloscope.




e. Perform the mode detection and observe the g-s waveform of power tube (P11NM60) with oscilloscope.


Observe the simple cutting, soft cutting, double pole cutting, doublepole typical coagulation and precise coagulation as indicated in the following diagram.




Waveform of mixed cutting 1 is as shown in the following figure:



Waveform of mixed cutting 2 is shown as follows:




Waveform of mixed cutting 2 is shown as follows:



f. Connect the voltage and current monitoring meter with power supply panel J1005.  
(the measuring span of voltmeter shall be more than DC 200V, and that of the ammeter shall be more than DC 5A)


g. Switch on the machine to enter the debugging status with scalpel head and pole plate subjected to no load. Once the simple cutting is initiated, the voltage value as displayed in the monitoring meter will increase with the increment of the AD value of PWM until it reaches 200V. The current value should not exceed 400mA.

3. Detection of Driving Signal from Single-tube Power Amplifying Circuit:
  - a. Pull out the short circuit plug of power supply panel J1005, and then switch on the machine to enter the debugging status for power setting.
  - b. Initiate the surface coagulation, observe the testing points TP2004 with oscilloscope and adjust the PR2008 so as to make sure that the pulse width is  $1.2\mu\text{s}$  (as indicated in the following figure), and the cycle is  $32\mu\text{s}$ .




Waveform: Cycle is  $32\mu\text{s}$ , and the pulse width is  $1.2\mu\text{s}$ .

- c. Initiate the point coagulation and observe the testing point TP2004 with oscilloscope (as indicated in the following figure).




Waveform: Cycle is  $16\mu\text{s}$ , and the pulse width is  $1.2\mu\text{s}$ .

- d. Detect the driving waveform of VMOS tube, initiate the surface coagulation, observe the g-s waveform of the power tube (450) with oscilloscope and adjust the PR2009 so as to make the falling edge of the waveform straight and guard against any small pulse on the negative terminal. See the following figure.



- e. Initiate the point coagulation and observe the g-s waveform of power tube as indicated in the following figure with oscilloscope.



- f. Connect the voltage and current monitoring meter with power supply panel J1005 (same as that of the full bridge).
- g. Switch on the machine to enter the debugging status with scalpel head and pole plate subjecting to no load. Once the simple cutting is initiated, the voltage value as displayed in the monitoring meter will increase with the increment of the AD value of PWM until it reaches 200V. The current value should not exceed 600mA.

## Output Panel 1

### Precise Operating Principle

The major function of output panel 1 is to produce an oscillation via NE555, which is used in turn to produce isolated initiation signal for manual control scalpel through the coil T5001, T5002 and self-restraint photo-coupler O5001-O5004.

### Interface

J5001: Connected with power panel J2009  
J5002: Connected with power panel J2008  
J5003: Connected with J6011 of output panel 2  
J5004: Connected with J6013 of output panel 2  
J5005: Connected with J6001 of output panel 2

### Functional Detection

1. Insert the manual control scalpel into the OUT1, and then switch on the machine to enter the operating status.
2. Initiate simple cutting and point coagulation respectively to activate the relevant start-up indicator lamps.
3. Insert the manual control scalpel into the OUT2.
4. Initiate simple cutting and point coagulation respectively to activate the relevant start-up indicator lamps.

## Output Panel 2

### Precise Operating Principle

The major function of output panel 2 is to make use of two independent resonant circuits to perform the full bridge and single tube power amplification so as to realize the uniform output via high voltage relay. The high voltage output part is equipped with current sampling circuit (T6006), high frequency creepage current sampling circuit (T6007), voltage sampling circuit (T6004), current sampling circuit (T6005) and pole plate resistance sampling circuit (T6003). All these will be transmitted to the power panel and the control panel for analysis.

### Interface

J6001: Connected with J5005 of output panel 1  
J6002, J6003: Connected with the 60<sup>th</sup> ring of large-size coil  
J6006: Connected with pole plate jack of the panel socket  
J6007, J6008: Connected with nipper jack of the panel socket  
J6009: Used to accept the foot control scalpel  
J6010: Connected with manual control scalpel 1 inserted into the panel socket  
J6011: Connected with J5003 of output panel 1  
J6012: Connected with manual control scalpel inserted into the panel socket  
J6013: Connected with J5004 of output panel 1  
J6014: Connected with pole plate switch of panel socket  
J6015: Connected with J2006 of power panel

### Functional Detection

1. Insert the manual control scalpel into OUT1, and then switch on the machine to enter the operating status.
2. Initiate the simple cutting and point coagulation respectively to activate power output.
3. Insert the manual control scalpel into the OUT2.



4. Initiate the simple cutting and point coagulation respectively to activate power output.
5. Insert the nipper and initiate the double pole cutting to activate the power output.

## Back Panel

### Precise Operating Principle

The major function of the back panel is to produce an oscillation via NE555 as well as isolated foot control scalpel start-up signal via coil T7002 and self-restraint photo-coupler O7002 and O7003. The back panel is also equipped with a LM386 vocalizing circuit.

### Interface

J7003: Connected with J2004 of power panel J7004:  
Connected with foot switch fixed to the shell

### Functional Detection

Switch on the machine to enter the operating status, and initiate the simple cutting and point coagulation concurrently with foot switch without imposing the power. At this point, the indicator lamp will flash accompanied by different start-up voice.



## Debugging of Integral Pole

This part includes the following four sections:

- Power Setting
- Regulation of Short Circuit Current
- Regulation of High Frequency Creepage Current
- Detection of Affixing Area






## Power Setting

1. Switch on the machine to enter the debugging status.
2. Press the double pole precise coagulation mode selector key or double pole typical coagulation mode selector key to select the appropriate mode. At this point, the indicator lamp for the selected mode shall flash. If the mode indicator lamp does not flash, it means that the machine is in double pole cutting mode. See the following figure.



3. At this point, the single pole cutting power set value display window will display by turns the set mode positioning point as well as the current or power value to be set.

See the following figure:




Press the simple cutting mode selector key to make switch between required current value and power value. The displayed value at the default start-up status is the current value.

See the following figure:



The single pole cutting power set value display window aims to display the positioning point of the set mode. See the following figure:



N-M

N: N refers to set mode (1-A).

1: Simple Cutting; 2: Mixed Cutting 1; 3: Mixed cutting 2; 4: Mixed Cutting 3; 5: Soft Coagulation;

6: Point coagulation; 7: Surface Coagulation                      8: Double Pole Cutting; 9: Precise Double Pole Coagulation;

A: Typical Double Pole Coagulation M:

Set Point (0-9)

4. Once the set mode is selected, the output will be switched into conduction with rated load. Except for the precise double pole coagulation, of which, the rated load is 200Ω, the rated load in other modes should be up to 500Ω.

5. The initial point in each mode is the zero point. At this point, the power value to be set as well as the AD value of PWM as displayed in single pole coagulation power set value display window is also zero.

See the following figure:



6. Initiate this mode until the displayed AD value of PWM flashes for 3 times, and then press the soft coagulation mode selector key (power setting key) to set the value at zero point.

7. Change the displayed value of set mode position point to 1 or the N-1. At this point, the set power value will also be changed to the set value of 1 automatically.

8. Press the single pole coagulation power setting increment/decrement key (increment/decrement key used to adjust the Ad value of PWM) to adjust the AD value of PWM so as to make sure that the output power is the same as that at the point 1.

See the following figure:



9. Initiate this mode until the displayed AD value of PWM flashes for 3 times, and then press the soft coagulation mode selector key (power setting key) to set the value at 1 point.



10. Set in turn all the other power points (0-9).
11. In case of error of set power point, it is applicable to press the single pole cutting power setting increment key or decrement key to reset it.




12. After all the 9 power points in this mode are set, select the next mode for setting.

## Regulation of Short Circuit Current

1. YR02144 aims to control the short circuit of hardware used for single pole cutting and coagulation. At the same time, the software can also be used to control the short circuit.
2. In normal operation, the software is frequently used to control the short circuit, and the hardware only plays a role of protection.
3. Switch on the machine to enter the operating status, and impose a load of 100Ω on the scalpel head and pole plate.
4. Initiate the maximum power for simple cutting, and adjust the potentiometer PR2002 on the power panel so as to limit the current. After that, slowly release the potentiometer until it is not unable to control the current.
5. Initiate the maximum power for mixed cutting 1, and adjust the potentiometer PR2006 on the power panel so as to limit the current. After that, slowly release the potentiometer until it is not unable to control the current.

## Regulation of High Frequency Creepage Current

1. Put the machine to be tested on the high frequency creepage current testing stand, and then lead the power supply to the 253V voltage regulator. Connection should be performed in line with the following figure:



1-Power Network 2-High Frequency Electric Scalpel 3- Table made of insulating materials 4-Surgery Electrode

5- Neutral Electrode 6-200Ω Measuring Electric Resistance 7-High Frequency Ammeter

2. Power should be set at the maximum value in each mode.
3. Adjust each potentiometer on the power panel in each mode so as to make sure that the high frequency creepage current is less than 150mA.  
Simple Cutting: PR2003 Mixed  
Cutting: PR2004 Surface  
Coagulation: PR2005 Point  
coagulation: PR2007

## Detection of Affixing Area

1. Test the pole plate socket on the machine panel and install the resistance box.
2. Switch on the machine to enter the operating status.
3. Set the pole plate resistance at 120  $\Omega$ , adjust the potentiometer PR1 on the regulating control panel so as to keep the machine from giving alarm. At this point, the 2 bars of affixing area indicator flash.
4. Reset the resistance at 121 $\Omega$  so as to make the machine give alarm and 1bar of the optical strip flash.
5. When the set resistance between pole plates is 9 $\Omega$ , the machine should give alarm. Otherwise, it is necessary to adjust the potentiometer PR2 on the regulating control panel to make the machine give alarm.
6. When the resistance between pole plates is 10 $\Omega$ , the machine shouldnot give alarm. Otherwise, it is also necessary to adjust the PR2 to keep the machine from giving alarm.
7. Adjust the potentiometer PR1 again to make the machine give alarmwhen the resistance is set at 121  $\Omega$  instead of 120  $\Omega$ .
8. Repeat the procedure 2, 3 and 4 until the machine does not give alarm when the resistance is set between 10-120  $\Omega$  other than other resistance values.
9. Perform verification in line with the following table:

Pole Plate Mode	Status	Requirement
Single -plate Mode	Pole plate not inserted	Bar lights but alarm lamp extinguish
	Insert pole plate	No alarm, all the bar lights extinguish
Double -plate Mode (to be connected with resistance box))	0-9 $\Omega$	Alarm, all the bar lights flash
	10 $\Omega$	No alarm, all the bar lights flash
	125 $\Omega$	No alarm, 2 bars of bar lights flash
	>125 $\Omega$	Alarm, 1 bar of bar lights flash or all the bar lights extinguish
	Set the resistance level of resistance box at 20 $\Omega$ , initiate the simple cutting and then set the resistance level of resistance box at 27 $\Omega$ ±1 $\Omega$ .	Alarm

## Fault List

This part includes the following two sections:

- Precise Fault List
- Fault Code Table



## Precise Fault List

Fault	Fault analysis	Cause
Display panel nixie light and indicator lamp do not flash	Power switch indicator lamp does not flash, main power source fault	Check whether the main power supply is available, whether the fuse is burnt or whether the power switch is normal.
	Power switch indicator lamp flashes. No 5V and $\pm 12V$ voltage output from power module	Power module is damaged.
	There is 5V and $\pm 12V$ voltage output from power module	panel or display panel line fault
Unable to initiate the machine	Check the accessories (manual control scalpel or foot switch)	Relevant accessories are damaged.
	There is start-up voice	Display panel line fault
	The manual control initiation but the foot control initiation is available	Back panel line fault
	The foot control initiation but the manual control initiation is available	Output panel 1 line fault
	Both manual and foot controlled initiation are not available.	CPU panel line fault
No voice	No square wave on the R7007 of the back panel	CPU panel line fault
	Waveform can be detected on the trumpet.	Back panel line fault
	Waveform can not be detected on the trumpet.	The trumpet is damaged
The selected pole plate mode is not in conformity with pole plate indicator lamp.	Check the pole plate socket fixed to the exterior shell.	The switching relay fixed to the socket of the pole plate is damaged.
	No switching signal	CPU pole plate line fault
When the single pole plate is used, the machine gives alarm.	Check accessories	Pole plate line is open circuit
	Check whether the affixing area transducer has been firmly installed.	Whether the installation of output panel T6003 is firm.
	Detect the internal fault of the machine	CPU Panel line fault
The affixing area indicator lamp flashes and the machine gives alarm even though the double pole plate has not been affixed to the patient's body	Check the wiring of soft pole plate.	The connecting line of soft pole plate is in short circuit or the holding seat has not been fixed yet.
	If it is still the case when the pole plate plug is pulled out.	CPU Panel line fault
When the affixing area is affixed to the patient's body, the	Check the wiring of soft pole plate.	The connecting line of soft pole plate is in open circuit or the holding seat is poor in contact.



affixing area indicator lamp extinguish and the machine give alarm.	Internal fault of the machine	CPU Panel line fault
When the affixing area is affixed to the patient's body, the affixing area indicator will display the obvious increment or decrement of affixing area.	Deviation of detected value	Readjust the affixing area parameter.
No power output	Check accessories (manual control scalpel, foot control scalpel, nipper or connecting line)	Relevant accessories are damaged
	Check the amplifying power	See Part 5-Power Supply Panel
	Check the driving waveform	See Part 5-Power Panel
	Check the resonant circuit of output panel.	See Part 5-Output Panel 2



## Fault Code Table

The fault codes are displayed respectively in double pole coagulation power setting display window, single pole cutting power setting display window and singlepole coagulation display window. Detail forms are shown as follows:

Double Pole coagulation  
Power Setting Display  
window

**Er**

Single Pole Power  
Setting Display window

xxx  
Fault Category

Single Pole Coagulation  
Power Setting Display  
Window

xxx  
Fault Number

Fault code			Fault Category	Fault Implication
Er	Cur	000	Cur Error marking of curve	AD value of power at the current point is not higher than the preceding one
Er	Cur	001		Extraordinarily high power curve offsetting value
Er	Cur	002		AD value of voltage at the current point is not higher than the preceding one
Er	Cur	003		AD value of current at the current point is not higher than the preceding one
Er	Cur	004		PWM value at the current point is not higher than the preceding one
Er	Cur	005		The PWM feedback to the CPU for reservation is different from that of the main CPU.
Er	Cur	006		The P feedback to the CPU for reservation is different from that of the main CPU.
Er	Cur	007		The V feedback to the CPU for reservation is different from that of the main CPU.
Er	Cur	008		The I feedback to the CPU for reservation is different from that of the main CPU.
Er	Cur	009		The VWS feedback to the CPU for reservation is different from that of the main CPU.
Er	Cur	00A		The AD value of power at the current point feedback to the CPU is not higher than that of the preceding one.
Er	Cur	00B		The AD value of voltage at the current point feedback to the CPU is not higher than that of the preceding one.
Er	Cur	00C		The AD value of current at the current point feedback to the CPU is not higher than that of the preceding one.
Er	Cur	00D		The AD value of independent voltage transducer curve at the current point is not higher than that of the preceding point.
Er	Cur	00E		The AD value of independent current transducer curve at the current point is not higher than that of the preceding point.
Er	Cur	00F		The AD value of independent power transducer curve at the current point is not higher than that of the preceding point.

Fault code			Fault Category	Fault Implication
Er	CAL	100	CAL Error in measurement and computation	Non-conformance curve and power mode saved in the RAM
Er	CAL	101		Extraordinarily high power curve offsetting value
Er	CAL	102		AD value is higher than 0x3FF PWM
Er	CAL	103		value is higher than 0x3FFError AD
Er	CAL	104		value of PWM
Er	CAL	105		Analog 12VPower Error
Er	USr	200	USr Error in keyboard operation	No single pole cutting mode mark No
Er	USr	201		double pole cutting mode mark
Er	USr	202		No single pole coagulation mode mark No
Er	USr	203		double pole coagulation mode mark
Er	USr	204		No foot control output/double pole mode indicator lamp markNo plus
Er	USr	205		and minus key code
Er	EEP	300	EEP EEPROM read-write error	No acknowledge signal
Er	EEP	301		Read-in is not available
Er	COF	400	COF Error of communication between main CPU and feedback CPU	The number of data transmitted exceeds the maximum numberspecified
Er	COF	401		Failure to receive error message within 500msData
Er	COF	402		verification error
Er	COF	403		The number of character strings received is incorrect.
Er	COF	404		Incorrect prefix and postfix of data string
Er	COF	405		Error in CRC verification of character string
Er	COF	406		“#esc5\n” is received, communication failure
Er	COF	407		“#esc5\n” is received, communication failure, error in CRC
Er	COF	408		verification by the receiver
Er	COF	409		“#esc5\n” is received, error in parity check by the receiver, communication failure
Er	COF	40A	Incorrect contents of replied order	
Er	COF	40A	Incorrect data contents replied	

Fault code			Fault Category	Fault Implication
Er	INT	500	INT Interruption error	Error in selection of AD channel number
Er	INT	501		Error in selection of channel number during the AD switching-over
Er	FPU	600	FPU 60× Incorrect marking of curve lineby feedback CPU	AD value of power at the current point is not higher than the preceding one
Er	FPU	601		Extraordinarily high power curve offsetting value
Er	FPU	602		AD value of voltage at the current point is not higher than the preceding one
Er	FPU	603		AD value of current at the current point is not higher than the preceding one
Er	FPU	604		PWM value at the current point is not higher than the preceding one
Er	FPU	610	FPU 61× Error in measurement and computation	Non-conformance curve and power mode saved in the RAM
Er	FPU	611		Extraordinarily high power curve offsetting value
Er	FPU	612		AD value is higher than 0x3FF
Er	FPU	613		PWM value is higher than 0x3FF
Er	FPU	614		Incorrect backup of important data
Er	FPU	620	FPU 62× Error in keyboard operation	No single pole cutting mode mark
Er	FPU	621		No double pole cutting mode mark
Er	FPU	622		No single pole coagulation mode mark
Er	FPU	623		No double pole coagulation mode mark
Er	FPU	624		No foot control output/double pole mode indicator lamp mark
Er	FPU	625		No plus and minus key code
Er	FPU	630	FPU 63× EEPROM read-write error	No acknowledge signal
Er	FPU	631		Read-in is not available
Er	FPU	640	FPU 64× Error communication with feedback CPU	The number of data transmitted exceeds the maximum numberspecified
Er	FPU	641		Failure to receive error message within 300ms
Er	FPU	642		Data verification error
Er	FPU	643		The number of character strings received is incorrect.

Fault code			Fault Category	Fault Implication
Er	FPU	644		Incorrect prefix and postfix of data string
Er	FPU	645		"#esc5\n" is received, communication failure
Er	FPU	646		Incorrect contents of replied order
Er	FPU	647		Incorrect data contents replied
Er	FPU	650	FPU 65× Interruption error	Error in selection of AD channel number
Er	FPU	651		Error in selection of channel number during the AD switching-over
Er	FPU	660	FPU 66× Control and detection error	Invalid current limitation
Er	FPU	661		Invalid power limitation
Er	FPU	662		Invalid termination of power current
Er	FPU	663		Invalid termination of power
Er	FPU	664		power mode error
Er	FPU	665		Single pole/double pole relay feedback input error
Er	FPU	666		Single pole/double pole relay feedback input error
Er	FPU	667		Power cutting/coagulation relay feedback input error
Er	FPU	668		Power cutting/coagulation relay feedback input error
Er	FPU	669		Power output channel 1 relay feedback input error
Er	FPU	66A		Power output channel 2 relay feedback input error
Er	FPU	66B		Power output relay feedback input error
Er	FPU	66C		Incorrect AD value of Vref/2
Er	FPU	66D		Analog 12Vpower error
Er	FPU	700	CCd Error in initiation of power	Irrelevant power is outputted.
Er	FPU	701		Tim-out in terminating power output of feedback CPU
Er	FPU	702		Feedback CPU fails to output the output relay control signal error in time.

Fault code			Fault Category	Fault Implication
Er	FPU	703		Feedback CPU fails to output the power mode error in time.
Er	FPU	704		Feedback CPU fails to output the PWM error in time.
Er	FPU	705		Feedback CPU fails to terminate the output power error in time.
Er	tAG	800	tAG Error in loading and saving the operating parameter array	The current operating parameters are not in conformity with the structural variable of operating parameter array.
Er	tAG	801		Invalid data showing the structural variable of operating parameter array
Er	UnE	900	UnE Error in integrating the main/feedback CPU software with hardware	Inconsistent software version No.WDT error
Er	UnE	901		In consistent AD switching value of 5V/2 channel of main/feedback CPU
Er	UnE	902		In consistent AD switching value of 12V/4 channel of main/feedback CPU
Er	UnE	903		In consistent AD switching value of voltage channel of main/feedback CPU
Er	UnE	904		In consistent AD switching value of current channel of main/feedback CPU
Er	UnE	905		In consistent AD switching value of power channel of main/feedback CPU
Er	UnE	906		In consistent AD switching value of mode channel of main/feedback CPU
Er	UnE	907		In consistent AD switching value of pole plate resistance channel of main/feedback CPU
Er	UnE	908		In consistent AD switching value of pole plate resistance channel of main/feedback CPU
Er	UnE	909		Inconsistent PWM value of main/feedback CPU


## Appendix

- Diagram for Control Panel
- Diagram for Display Panel
- Diagram for Power Supply Panel
- Diagram for Power Panel
- Diagram for Output Panel 1
- Diagram for Output Panel 2
- Diagram for Back Panel






# Display Panel







# Power Panel




# Output Panel 1



## Output Panel 2



# Back Panel





## YR02144 Service parts List

NO.	Part Number	Name	PCS
1	7PB-YR02144DS-002	Display board	1
2	7PB-YR02144CB-002	Control board	1
3	7PB-YR02144PW-001	Power board	1
4	7PB-YR02144OP1-001	Output board 1	1
5	7PB-YR02144OP2-001	Output board 2	1
6	7PB-YR02144DY-002	Power supply board	1
7	7PB-YR02144HM-002	Back Panel board	1
8	7DY-MPS3012-001	12V Power Module (MPS-30-12)	1
9	7DY-MPT45B-001	5V ±12V Power Module (MPT-45B)	1
10	1VM-840-001	VMOS 840	4
11	1VM-W12NK9-001	VMOS W12NK90Z	2
12	1VM-P11NK6-001	VMOS P11NK60	4
13	1DP-R311-001	Fuse seat R311,250VAC,10A	2
14	1FS-T6.3A250-001	FUSE 6.3A	4
15	1DP-10A250-001	Power jack	1
16	1AS-PAGM04-001	Petal switch seat	1
17	YR02144-36	Radiator A	1
18	YR02144-37	Radiator B	1
19	1RX-471100-001	Power Resistor	1
20	2JT-D3H3, 74-001	Insulating screw casing	10
21	6GY-015012-001	Insulating gaskets	8
22	6GY-026020-001	Big Insulating gaskets	2
23	7HT-00T1-001	Power coil	1
24	1WC-10A110-SW1	Power line	1
25	1SW-RF1004-G01	Main power switch	1
26	1LT-FN2060-1006	Filter	1