




Model Series YR02937 to YR02944
Intelligent Hot Plate Magnetic Stirrer
Instruction Manual





INSTALLATION

Upon receiving the Hotplate-Stirrer, check to ensure that no damage has occurred in shipment. It is important that any damage that occurred in transport is detected at the time of unpacking. If you do find such damage the carrier must be notified immediately. After unpacking, place Hotplate-Stirrer on a level bench or table, away from explosive vapors. Ensure that the surface on which the unit is placed will withstand typical heat produced by the unit and place the unit a minimum of 30 cm from vertical surfaces. Always place the unit on a sturdy work surface. The Hotplate-Stirrer can be plugged into a properly grounded outlet. The 110V unit plugs into a 110-volt, 50/60 Hz source. The 220V unit plugs into a 220-volt, 50/60 Hz source.

MAINTENANCE & SERVICING

Hotplate-Stirrer is built for long, trouble-free, dependable service. It needs no user maintenance beyond keeping the surfaces clean. The unit should be given the care normally required for any electrical appliance. Avoid wetting or unnecessary exposure to fumes. Spills should be removed promptly after the unit has cooled down. Do not use a cleaning agent or solvent on the front panel or top plate which is abrasive or harmful to plastics, nor one which is flammable. Always ensure the power is disconnected from the unit prior to any cleaning.

INTENDED USE

These Hotplate-Stirrers are intended for general laboratory use.

ENVIRONMENTAL CONDITIONS

Operating Conditions: Indoor use only.

Temperature: 5 to 40°C (41 to 104°F)

Humidity: 20% to 80% relative humidity, non-condensing
Altitude: 0 to 6,562 ft
(2000 M) above sea level

Non-Operating Storage:

Temperature: -20 to 65°C (-4 to 149°F)

Humidity: 20% to 80% relative humidity, non-condensing

SAFETY INSTRUCTIONS

Please read the entire instruction manual before operating Hotplate-Stirrer.

WARNING! DO NOT use the Hotplate-Stirrer in a hazardous atmosphere

or with hazardous materials for which the unit was not designed. Also, the user should be aware that the protection provided by the equipment may be impaired if used with accessories not provided or recommended by the manufacturer or used in a manner not specified by the manufacturer. Always operate unit on a level surface for best performance and maximum safety.

DO NOT lift unit by the top plate.

CAUTION! To avoid electrical shock, completely cut off power to the unit by disconnecting the power cord from the unit or unplug from the wall outlet. Disconnect unit from the power supply prior to maintenance and servicing. Spills should be removed promptly after the unit has cooled down.

DO NOT immerse the unit for cleaning. Alkalis spills, hydrofluoric acid or phosphoric acid spills may damage the unit and lead to thermal failure.

CAUTION! The top plate can reach 350°C, **DO NOT** touch the heated surface. Always use caution. Keep the unit away from explosive vapors and clear of papers, drapery, and other flammable materials. Keep the power cord away from the heater plate. **DO NOT** operate the unit at high temperatures without a vessel/sample on the top plate. **WARNING!** Units are **NOT** explosion proof. Use caution when heating volatile materials. **DO NOT** operate the unit if it shows signs of electrical or mechanical damage. Earth Ground - Protective Conductor Terminal Alternating Current



14X14cm HOTPLATE-STIRRER SPECIFICATIONS

Top plate dimensions (L x W): 140 x 140mm

Electrical (50/60 Hz):

	Hotplate	Stirrer	Hotplate-Stirrer
220V/110v:	450 watts	450watts	450 watts

Fuses: 5mm x 20mm, 5-amp quick acting **Temperature**

range: ambient +5°C to 400°C**Max**

Temp for top plate surface: 325°C **Temperature stability:**

+/-3%

Speed range: 100 to 1800rpm

Speed stability: +/-2%

Capacity: 50- 500mL, gross weight should not exceed 20lbs

Controls: see diagram

** Below 100°C +/-2°C. Environmental and sample conditions permitting.*

NOTE: *On all units, the Max. temperature setting on the display is 325°C for hot*

plate top

18X18cm HOTPLATE-STIRRER SPECIFICATIONS

Top plate dimensions (L x W): 180 x 180mm

Electrical (50/60 Hz):

	Hotplate	Stirrer	Hotplate-Stirrer
220V/110V:	550watts	550 watts	550 watts

Fuses: 5mm x 20mm, 5-amp quick acting **Temperature**

range: ambient +5°C to 400°C**Max**

Temp for top plate surface: 325°C **Temperature stability:**

+/-3%

Speed range: 100 to 1800rpm

Speed stability: +/-2%

Capacity: 50-2000mL, gross weight should not exceed 20lbs

Controls: see diagram

** Below 100°C +/-2°C. Environmental and sample conditions permitting.*

NOTE: *On all units, the Max. temperature setting on the display is 325°C for hot*

plate top



23X23cm HOTPLATE-STIRRER SPECIFICATIONS

Top plate dimensions (L x W): 230 x 230mm

Electrical (50/60 Hz):	Hotplate	Stirrer	Hotplate-Stirrer
220V/110V:	750 watts	750 watts	750watts

Fuses: 5mm x 20mm, 5-amp quick acting Temperature

range: ambient +5°C to 400°CMax

Temp for top plate surface: 325°C Temperature stability:

+/-3%

Speed range: 100 to 1800rpm

Speed stability: +/-2%

Capacity: 50-10000mL, gross weight should not exceed 20lbs

Controls: see diagram

* Below 100°C +/-2°C. Environmental and sample conditions permitting.

NOTE: On all units, the Max. temperature setting on the display is 325°C for hot plate top



28X28cm HOTPLATE-STIRRER SPECIFICATIONS

Top plate dimensions (L x W): 280 x 280mm

Electrical (50/60 Hz):

	Hotplate	Stirrer	Hotplate-Stirrer
220V/110V:	1000 watts	1000 watts	1000 watts

Fuses: 5mm x 20mm, 5-amp quick acting **Temperature**

range: ambient +5°C to 400°C**Max**

Temp for top plate surface: 325°C **Temperature stability:**

+/-3%

Speed range: 100 to 1800rpm

Speed stability: +/-2%

Capacity: 50-13000mL, gross weight should not exceed 20lbs

Controls: see diagram

* Below 100°C +/-2°C. Environmental and sample conditions permitting.

NOTE: On all units, the Max. temperature setting on the display is 325°C for hot

plate top

HEATING OPERATING INSTRUCTIONS

The Hotplates and Hotplate-Stirrers have a micro-processor-controlled heater that is designed to bring samples to temperature quickly and accurately.

1. Getting ready:

a. Fix the support rod into the screw hole (at the back of the instrument); Adjust the height of the cross clip and fasten the glassware tightly with clip.



b. Place a vessel with solution and the appropriate accessories in the center of the top plate. This is important because the vessel should be over the hottest part of the top plate.

c. Put the inner thermocouple on the machine.




(Inner Thermocouple)

- D. R/min windows display the stir speed value.
- E. Set button: to set temperature value and stir speed
- F. Self-tuning Button / Translocated button
- G. Reduce stir speed and temperature Button
- H. Add Stir speed and temperature button. Stir Power Button


• • How to set Temperature value


- a. Turn on the power button in the lateral (as the following photos).


2. Operating Instruction for Temperature Adjustment






- A. PV windows displays the measured temperature value of the plate
- B. SV windows displays the settled temperature value.
- C. Heating Power Button





- b. Press  button, PV windows display the ambient temperature, and SV display the settled temperature last time. (Notice: The hotplate stirrer can reach the settled temperature automatically, if you don't set the temperature you need.)




c. Press  button 3 seconds, The last figure flashes in the SV windows located display

d. Press  or  to for more than 2 seconds quickly to set the temperature you required. The SV windows displays the changed value, PV windows display the ambient temperature.

e. Press **Self-tuning button / translocated button** 2 seconds again when the **1** temperature finished if you need higher temperature.


Then **2** flashes, press  or  to set the temperature you required for


tens digit. The SV windows displays the changed value, PV windows display the

ambient temperature. (Note: Press  again if you don't need the tens digit temperature, then settled temperature is finished)

F. Press **Self-tuning button / translocated button** 2 seconds again when the **2** temperature finished if you need higher temperature.

Then **3** flashes, press  or  to set the temperature you required for hundreds digit. The SV windows displays the changed value, PV windows display the

ambient temperature. (Note: Press  again if you don't need the hundreds digit temperature, then settled temperature is finished)

1 G: Press  again when three digits finished their set. Then The SV windows displays the settled value, PV windows display the ambient temperature.

• **How to use the self-tuning function.**

1. Turn off the power button in the lateral when the PV windows temperature is increasing higher than settled temperature.
2. Wait for the machine temperature return to ambient temperature.
3. Press self-tuning 8 seconds until your settled temperature flashes.
4. Self-tuning function is finished when the SV window stop flashing, then temperature value of PV windows is the same as SV windows.



• **Measured the liquid temperature adjustment**

a. Pull out the inner thermocouple at the back of the product by hand. The PV windows displays HHHH, SV windows displays the settled temperature last time. (As following)



(Outer thermocouple)








(Inner Thermocouple).



- b. Connect the outer thermocouple 4 holes with the hotplate stirrer 4 pins
PV windows display ambient temperature, and SV display the settled temperature.
- c. Put another stainless-steel part of outer thermocouple length 20mm- 30mm into the liquid you measured

PV windows display the changed temperature and SV display the settled temperature. d. Burn out protection Function, when outer thermocouple disconnected, window display "hhhh", green light off, hotplate stirrer stops working. When this occurred check it before using.

e. Pull out the outer thermocouple and put on the inner thermocouple again when you don't need to operate the temperature.

• **Stirring Speed Operation (Two conditions)** 1. **When the Heating Function is working**

- a. Press the  button for stirring on. (Notice: The r/min windows display the changed stir speed until reaching the settled stir speed if you don't set again.)
- b. Press the set button  2 seconds, The last figure **1** flashes in the SV windows
- c. Press the set button  2 seconds again, The last figure **4** flashes in the r/ min windows
- d. Press  or  to for more than 2 seconds quickly to set the stir speed you required. The r/min windows display the settled stir speed. The SV windows displays the changed value, PV windows display the ambient temperature.
- e. Press **Self-tuning button / translocated button** 2 seconds again when the **4-stir** speed finished if you need higher speed.

Then **5** flashes, press  or  to set the temperature you required for




tens digit.


F. You can set **6,7** like **4,5**



e. Press  button 2 seconds again when the settled stir finished.

2. When the Heating Function isn't working.



A. Turn on the power button in the lateral (as the following photos).

b. Press the  button for stirring on. The r/min window display the stir speed
(Notice: The r/min windows display the changed stir speed until reaching the settled stir speed if you don't set again.)


c Press the set button  2 seconds again, The last figure **4** flashes in the r/ min windows

d. Press  or  to for more than 2 seconds quickly to set the stir speed you required. The r/min windows display the settled stir speed. The SV windows displays the changed value, PV windows display the ambient temperature.


E. Press **Self-tuning button / translocated button** 2 seconds again when the **4-stir** speed finished if you need higher speed.

Then **5** flashes, press  or  to set the temperature you required fortens digit.

F. You can set **6,7** like **4,5**.

G. Press  button 2 seconds again when the settled stir finished.Note:
Then the r/min windows display the changed stir speed. And the figure can be stored next time if you don't set the stir speed again.



Front Panel Parts





- A. Heating Power Button
- B. Stir Power Button
- C. Set temperature and stir speed button
- D. Self-tuning Button
- E. Reduce stir speed and temperature Button
- F. Add Stir speed and temperature button
- G. PV windows displays the surface temperature value of the plate
- H. SV windows displays the settled temperature, stir speed value. R/min windows display the stir speed value.

• Standard accessories

Photos	Item	Quantity
	Can be Removed Holder 1. Length 280mm 2. Diameter 10mm 3. Material: stainless steel.	1 PC
	Clamp with screw	1 PC

	Inner Thermocouple 4 -spin	1 PC
	Outer Thermocouple 4-spins Stainless Steel Length: 235mm Wire length: 1000mm	1 PC
	•Designed specifically for round bottom flasks •PTFE-coated Length: 20mm Diameter: 10 mm	1 PC

• Optional accessories

- c. One ps clip

WARNINGS

- Use an electrical outlet that operates with a fuse or circuit breaker and a ground fault interrupt circuit (GFCI)
- Disconnect the mantle from the power supply prior to maintenance and servicing to avoid personal injury:
- Don't use in the presence of flammable or combustible materials, fire or explosion may result.
- Refer servicing to qualified personal only.
- Flask or vessels should be handled with care. Appropriate safety clothing, glasses,



<https://kalstein.eu/>


gloves, and coats should be worn when operating all mantles.

CAUTION! To avoid electrical shock, completely cut off power to the unit by disconnecting the power cord from the unit or unplug from the wall outlet. Disconnect unit from the power supply prior to maintenance and servicing.

Spills should be removed promptly after the unit has cooled down.

CAUTION! The heating mantle can reach 400°C, **DO NOT** touch the heated surface. Always use caution. Keep the unit away from explosive vapors and clear of papers, drapery, and other flammable materials. Keep the power cord away from the heating mantle.

DO NOT operate the unit at high temperatures without a vessel/sample into the heating mantle.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tlf: +33 467158849 / +33 680760710 / +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S

