



Steam Sterilizer
Model Series YR03391
Instruction Manual

Thank you very much for purchasing our Kalstein's Steam Sterilizer Model Series YR03391

Please read the "Operating Instructions" and "Warranty" before operating this unit to assure proper operation. After reading these documents, be sure to store them securely together with the "Warranty" at a hand place for future reference.



Warning: Before operating the unit, be sure to read carefully and fully understand important warnings in the operating instructions.



OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.






Thank you for choosing our steam sterilizers.

Prior to operating this instrument, please read the operations manual carefully and follow all installation instructions.

IMPORTANT NOTICE:

If you can't open the door, please unlock the door according to the instruction "How to open the door in the case of power outage" in the manual.

Need Maintenance



If this picture appears on the screen when power on, please call your dealer or local service maintenance. Your steam sterilizer needs a regular maintenance.



Table of Contents

| | |
|------------------------------------------------------------|-----------|
| 1. General | 1 |
| 1.1 Scope of manual | 1 |
| 1.2 Intended use | 1 |
| 1.3 General Safety Instructions | 1 |
| 1.4 Standards and directives | 1 |
| 1.5 Symbols | 1 |
| 2. Description of the autoclave | 2 |
| 2.1 Autoclave views | 2 |
| 2.2 Control panel | 2 |
| 2.3 Technical Specifications | 3 |
| 2.4 Packing content | 3 |
| 3. Installation | 4 |
| 3.1 General conditions | 4 |
| 3.2 Power supply connection | 4 |
| 4. Setup | 4 |
| 4.1 Basic set | 5 |
| 4.2 Fill the distilled Water Tank | 5 |
| 4.3 Preparation of the sterilization materials | 5 |
| 5. Operation | 6 |
| 5.1 Select the program | 6 |
| 5.2 Running the sterilization program | 6 |
| 5.3 Start the sterilization program | 6 |
| 5.4 End of cycle | 6 |
| 5.5 Manual abort the program | 7 |
| 5.6 Record of the cycle | 7 |
| 5.7 Printer | 7 |
| 5.8 Report | 7 |
| 5.9 Label | 9 |
| 5.10 About device | 9 |
| 6. Advance settings | 9 |
| 6.1 Parameter | 9 |
| 6.2 Unit | 10 |
| 6.3 Preheat | 10 |
| 6.4 Expiry date | 10 |
| 6.5 Water quality | 10 |
| 6.6 Last error | 10 |
| 6.7 Factory reset | 10 |
| 7. Maintenance | 11 |
| 7.1 Clean the distilled water tank | 11 |
| 7.2 Clean Chamber, door seal, trays and tray Rack | 11 |
| 7.3 Door adjustment | 11 |
| 7.4 Replacement of the door seal | 12 |
| 7.5 The drain valve | 12 |
| 7.6 How to open the door in the case of power outage | 12 |
| 8. Error codes | 13 |
| 9. Transportation and storage | 14 |
| 10. Safety devices | 14 |
| Appendix 1 - Water properties / Characteristics | 15 |
| Appendix 2 - Diagrams of the sterilization programs | 16 |
| Appendix 3 - Wiring diagram | 21 |
| Appendix 4 - Hydraulic diagram | 22 |



1. General

1.1 Scope of Manual

This manual contains information concerning the installation, operation and maintenance of the steam sterilizers. To ensure proper performance of the autoclave, the instructions given in this manual should be thoroughly understood and followed.

Keep the manual near to the sterilization in a readily accessible location for future reference.

1.2 Intended Use

The device designed for total elimination and/or inactivation of microorganisms from medical devices and related products, placed in sterilization wraps/packaging, using pressurized steam (i.e., moist heat) as the sterilizing agent; it is used for products non-sensitive to high temperature, water, or steam.

1.3 General Safety Instructions

- Read and understand this manual before attempting to install or operate the sterilizer.
- Make sure that all the installation conditions are fully complied with.
- Ensure that the supply voltage agrees with the supply voltage specified on the type plate of the sterilizer.
- This appliance must be grounded. Connect only to a properly grounded outlet.
- Do not cover or block any openings on this appliance.
- Use this appliance only for its intended use as described in this manual.
- Do not exceed the maximum weight limit of the loads specified in this manual.
- Do not operate this appliance if it has a damaged cord or plug, if it is not working properly, or if it has been damaged or dropped.
- Never must put into the sterilizer in flammables or explosives products.
- The sterilizer may not be operated in areas in which gas or any other explosive volatile substance is present.
- Installation and repair work should only be performed by authorized service technician. Work by unqualified persons could be dangerous and void warranty.

1.4 Standards and directives

The steam sterilizers were designed and have been produced in conformity with the following directives and standards:

Directives:

2014/68/EC Pressure equipment.

93/42/EEC Medical devices (Class II b).

Standards:

EN13060 Relative to small steam sterilizers.

EN61010-1 Safety regulations for laboratory devices-Part 1:General regulations.

EN61010-2-040 Safety regulations specific to sterilizers used in the processing of medical material.

EN61326-1 Electromagnetic compatibility regulations for laboratory devices.

1.5 Symbols

For safe operation, please pay close attention to the alert symbols below which can be found in the sterilizer or throughout this manual.



This symbol represents an electrical caution - ground protection.



Hot Surface

This symbol represents a warning of a potential hot surface.





Important safety information.

This symbol represents a warning for extra caution.

2. Description of the sterilizer




2.1 Sterilizer views



| | |
|----|---------------------------|
| | Date of manufacture |
| SN | Manufacture number |
| | Manufacturer Name |
| | Manufacturer Address |
| | EC-Representative Name |
| | EC-Representative Address |

- | | | |
|-------------------------|--------------------------------------------|---------------------------|
| 1. Door | 9. Control Panel | 17 Bacteriological filter |
| 2. Seal | 10. Printer (Optional) | 18. Circuit breaker |
| 3. Door lock | 11. USB port | 19. Power socket |
| 4. Door handle | 12. Main power switch | 20. Used water tank vent |
| 5. Chamber | 13. Drain connector (Distilled water tank) | 21. Condenser vent |
| 6. Door spindle | 14. Drain connector (Used water tank) | 22. Rating plate |
| 7. Distilled water tank | 15. Safety valve | |
| 8. LCD | 16. Label printer port (Optional) | |

2.2 Control panel



| | |
|----|-----------------------------|
| 1 | Temperature of the cycle |
| 2 | Program |
| 3 | Holding time |
| 4 | Pressure |
| 5 | Fill distilled water tank |
| | The quality of water is bad |
| | The water pump is working |
| 6 | Drain used water tank |
| 7 | Printer is connected |
| 8 | Door opens |
| | Door is closed |
| | Door is locked |
| 9 | Time |
| 10 | Current temperature |
| 11 | Up button |
| 12 | Down button |
| 13 | Enter button |
| 14 | Menu button |






2.3 Technical specifications

| Item | YR03391 | YR03392 | YR03393 | YR03394 |
|------|---------|---------|---------|---------|
|------|---------|---------|---------|---------|



| | | | | |
|--------------------------------------|----------------------------------------|-------------|-------------|-------------|
| Chamber(mm) | φ170 x 320 | φ200 x 360 | φ247 x 350 | φ247 x 450 |
| Rated Voltage | 110-130 VAC; 60 Hz, 220-240 VAC; 50 Hz | | | |
| Circuit breaker | F16A, F20A/400V | | | |
| Nominal power | 1600VA | | | |
| Sterilization temperature | 121°C /134°C | | | |
| Capacity of the distilled water tank | Approx 2.5 L (Water at level Max) | | | |
| | Approx 0.5 L (Water at level Min) | | | |
| Operation temperature | 5°C-40°C | | | |
| Operation relative Humidity | Max. 80%, non condensing | | | |
| Overall dimensions WxHxD(mm) | 420x370x525 | 420x370x595 | 490x455x600 | 490x455x690 |
| Net weight | 35kg | 38kg | 45 kg | 50 kg |
| Max. Noise level | <70 dB | | | |
| Atmospheric pressure | 76 kPa - 106 kPa | | | |

2.4 Packing content

| No. | Accessory | Quantity |
|-----|------------------------------------------------------------------------------------------------------------|-------------|
| 1 | Steam sterilizer | 1 |
| 2 | Instrument tray  | 3 For 18/23 |
| | | 2 For 8/12 |
| 3 | Instrument tray rack  | 1 |
| 4 | Instrument tray handle  | 1 |
| 5 | Door adjustment tool  | 1 |
| 6 | Draining hose  | 2 |
| 7 | Instructions manual | 1 |
| 8 | Door seal | 1 |



3. Installation

3.1 General conditions

Position the device on a plane surface with minimum capacity 60 kgs.

The sterilizer should be placed on a level worktable.

Improper water level in the chamber could cause a sterilizer malfunction.

Leave at least 10cm between the device rear part and the wall. The clearance required to open the door is 40cm.

Position the autoclave at such a height as to make it possible for the operator to check the whole sterilization chamber and carry out the normal cleaning operations.

The room where the device is installed must be enough ventilated.

Do not install the device near washing basins, taps, etc. where it is likely to be splashed.

Do not lean on the door when it is opened.


Do not place trays , papers, fluid containers, etc. on the sterilizer.

3.2 Power supply connection

Check the label on back panel o sterilizer to verify voltage rating for the unit. Failure to connect the autoclave to an appropriate power supply could result in damage to the unit, and electrical shock to personnel.

Plug power cord into a properly polarized and grounded receptacle rated. A dedicated circuit only used for the sterilizer is recommended.

Never connect the device pin to reductions of any type.



4 Setup

Open the door and remove all of the inner contents for unpacking.

Connect the power cord to an outlet of the appropriate voltage.

Turn on the main power switch on the right side. After switching on, the machine turns on the LCD and shows the door position, water level, working program, date, time and etc.

Note: The control panel will be locked for the initial 10 seconds after powering up for system initialization.


Notice: Before using the sterilizer or at any time the low water level icon  blinks, fill the distilled water tank with distilled water.



4.1 Basic Set

The "Basic Set" Menu permits to set the following options:

*Date *Time *Language

Select the "Basic Set" from the main menu and tap 

Select the item by taping the  button. The unit you selected will be lighted.

Adjust the value by taping   on. Tap  to the next item.

Tap  button to save and exit.

Program
Basic Set
Report
Label

DATE:23-04-17
TIME:09:10:08
LANGUAGE:ENG
Counter:12

Note: The Counter (cycle No) can not be set by the operator.





Abbreviation of language options

| | | | | | | | |
|-----|---------|-----|------------|-----|------------|-----|----------|
| CHN | Chinese | ENG | English | DEU | German | ESP | Spanish |
| PL | Polish | FR | French | HUN | Hungary | ROM | Romanian |
| NL | Dutch | LTU | Lithuanian | LAT | Latvian | CZE | Czech |
| ITA | Italian | RUS | Russian | PT | Portuguese | HR | Croatian |

4.2 Fill the distilled water tank

Ensure that the drain valve is closed.

Tap the button and open the water tank cover.



Use only high quality distilled water. (see Appendix 1)

4.3 Preparation of sterilization materials

For the most effective sterilization and to preserve the sample, please follow below:

*Clean instruments immediately after use.

*Treat the instruments by ultrasound cleaner.

*Residual chemicals left over after cleaning and disinfecting process may damage and corrode parts of the autoclave, always rinse off the instruments using distilled water.

*Follow instrument manufacturer's guidelines and recommendations for handling and cleaning instruments prior to sterilization.

*Check the manufacturer's instructions as to proper procedure for sterilizing of each item.

*Arrange the samples of different materials on different trays or with at least 3cm of space between them.

*Clean and dry instruments thoroughly before placing them into tray.

*Always insert a sterilization paper or cloth between the tray and sample to avoid direct contact.

*Arrange the containers (glasses, cups, test-tubes, etc) on one side or inverted position, avoiding possible water stagnation.

*Don't stack the trays one above the other or put them in direct contact with the walls of the sterilization chamber.

*Always use the instrument tray handle.

*Wrap the samples one by one or, if more tools have to be set in the same bag, verify that these are made of the same material.


*Don't use metallic clips, pins or other, as this jeopardizes the maintenance of the autoclave.

*Don't overload the trays over the stated limit (see appendix 2).



5. Operation

5.1 Select the program

Tap  button to the main menu, select "Program". You will see the available sterilization programs. See Appendix 2.




Select the program by   buttons and confirm it by taping . If you don't want to select a program you may Tap  button to exit .


5.2 Running the sterilization program.


After selecting program, the materials to be sterilized can now be placed on the tray, placed inside the chamber by the tray handle.

After the instruments are loaded, you may close the door.



5.3 Start the sterilization program.


After the  button is taped, the stage and the status of the current cycle will appear on the display. The sterilizers will perform the program automatically. (see appendix 2).



Total time or count down until completion

5.4 End of cycle


After cycle is completed, the printer will be activated and print out a report of the cycle (if the optional printer has been connected) or save the report in the USB drive (optional).



Caution: Always use the tray handle to load or unload the tray into the autoclave. Failure to do so can result in burning.




Note: If the power shut off during the cycle is working, the screen will show a special picture when power on again.





5.5 Manual abort of the program

It is possible to interrupt a started cycle prematurely. If you need to interrupt a cycle and remove the items urgently, you may hold the button for 3 sec  during the drying time to skip the dry cycle.

If you interrupt a cycle before it reaches the "Drying" step, the items inside the autoclave must be considered not be sterilized.



If you need to interrupt a cycle after the holding time of the sterilization cycle and during the drying step, the items inside the autoclave can be considered sterilized.



Caution: Depending on the status of the Cycle, steam can escape from the sterilization chamber when you open the door.

5.6 Report of the cycle

USB Flash memory (Optional)

A USB drive can be used as a method of storing a report of the cycle. To do so, insert the USB drive into the slot located on the service door of the sterilizer.

The information will automatically output directly to the USB drive after the cycle has completed. The name of the file is determined by the serial number of the machine and the cycle number.

For example:

The serial number is A12345. The cycle number is 00012.

The file name in the USB stick is A12345_00012E00.txt.


The last three numbers represent error code.

For example, E00 means no error. E01 means error E01.

5.7 Printer (Optional)

If installed you can see the Icon in the screen stop flashing.

The printer (Optional) will print a report of the cycle that just ended. At the end of each cycle the printer will print out a report of the cycle.



5.8 Report

Internal Memory

In this menu you can find the reports stored in the internal memory of the sterilizer

5.8.1 Select "Report" from the main menu and Tap  you will see the list of records.


5.8.2 Select the record by taping .



5.8.3 Tap  button.

Then the record is printed (if the optional printer has been connected) the report or/and is saved in the USB drive (optional).


Note: It can save in the internal memory only the reports of the last 20 cycles.

Note: The storage system is based on the principle of "first In-first-out".

Tap  button to exit.



When viewing printed data records, refer to the diagram below:



```
=====
```

Program: Wrapped
Temperature: 134 C
Pressure: 204 kPa
Drying Time: 2.0 MIN
Holding Time: 04.0 MIN

```
-----
```

| | Time | Temp. | Pressure |
|------------------|----------|-------|----------|
| | HH:MM:SS | C | kPa |
| Start | 11:38:02 | 028.2 | 002.5 |
| T1: | 11:50:46 | 115.9 | 081.0 |
| T2: | 11:51:17 | 107.2 | 020.6 |
| T3: | 11:51:29 | 111.0 | 058.8 |
| T4: | 11:51:52 | 107.3 | 020.1 |
| T5: | 11:52:26 | 112.8 | 069.3 |
| T6: | 11:52:52 | 107.5 | 020.1 |
| TS: | 11:57:47 | 134.5 | 211.9 |
| Max Temperature: | 135.4 | | |
| Min Temperature: | 134.2 | | |
| Max Pressure: | 221.1 | | |
| Min Pressure: | 210.7 | | |
| T7: | 12:01:48 | 135.2 | 218.8 |
| T8: | 12:08:01 | 100.6 | 000.1 |
| T9: | 12:17:45 | 098.7 | 001.1 |
| End | 12:17:47 | 098.7 | 001.6 |

```
-----
```

Cycle No.: 00022
Ster. Value: Success
Water quality: OFF, 000
Date: 20-05-2019
SN:A00001
Operator:
3BN23D 11100010V2.9.1.3
0000

```
=====
```



5.9 Label printer (Optional)

5.9.1 Connect the label printer and switch on.

5.9.2 Select "Label", Tap

5.9.3 Select the records by taping button.





5.9.4 Tap the button to the interface for setting the quantity.

5.9.5 Change the quantity of the label that you want to print by taping button.

The range is 1-99.

5.9.6 Tap button to print.

5.9.7 Tap button to go back.




5.10 About device

5.10.1 Select "About device", Tap .

5.10.2 It will show the version of the firmware and the serial number.

5.10.3 Tap button to go back.



6. Advance setting

The advance setting interface permits to set the following

Options:

***Parameter:** Permits to change the time of Holding time and Drying time.

***Unit:** Permits to change the unit of measure temperature and pressure.


***Preheat:** This option allows you to maintain the required temperature in the sterilization chamber and the steam generator to start a new cycle immediately for the next 60 minutes so to run a new cycle should expect preheating time from 3 to 5 minutes upon start.

If the option is disabled (OFF) once the sterilization cycle end immediately the Autoclave is no longer heated, so when you start a new cycle should expect a preheating time between 3 to 5 minutes upon start.

Note: To maintain the temperature for longer time it is recommended that after each cycle hold the door closed.

Enter the setting

Select the "Setup", tap the button. Input the password 1111. Tap the button to change the number, tap the button to change the position. Tap the button will enter setting interface after inputting the 4 passwords.



6.1 Parameter

Tap button to select Parameter. Tap button to enter the menu.

Select the program that you need to change by taping button. Tap button to enter the setting.



Adjust the parameter by taping the button.

Drying time: 0-60 minutes

121°C holding time: 1-60 minutes

134°C holding Time: 1-20 minutes

After you finish adjusting the parameter Tap button to save and return to the above menu.





6.2 Unit

Select the "Unit" to adjust the unit of temperature and pressure.

Tap button to enter the menu.

Select the item by tapping ton.

Change the unit by tapping on.

Pressure: kPa/bar/psi

Temperature: °C/°F

Tap button to save and exit.

Pressure: kPa
Temperature: C

6.3 Preheat

Select the "Preheat" to adjust the Preheat.

Tap button to enter the menu.

If you don't want to preheat after switch on, set the value to Off by tapping button.

Tap button to save and exit.

Preheat: On

6.4 Expiry date

Select the "Expiry date" to enter the setting interface.

Tap button to enter the menu.

Adjust the value by tapping the button. The range is 1-12

Tap button to save and exit.

Month: 03

6.5 Water quality

Select the "Water quality" to change the function.

Tap button to enter the menu.

If you don't want to have a alarm if the quality of water is bad, set the value to Off by tapping button.

Tap button to save and exit.

Water quality
On

6.6 Last error

Select the "Last error" to adjust the Preheat.

Tap button to enter the interface.

It will record the parameters of the sensors when the alarm appears.

Tap button to save and exit.

LAST ERROR: E30
2019-06-06 13:40
PC:13 ST:03 CN:00011
Pressure: 101kPa
T1:153.9 T3: 093.2°C
T2:028.1 T4: 220.5°C

6.7 Factory reset

Select the "Preheat" to adjust the Preheat.

Tap button to enter the interface.

Change Yes/No by tapping ton.

Tap button to confirm and exit..

The value of holding time and drying time will restore the default value if you confirm "Yes".

Factory reset
Yes



7.Maintenance

To assure proper operation and maximum autoclave life, carefully follow all recommendations for periodic maintenance. One of the MOST important steps you can take to prevent problems with your sterilizer is to use ONLY distilled water.

| Frequency | Number of cycles | Maintenance operation |
|----------------|------------------|-----------------------------------------------------------------|
| Monthly | 50 | Clean the door seal |
| | | Clean the filter inside the chamber and in the clean water tank |
| | | Clean the chamber the trays and the rack |
| | | Clean the external surface |
| Every 3 months | 200 | Clean the distilled water tank |
| Every year | 800 | Replace the door seal |


7.1 Clean the distilled water tank

Disconnect the main cable.

Drain the tank completely using the drain tube and leave it connected into the connector in a open position.

Clean the internal surface with a soft sponge and a small soft brush for the areas that are difficult to reach using and a mild soap.

Remove the filter and clean it with a small soft brush and mild soap, rinse it with distilled water and put it in to the position.



back

7.2 Clean Chamber, door seal, trays and tray Rack.


Remove the trays and tray rack from the chamber.

Clean trays, rack and inside of chamber with mild soap.

Rinse the trays, rack and inside of chamber with a smooth cloth and distilled water.

Examine door seal for possible damage.

Clean door seal and mating surfaces with a damp cloth.



Note: Do not use bleaching agents or any abrasive materials / substances in chamber. Failure to comply may result in damage to the chamber and/or other components.



Caution: To prevent burns, let unit to cool before cleaning gaskets and touch the surface.


7.3 Door adjustment

Under normal circumstances the chamber door does not require adjustments. However, if the seal fails (resulting in steam leaking from the front of the chamber), you may use the spanner tool to tighten the seal.


Open the door.

Insert the spanner tool in the gap beneath the plastic cover; use the spanner to grip the adjusting nut. the nut counter clockwise as the figure below. This will tighten the sealing plate.

Turn the nut until the sealing plate is tight. If the door knob is too tight, you may also turn the nut clockwise to loosen it.



door



Turn



Caution: Never adjust the chamber door while the door is closed.




7.4 Replacement of the door seal

Open the chamber door. Remove the door seal ring carefully by hand.

Clean the door seal ring carefully with a smooth cloth with distilled water.


Moisten the new seal with distilled water.

Insert the new seal and tap in sequence as follows:




Caution: Please ensure the chamber and the door are cooled prior to replacing the seal ring.

7.5 The drain valve





7.6 How to open the door in the case of power outage





8. Error codes

| Code | Description | Proposed solution |
|------|----------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------|
| E1 | Steam generator temperature sensor error | Power off & run a new cycle Contact your supplier if error persists. |
| E2 | Inner temperature sensor error | Power off & run a new cycle Contact your supplier if error persists. |
| E3 | Temperature sensor of the chamber wall error | Carefully ensure that the chamber wall is heated and contact your supplier |
| E5 | Fail to release the pressure | Power off & run a new cycle Contact your supplier if error persists. |
| E6 | Door lock problem during the cycle | Make sure you had closed the door properly. check the door switch |
| E7 | The pressure is too lower during holding time. | Contact your supplier if error persists. |
| E8 | The pressure is too high during holding time. | Contact your supplier if error persists. |
| E9 | Failure to hold temperature | Ensure the distilled tank isn't empty. Check the inner temperature sensor. Check somewhere for leaking. |
| E10 | The door locking system doesn't work. | The electromagnet of locking system doesn't work. The switch of locking system doesn't work. |
| E11 | Failure to preheat the steam generator | Power off & run a new cycle Contact your supplier if error persists. |
| E12 | Failure to preheat the chamber | Power off & run a new cycle Contact your supplier if error persists. |
| E16 | The pressure doesn't reach 0 in 5 minutes after drying period. | Contact your supplier if error persist |
| E18 | The filling water pump working time is overime | Check the water pump or Contact your supplier if error persist |
| N20 | Program manually interrupted | Holding the  button for 3 seconds |
| E24 | It takes too long time to enter the next status. | Check somewhere leaking. Or contact your supplier if error persists. |
| E28 | The pressure is overload. | Power off and contact your supplier if error persists. |
| N29 | Power failure during working. | A notification message. |
| E34 | The pressure is higher than 50kPa during drying. | The solenoid valves are blocked. |

Caution: You may cancel the voice of alarm by taping any button. And cancel the alarm by holding the  button for 3 seconds after you repair it. Then swith off and switch on again.



9. Transportation and storage

9.1 Switch off the sterilizer before transportation or storage.

9.2 Pull out the plug. Let the machine cool down.

9.3 Drain the distilled water tank and the used water tank.

Condition for transport and storage

Temperature: -20°C ~ +50°C

Relative humidity: $\leq 85\%$

Atmospheric pressure: 50kPa~ 106kPa.

10. Safety devices

1.Main fuses: Protection the instrument against possible failures of the heating resistor.

Action: Interruption of the electric power supply.

2.Thermal cutouts on the main transformer winding: protection against possible short circuit and main transformer primary winding overheating

Action: Temporary interruption of winding.

3.Safety valve: Protection against possible sterilization chamber over-pressure.

Action: Release of the steam and restoration of the safety pressure.

4.Safety micro-switch for the door status: Comparison for the correct closing position of the door.

Action: Signal of the wrong position of the door

5.Thermostat on chamber heating resistors: Protection for possible over heating of the chamber heating resistors.

Action: Interruption of the power supply of the chamber resistors.

6.Thermostat on steam generator heating resistors: Protection for possible overheating of the steam generator heating resistors.

Action: Interruption of the power supply of the steam generator resistors.

7.Door safety lock: Protection against accidental opening of the door.

Action: Impediment of the accidental opening if the door during the program.

8.Self-leveling hydraulic system: Hydraulic system for the natural pressure leveling in case of manual cycle interruption, alarm or black-out.

Action: Automatic restoration of the atmospheric pressure inside chamber.



Appendix 1

Water properties / Characteristics

| Description | Feed water | Condensate |
|-------------------------------|-------------------------|-------------------------|
| Evaporate residue | $\leq 10\text{mg/l}$ | $\leq 1.0\text{mg/kg}$ |
| Silicium oxide SiO_2 | $\leq 1\text{mg/l}$ | $\leq 1.0\text{mg/kg}$ |
| Iron | $\leq 0.2\text{mg/l}$ | $\leq 0.1\text{mg/kg}$ |
| Cadmium | $\leq 0.005\text{mg/l}$ | $\leq 0.05\text{mg/kg}$ |
| Lead | $\leq 0.05\text{mg/l}$ | $\leq 0.1\text{mg/kg}$ |
| Rest of heavy metals | $\leq 0.1\text{mg/l}$ | $\leq 0.1\text{mg/kg}$ |
| Chloride | $\leq 2\text{mg/l}$ | $\leq 0.1\text{mg/kg}$ |
| Phosphates | $\leq 0.5\text{ mg/l}$ | $\leq 0.1\text{mg/kg}$ |
| Conductivity | $\leq 15\mu\text{s/cm}$ | $\leq 3\mu\text{s/cm}$ |
| PH Value | 5 – 7.5 | 5-7 |
| Appearance | Colorless, clean | Colorless, clean |
| Hardness | 0.02 mmol/l | 0.02 mmol/l |



Appendix 2

Diagrams of the sterilization programs(YR03394)

| Programs | Te | a) kP | Ho | Total | Type | Ma | Max |
|-------------------|-----|----------|----|-------|-----------------------------------------|------|------|
| | | Pre | | | | | |
| SOLID | 134 | 210 | 4 | 17-30 | Unwrapped solid material | 5.00 | 1.50 |
| | 121 | 110 | 20 | 30-40 | | | |
| WRAPPED | 134 | 210 | 4 | 20-30 | Unwrapped solid material | 5.00 | 1.50 |
| | 121 | 110 | 20 | 30-35 | Single-wrapped solid or hollow material | 4.00 | 1.20 |
| TEXTILE | 134 | 210 | 8 | 24-34 | Unwrapped porous material | 1.25 | 0.40 |
| | | | | | Single-wrapped porous material | 1.10 | 0.30 |
| | | | | | Dual-wrapped porous material | 0.75 | 0.25 |
| | 121 | 110 | 30 | 40-45 | Single-wrapped hollow material | 4.00 | 1.25 |
| | | | | | Dual-wrapped solid and hollow material | 2.00 | 0.60 |
| PRION | 134 | 210 | 18 | 34-44 | Unwrapped porous material | 1.25 | 0.40 |
| | | | | | Single-wrapped porous material | 1.10 | 0.30 |
| | | | | | Dual-wrapped porous material | 0.75 | 0.25 |
| | | | | | Single-wrapped hollow material | 4.00 | 1.25 |
| | | | | | Dual-wrapped solid and hollow material | 2.00 | 0.60 |
| LIQUID (optional) | 121 | 110 | 30 | 45-55 | Liquid | 1.20 | 0.40 |
| Drying (optional) | — | — | — | 1-20 | — | — | — |



Diagrams of the sterilization programs(YR03393)

| Programs | Te | Pre | Ho | Total | Type | Ma | Max |
|-------------------|-----|-----|----|-------|-----------------------------------------|------|------|
| | | | | | | | |
| SOLID | 134 | 210 | 4 | 17-27 | Unwrapped solid material | 4.50 | 1.20 |
| | 121 | 110 | 20 | 27-37 | | | |
| WRAPPED | 134 | 210 | 4 | 22-32 | Unwrapped solid material | 4.50 | 1.20 |
| | 121 | 110 | 20 | 32-42 | Single-wrapped solid or hollow material | 3.50 | 1.10 |
| TEXTILE | 134 | 210 | 8 | 26-36 | Unwrapped porous material | 1.00 | 0.30 |
| | | | | | Single-wrapped porous material | 0.80 | 0.25 |
| | | | | | Dual-wrapped porous material | 0.60 | 0.20 |
| | 121 | 110 | 30 | 42-52 | Single-wrapped hollow material | 3.50 | 1.00 |
| | | | | | Dual-wrapped solid and hollow material | 1.50 | 0.50 |
| PRION | 134 | 210 | 18 | 36-46 | Unwrapped porous material | 1.00 | 0.30 |
| | | | | | Single-wrapped porous material | 0.80 | 0.25 |
| | | | | | Dual-wrapped porous material | 0.60 | 0.20 |
| | | | | | Single-wrapped hollow material | 3.50 | 1.00 |
| | | | | | Dual-wrapped solid and hollow material | 1.50 | 0.50 |
| LIQUID (optional) | 121 | 110 | 30 | 43-53 | Liquid | 1.00 | 0.30 |
| Drying (optional) | — | — | — | 1-20 | — | — | — |



Diagrams of the sterilization programs(YR03392)

| Programs | Temp | a) kPa | Holds | Total | Type | Max | Max |
|-------------------|------|-----------|-------|-------|-----------------------------------------|------|------|
| | | Pre | | | | | |
| SOLID | 134 | 210 | 4 | 15-25 | Unwrapped solid material | 3.00 | 1.00 |
| | 121 | 110 | 20 | 25-36 | | | |
| WRAPPED | 134 | 210 | 4 | 15-25 | Unwrapped solid material | 3.00 | 1.00 |
| | 121 | 110 | 20 | 27-30 | Single-wrapped solid or hollow material | 2.50 | 0.80 |
| TEXTILE | 134 | 210 | 8 | 20-30 | Unwrapped porous material | 0.75 | 0.25 |
| | | | | | Single-wrapped porous material | 0.50 | 0.15 |
| | 121 | 110 | 30 | 37-40 | Dual-wrapped porous material | 0.30 | 0.10 |
| | | | | | Single-wrapped hollow material | 2.50 | 0.80 |
| | | | | | Dual-wrapped solid and hollow material | 1.20 | 0.40 |
| PRION | 134 | 210 | 18 | 30-40 | Unwrapped porous material | 0.75 | 0.25 |
| | | | | | Single-wrapped porous material | 0.50 | 0.15 |
| | | | | | Dual-wrapped porous material | 0.30 | 0.10 |
| | | | | | Single-wrapped hollow material | 2.50 | 0.80 |
| | | | | | Dual-wrapped solid and hollow material | 1.20 | 0.40 |
| LIQUID (optional) | 121 | 110 | 30 | 40-50 | Liquid | 0.80 | 0.25 |
| Drying (optional) | — | — | — | 1-20 | — | — | — |







Diagrams of the sterilization programs(YR03391)

| Programs | Te | a) kP Pre | Ho | Total | Type | Ma | Max |
|----------------------|-----|-----------------|----|-------|-----------------------------------------|------|------|
| | | | | | | | |
| SOLID | 134 | 210 | 4 | 15-25 | Unwrapped solid material | 2.00 | 0.60 |
| | 121 | 110 | 20 | 25-35 | | | |
| WRAPPED | 134 | 210 | 4 | 15-25 | Unwrapped solid material | 2.00 | 0.60 |
| | 121 | 110 | 20 | 25-30 | Single-wrapped solid or hollow material | 1.50 | 0.50 |
| TEXTILE | 134 | 210 | 8 | 19-29 | Unwrapped porous material | 0.50 | 0.15 |
| | | | | | Single-wrapped porous material | 0.35 | 0.10 |
| | | | | | Dual-wrapped porous material | 0.25 | 0.10 |
| | 121 | 110 | 30 | 35-40 | Single-wrapped hollow material | 1.50 | 0.30 |
| | | | | | Dual-wrapped solid and hollow material | 1.00 | 0.30 |
| PRION | 134 | 210 | 18 | 29-39 | Unwrapped porous material | 0.50 | 0.15 |
| | | | | | Single-wrapped porous material | 0.35 | 0.10 |
| | | | | | Dual-wrapped porous material | 0.25 | 0.10 |
| | | | | | Single-wrapped hollow material | 1.50 | 0.30 |
| | | | | | Dual-wrapped solid and hollow material | 1.00 | 0.30 |
| LIQUID (Optional) | 121 | 110 | 30 | 25-35 | Liquid | 0.60 | 0.20 |
| Drying (Optional) | — | — | — | 1-20 | — | — | — |

The time required for sterilizer to be ready for routine use after the power is switched is less than 5 minutes.

The max. Temperature of the 134°C sterilization cycle is 137°C


The max. Temperature of the 121°C sterilization cycle is 124°C





Appendix 3

Wiring diagram




- TP1: Steam generator temperature sensor
- TP2: Inner temperature sensor 1
- TP3: Temperature sensor of chamber wall
- V1: Air release valve
- V4: Water release valve



Appendix 4

Hydraulic diagram



V1: Air release valve

V4: Water release valve