









**Proforma invoice**
**Invoice No: SFSMI20230817**
**Date: April 17th, 2023**
**Customer:** Optimum Business Center 450 Rue Baden Powell,  
 • 34000 Montpellier, France.  
**Email:** Tlf: +33 467158849 / +33 680760710 / +33 663810023  
**Tel:** https://kalstein.eu  
**Add:** KALSTEIN FRANCE, S. A.

**From: Patty xiao**  
**Email: sfslab\_3@foxmail.com**  
**Tel: +(86) 020-2988 1968 Fax: +(86) 020-2988 1968**  
**Add: No.2,Shijiao commercial ST,Huadong town,Huadu district,GZ ,China**

Item number	Description	Image	Specifications	Dimension (mm)	Quantity (pc)	Unit price (USD)	Total Price (USD)
YR06013	Side workbench		Cabinet in full stainless steel 304, both the top and the cabinet. 	2000*750*800	1	\$885	\$885
YR06015	Conner workbench		Cabinet in full stainless steel 304, both the top and the cabinet.	1000*1000*800	1	\$500	\$500
YR06017	SS Faucet and SS sink		SS triple way faucet and SS sink -304# stainless steel , 3mm thick .		1	\$142	\$142
YR06024	Balance bench (scale)		Frame: 1.0mm cold roll steel . Epoxy resin powder sprayed for steel frame. Worktop: 12.7 mm thickness Chemical resistant phenolic resin, put the balance position with 18 mm thick marble countertop, shock-reduction design, to ensure the accuracy use of laboratory equipment. Adjustable caster: Use shock absorption casters of our own design. Two sockets with anti-splash cover (two & three plugs) Built-in three-stage shock absorption spring under main table: Highest accuracy millionth	900*600*800	1	\$300	\$300
YR06023	Drip rack		Drip rack- high density PP Type: single side	550*400*120	1	\$25	\$25
YR05868	F9 -BT-W1800 fume hood		1. Worktop : no worktop 2. Shell: 1.0mm thick steel covered with a grayish epoxy-polyester paint. 3. The inner lining board: 5mm thick compact panel (anti-corrosion, anti-acid) 4. Sliding door: 5mm thick tempered glass 5. Lighting: LED lamp. 6. Power: 110V with American Plug , Socket: 110V Two 10A multifunctional anti- splash sockets. 7. Fan: 110w axial flow fan with air velocity regulator. Surface wind speed: 0.3~0.6m/s 8. The exhaust pipe: $\varnothing$ 250mm 3meter soft exhaust air pipe. 10. All knock down packing, with hard carton.	1800*800*1500	1	\$532	\$532

YR05822	F9S-BT- W1800 SS fume hood		<p>1. Worktop : no worktop  2. Material: 304# stainless steel 1.0mm thick .  4. Sliding door: 5mm thick tempered glass  5. Lighting: LED lamp.  6. Power: 110V with American Plug , Socket: 110V Two  10A multifunctional anti- splash sockets.  7. Fan: Centrifugal  Blower with air velocity regulator.  Surface wind speed: 0.3~0.6m/s  8. The exhaust pipe: <math>\varnothing</math> 250mm 3meter soft exhaust air  pipe.  10. SS water tap one inner wall side .</p>	1800*800*1500	1	\$1,482	\$1,482
---------	-------------------------------------	--	---	---------------	---	---------	---------

<b>Total amount</b>							<b>\$3,866</b>
---------------------	--	--	--	--	--	--	----------------

**Notice:**  
1.Price term: EX-FACTORY .  
2.Payment terms: 100% in advance by T/T.  
3.Production time: 25 working days.  
4. Payment currency:USD  
Beneficiary account number: 101700 07047693  
Swift code: CHASHKHH  
Beneficiary country/region: Hong Kong, China  
Beneficiary name: Guangzhou Scienfocus Lab Equipment Co., Ltd.  
Beneficiary address: 26/F TOWER ONE TIME SQUARE 1 MATHESON STREET CAUSEWAY BAY HK  
Beneficiary bank: JPMorgan Chase Bank N.A., Hong Kong Branch  
Beneficiary bank address: Chater House, 8 Connaught Road Central, Hong Kong  
Bank code: 007  
Branch code: 863  
Remark: 157753842001026580 . Please include the order number

GUANGZHOU SCIENFOCUS LAB EQUIPMENT CO. LTD