




REAGENTS & CONSUMABLES SECTOR REAGENT



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Reagents
Clinical Chemistry Reagents
Model YRA1 // YRA14



Product Description:

Laboratory chemical reagents are chemical substances used in laboratories to make determinations of the content of other substances, the distribution of the constituents of the sample, to perform analysis, synthesis and finally to be able to do research and scientific work.

The blood lipid elements: TCH, TG, HDL-C LDL-C , APOA1, the APOB, Lp(a), HCY, TIR.

The elements of liver function: ALP, γ -GT, AST, ALT, TP, ALB, TBIL, DBIL, CHE, TBA, safe, ADA, 5'-NT, MAO, AFP, AMM, legs, ALC, GLDH.

Myocardium Zymogram elements: LDH, α -HBDH, CK-MB, CK, MB, ACE, LDHI, TNI, IMA.

The elements of renal function: UREA, UA, CREA, MALB, Cys-C, BMG, NAG, UTP, α 1-MG, RBP.

Glycometabolism elements: GLU, SGP, HbA1c, D3H, LAC.

Pancreas elements: α -Amy, LPS.

Pepsinogen elements: PGI, PGI.

Blood gas electrolyte elements: CL, Ca, P, Mg, CO₂, Na, K, Cu, Fe, Zn.

Rheumatism elements: RF, ASO, CRP.

The special protein elements: PALB, TRF, G6PD, DD3, C, C4, IgA, IgG, IgM, UIBC, CP.

Inspection of Hemagglutination topic: FB.



#Let'sgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 2 Rue Jean Lantier • 75001 Paris •
 +33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® Reagents
Clinical Chemistry Reagents
Model YRA1 // YRA14



Technical Description:

Model	Specification (713)/Box	Package/Box	Method	Reagent test
YRA1	R1:48mlx2 R2:12mlx2	120ml	IFCC Recommendation Method	ALT
YRA2	R1:48mlx2 R2:12mlx2	120ml	IFCC Recommendation Method	AST
YRA3	R1:48mlx2 R2:16mlx2	128ml	SZASZ Method	GGT
YRA4	R:20mlx8	160ml	Bromocresol Green Method	ALB
YRA5	R1:48mlx2 R2:12mlx2	120ml	Diazo 2, 4-2 chloro aniline colorimetric	TBIL
YRA6	R1:48mlx2 R2:12mlx2	120ml	Vanadate oxidation method	TBIL
YRA7	R1:48mlx2 R2:16mlx2	128ml	Enzymatic Cycling Assay	TBA
YRA8	R1:48mlx2 R2:16mlx2	128ml	Urease – glutamate dehydrogenase enzyme hydrolysis	UREA
YRA9	R1:48mlx2 R2:16mlx2	128ml	Enzymatic Assay Method	CREA
YRA10	R1:48mlx2 R2:16mlx2	128ml	Enzymatic Assay Method	K
YRA11	R1:48mlx2 R2:16mlx2	128ml	Enzymatic Assay Method	Na
YRA12	R:20mlx8	160ml	Mercuric Thiocyanate Method	Cl
YRA13	R1:48mlx2 R2:16mlx2	120ml	IFCC Recommendation Method	CK-MB
YRA14	R1:20mlx2 R2:20mlx2	80ml	Immuno Turbidimetric	CRP



#Let'sgivemore ❤️
 Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 2 Rue Jean Lantier • 75001 Paris •
 +33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:

 **Fundación
Jacinto Convit**
<https://www.jacintoconvit.org/>

 **ONE
TREE
PLANTED**

<https://onetreepanted.org>

 **Humatem**
Solidarité et coopération
biomédicale
<https://www.humatem.org/>


**MANIAPURE
FOUNDATION USA**
<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed

The 3 articles of the week!

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW

How does the climate affect your microscope?

Kalstein Blog

Article of the week:
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW

What is the importance of the laboratory refrigerator?

Kalstein Blog

Article of the week:
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER

2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89

Cel: +33 (0) 6 80 76 07 10

sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

