




REAGENTS & CONSUMABLES SECTOR REAGENT



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Reagents
Hematology Reagent
Model YRA45 // YRA71




Product Description:

Laboratory chemical reagents are chemical substances used in laboratories to determine the content of other substances, the distribution of sample components, to perform analysis, synthesis and finally to carry out research and scientific work.

This product is the cleaning solution of the reaction system in the detection process of the blood cell analyzer, so as to facilitate the in vitro detection of blood cells, hemoglobin, etc.

It is used to destroy red blood cells, dissolve hemoglobin, and maintain the shape of the cells to be analyzed before blood cell analysis, thereby facilitating cell classification and counting or quantitative determination of hemoglobin.



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 2 Rue Jean Lantier • 75001 Paris •
 +33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® Reagents

Hematology Reagent

Model YRA45 // YRA71



Technical Description:

YRA45	Diluent	20L
YRA46	Lyse	500ml
YRA47	Rinse	5L, 20L
YRA48	EZ Cleaner	100ml
YRA49	Probe Cleaner	50ml
YRA50	Diluent	20L
YRA51	EO-I Lyse	1L
YRA52	EO-II Lyse	500ml
YRA53	LH Lyse	500ml
YRA54	Cleaner	1L
YRA55	Probe Cleaner	50ml
YRA56	Diluent	20L
YRA57	EO-I Lyse	1L
YRA58	EO-II Lyse	500ml, 200ml
YRA59	LH Lyse	500ml
YRA60	Probe Cleaner	50ml
YRA61	Diluent	20L
YRA62	52 DIFF Lyse	500ml
YRA63	52 LH Lyse	100ml
YRA64	Probe Cleaner	50ml
YRA65	Diluent	20L
YRA66	EO-I Lyser	1L
YRA67	EO-II Lyser	500ml
YRA68	LH Lyser	500ml
YRA69	LBA Lyser	1L
YRA70	Cleaner	1L
YRA71	Probe Cleaner	50ml



#Let'sgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale


<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed

The 3 articles of the week!

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW

How does the climate affect your microscope?

Kalstein Blog

Article of the week:
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW

What is the importance of the laboratory refrigerator?

Kalstein Blog

Article of the week:
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER

2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89

Cel: +33 (0) 6 80 76 07 10

sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

