




Kalstein®

A different accompaniment, at your service

LABORATORY LINES

BIOSAFETY CABINETS AND HOODS





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Biosafety Cabinets and Hoods
Horizontal Laminar Flow Cabinet
Model YR05724



Product Description:

Desktop Horizontal Laminar Flow Cabinet is a local purification unit which is used in a clean aseptic and dust-free environment in the fields of medicine and health, biological engineering, scientific experiments, etc.

Product Features :

- ✓ The whole equipment adopts cold-rolled plate welding and the surface is electrostatically sprayed.
- ✓ The air supply mode of the equipment is divided into vertical air supply and horizontal air supply, almost closed glass damper, which is easy to operate.
- ✓ Small and exquisite, it can be freely placed on the general workbench for operation, which is convenient for use in small studios.
- ✓ The ventilation system is controlled by remote switch with double fast and slow speed to ensure the wind speed in the working area is always in an ideal state.

Technical Specifications :

Model	YR05724
Cleanliness level	Class 100@≥0.5μM (US Federal 209E)
Number of colonies	≤0.5pcs/dish·hour (Φ90mm culture plate)
Average Wind Speed	0.25~0.45m/s (Fast and slow)
Noise	≤62dB (A)
Vibration half peak	≤5μM (x、y、z direction)
Illuminance	≥300Lx
Power Supply	AC single-phase 220V/50Hz
Working area	690×350×490
W1×D1×H1	
Overall dimension	850×530×850
W×D×H	
High efficiency filter specification and quantity	670×455×38×1
Specification and quantity of fluorescent lamp/ultraviolet lamp	12W×1/20W×1



#Let'sgivemore ❤️
 Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 Optimum Business Center 450 Rue Baden Powell,
 • 34000 Montpellier, France.
 Tel: +33 467158849 / +33 680760710/ +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S.





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale

<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed
The 3 articles of the week!
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW
How does the climate affect your microscope?
microscope is an instrument that allows observing and measuring small objects. This is a scientific instrument that is used in laboratories, schools and research centers.
Kalstein Blog
Article of the week:
How does the climate affect your microscope?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW
What is the importance of the laboratory refrigerator?
Refrigerator is a device that is one of the most essential equipment that is used in laboratories, in a controlled environment (refrigerated space), without frost and humidity.
Kalstein Blog
Article of the week:
What is the importance of the laboratory refrigerator?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER
2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89
Cel: +33 (0) 6 80 76 07 10
sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

