




## LABORATORY SECTOR

HOTPLATE AND  
HOTPLATE STIRRER



**Kalstein**®

A different accompaniment, at your service





# Kalstein®

A different accompaniment, at your service

OUR SERVICES

## Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



### Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



### Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



### #Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



### Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



### Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



### Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



## KALSTEIN® Hotplate and Hotplate Stirrer

### Hotplate Magnetic Stirrer

Model YR02955



#### Product Description:

- ✓ 6-position Hotplate Magnetic stirrer
- ✓ Digital Display, Temperature: Rt -380°C, Mixing speed :100-1800 RPM
- ✓ Ceramic coated aluminum plate for corrosion resistance
- ✓ Specially designed housing prevents spill from damaging the electronic components inside
- ✓ Adjustable rotation speed and temperature
- ✓ OEM service

<b>Model</b>	YR02955
<b>Top Plate material</b>	White Ceramics
<b>Stir Range</b>	100~1800rpm
<b>Stir Speed Accuracy</b>	±1 rpm
<b>Top Ceramic Dimension for each position (MM)</b>	180 × 180
<b>Working Time</b>	Continue
<b>Position</b>	6
<b>Stir Speed Display</b>	Digital Display
<b>Max Temp for top Plate</b>	325°C
<b>Temperature Range for each position</b>	ambient to 400degree
<b>Temperature Accuracy</b>	±1 degree
<b>Temperature Control</b>	PID
<b>Voltage</b>	160-240V, 110-120V
<b>Temperature Display</b>	Digital display
<b>Outer Thermo</b>	to measure and control the liquid temperature value together with the front panel
<b>Inner Thermo .</b>	control and Measure the hotplate stirrer temperature
<b>Motor Power for each position</b>	40W, DC20V
<b>Product power for each position</b>	550W
<b>Stir Capacity for each position</b>	50ml-5000ml



#Letsgivemore ❤️

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
 2 Rue Jean Lantier • 75001 Paris •  
 +33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
**KALSTEIN FRANCE, S. A. S**







 **Kalstein**<sup>®</sup>  
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE  
ACQUISITION  
OF A KALSTEIN  
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación  
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



**Humatem**

Solidarité et coopération  
biomédicale

<https://www.humatem.org/>



**MANIAPURE  
FOUNDATION USA**

<https://www.maniapurefoundation.org/>



**Articles Kalstein**  
(Social Networks and Official Websites)



**Kalstein Blog**  
Not just be informed, be fully informed  
**The 3 articles of the week!**  
Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**1**  
Read this ARTICLE NOW  
**How does the climate affect your microscope?**  
Microscopes are an essential tool for scientific research, allowing researchers to observe and analyze small objects. This is why it is important to take care of your microscope and protect it from environmental factors such as temperature and humidity.

**Kalstein Blog**  
**Article of the week:**  
How does the climate affect your microscope?  
Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**2**  
Read this ARTICLE NOW  
**What is the importance of the laboratory refrigerator?**  
Refrigerators in a laboratory are one of the most important pieces of equipment that you need to maintain in a controlled environment (refrigerated space), without fluctuations in temperature.

**Kalstein Blog**  
**Article of the week:**  
What is the importance of the laboratory refrigerator?  
Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>

**Kalstein Training Videos**  
(Youtube Channel and Official Websites)



**Kalstein France - Jacinto Convit Foundation / Alliances**  
23 views • 3 years ago



**Kalstein's Microscopes**  
445 views • 1 week ago



PARIS - FRANCE  
**HEADQUARTER**  
2 Rue Jean Lantier  
75001 Paris - France

**Tlf:** +33 (0) 01 78 95 87 89  
**Cel:** +33 (0) 6 80 76 07 10  
[sales@kalstein.eu](mailto:sales@kalstein.eu)

<https://kalstein.eu/>  
<https://www.kalstein.fr/>

