




LABORATORY LINES

VORTEX MIXER



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Vortex Mixer

Multi-Plate Vortex Mixer

Model YR06044 // YR06047



Product Description:

- ✓ Digital display available; Stepless speed regulation; Easy touch/run switch.
- ✓ Ultra low stable operation of brushless body/motor with Power.
- ✓ With various accessories, suitable for centrifuge tube, blood tube, PCR plate, cell culture plate and so on.
- ✓ It can mix whatever you want, such as nail polishes, eyelash adhesives, paints, glues, tattoo ink, pigments and other resistant liquids.
- ✓ It applies to nail salons, beauty salons, artist studios, tattoo stores, schools, enterprises and home use.

Technical Specifications :

| Model | YR06044 | YR06045 | YR06046 | YR06047 |
|------------------------|--------------------------------------|---------|----------------------|---------|
| Operating mode | Touch / Run | | | |
| Control mode | Konb speed adjustment | | PWM speed adjustment | |
| Display | NO | | LCD | |
| Motor | Brushless motor | | | |
| Amplitude | 4.5MM circular oscillation | | | |
| Capacity | 50ml*4 / 15ml*9 / 5-10ml*12 / 2ml*24 | | | |
| Speed range | Touch: 0-3000 rpm | | Touch: 0-3000 rpm | |
| | Run: 200-3000 rpm | | Run: 100-3000 rpm | |
| Timing / Time | 1min – 8h | | | |
| Shell protection level | IP21 | | | |
| Voltage | AC 100-240V / DC24V | | | |
| Power | 20W | | | |
| Temperature | 5 – 40°C, 80%RH | | | |



#Let'sgive more ♥

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél: +33 467 158849 / +33 690760710 / +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A.S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale

<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed
The 3 articles of the week!
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW
How does the climate affect your microscope?
Microscopes are instruments used for observing microscopic objects. This is a sensitive device that requires a controlled environment to maintain its performance.
Kalstein Blog
Article of the week:
How does the climate affect your microscope?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW
What is the importance of the laboratory refrigerator?
Refrigeration is a key element of the most advanced equipment that is used in laboratories. It is a controlled environment (refrigerated space), without fluctuations in temperature.
Kalstein Blog
Article of the week:
What is the importance of the laboratory refrigerator?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Alliances
23 views • 3 years ago



Kalstein's Microscopes
445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER
2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89
Cel: +33 (0) 6 80 76 07 10
sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

