




LABORATORY SECTOR

SHAKERS



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Shakers

Orbital Shaker

Model YR00S240C // YR0RS-320C



Features:

- ✓ Large-diameter barbital oscillation, efficient mixing.
- ✓ Metal base, vibration stability. Aluminum alloy platform with oxidation coating, good corrosion resistance.
- ✓ Platform with quick release design, effectively improves the efficiency replacement.
- ✓ A variety of clamps can be selected to meet the experimental needs of various vessels.
- ✓ DC brushless motor, quiet operation, maintenance-free, long service life.
- ✓ Speed regulation function, timing function, LED display, easy to set and repeat experimental conditions, precise control of the experimental process.



#Letsgivemore 

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® Shakers

Orbital Shaker

Model YR00S240C // YR0RS-320C



Specifications:

Model	YR00S-240C	YR00S-350C	YR00S-350C Plus	YR0RS-320C
Shaking mode	Orbital			Rocking
Amplitude	Φ30mm	Φ30mm	Φ6mm	Max 8°
Speed range	30-250rpm	30-250rpm	30-650rpm	1-80rpm
Max loading	2Kg	5Kg		2kg
Max capacity	250ml*5	1000ml*5/500ml*9/250ml*16		/
Control	Microcomputer PID control			
Display	LED			
Work mode	continuous/timing			
Timer range	0-99h99min			
Standard clamp	platform with rubber pad			
Platform size	240*240mm	350*350mm		320*320mm
Overall size(mm)	240*316*109mm	485*330*140mm		320*446*109mm
Motor power	50W			40W
Ambient temperature	5-45°C			
Power supply	100-240V, 50Hz			

Optional Accessories:

		
standard platform	rubber bands fixed platform	spring rack
		
platform with bars	double layers platform	triangular flask clamp



#Letsgivemore 

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® Shaker

Open Air Shaker

Model YR00S-420C // YR00S-520C



Features:

- ✓ Orbital oscillation.
- ✓ Adjustable speed, timing function, LED display.
- ✓ Metal base, vibration stability. Aluminum alloy platform with oxidation coating, good corrosion resistance.
- ✓ DC brushless motor, quiet operation, maintenance-free, long service life.
- ✓ A variety of clamps can be selected to meet the experimental needs of various vessels.



#Letsgivingore ❤️

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® Shaker

Open Air Shaker

Model YR00S-420C // YR00S-520C



Specifications:

Model	YR00S-420C	YR0LS-420C	YR00S-520C
Shaking mode	Orbital	Linear	Orbital
Amplitude	Φ20mm	20mm	Φ20mm
Speed range	30-350rpm	30-350rpm	30-450rpm
Max loading	15Kg		25Kg
Max capacity	2000ml*5/1000ml*9/500ml*16/250ml*32		2000ml*9/1000ml*16/500ml*25/250ml*50
Control	Microcomputer PID control		
Display	LED digital display		
Work mode	continuous/timing		
Timer range	0-99h99min		
Platform size	420*420mm		520*520mm
Standard clamp	platform with rubber pad		
Overall size(mm)	450*476*169mm		590*650*190mm
Motor power	120W	120W	220W
Ambient temperature	5-45°C		
Power supply	100-240V, 50Hz		

Optional Accessories:

			
standard platform	rubber bands fixed platform	spring rack	platform with single layer bars / double layer bar
			
triangular flask clamp	double/triple layer platform	sticky platform	centrifugal tube platform



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale

<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed

The 3 articles of the week!

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW

How does the climate affect your microscope?

Kalstein Blog

Article of the week:
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW

What is the importance of the laboratory refrigerator?

Kalstein Blog

Article of the week:
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER
2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89
Cel: +33 (0) 6 80 76 07 10
sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

