



LABORATORY LINES

CENTRIFUGES



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



Portable Centrifuge

Model YR142-6



YR142-6 Portable Centrifuge

Kalstein's Low Speed Centrifuge is portable design and has low noise. It is popular with clinic, research and university.

Technical performance

- . Microcomputer control, brushless motor.
- . It is programmable and has LED display.
- . It has speed/RCF switch easily.
- . There are lid safety device, speed safety device errors detection.
- . ABS housing make the clean easy.
- . Silastic sealed ring comfort to GMP.US FDA certified.

Low-speed centrifuge installation requirements

- ✓ Low-speed centrifuge uses brushless DC motor,

maintenance-free.

- ✓ One-piece electronic lock: to avoid injury to the experimental environment and experimental personnel when the centrifuge is running.
- ✓ Intelligent control, touch panel, intuitive digital display, can directly set parameters such as rotation speed and time, and is convenient to use.
- ✓ The low-speed centrifuge body adopts high-quality steel structure, stainless steel centrifugal cavity, durable, easy to clean.
- ✓ The main application scope of low-speed centrifuges is: the new drug research and development unit tests the heart of SD rat blood, the drug purification, and the powder in the centrifugal liquid of the Chinese medicine decoction pieces.



#Letsgivingore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél: +33 467158849 / +33 680760710/ +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S



Portable Centrifuge Model YR142-6



Specification

Max Speed	4000r/min
Max RCF	1980xg
Max Capacity	15mlx6
Speed Accuracy	± 30r/min
Time Setting	1min~99min
Noise	<65dB(A)
Power supply	110V/60Hz, 220V/50Hz
Dimension(LxWxH)	310mmx260mmx220mm
Packing size(LxWxH)	370mmx320mmx280mm
Net weight	5 KG

Rotor

No.1



Rotor Capacity: 7mlx12 10mlx12
Max Speed: 4000r/min
Max RCF: 1980xg

No.2



Rotor Capacity: 15mlx6 15mlx8
Max Speed: 4000r/min
Max RCF: 1980xg



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél: +33 467158849 / +33 680760710/ +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:

 **Fundación
Jacinto Convit**
<https://www.jacintoconvit.org/>

 **ONE
TREE
PLANTED**

<https://onetreepanted.org>

 **Humatem**
Solidarité et coopération
biomédicale
<https://www.humatem.org/>

 **MANIAPURE
FOUNDATION USA**
<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed

The 3 articles of the week!

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

1
Read this ARTICLE NOW



How does the climate affect your microscope?

microscope is an instrument that allows observing and measuring objects that are too small to be seen by the naked eye. This is achieved through a complex system of lenses and mirrors, which form

Kalstein Blog

Article of the week:
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

2
Read this ARTICLE NOW



What is the importance of the laboratory refrigerator?

refrigerator is a device that is one of the most common laboratory equipment. Its purpose is to maintain, in a controlled environment (refrigerated space), various fluids and substances

Kalstein Blog

Article of the week:
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



MONTPELLIER - FRANCE
HEADQUARTER
Optimum Business Center,
450 Rue Baden Powell,
34000 Montpellier - France.

Fax: +33 (0) 4 67 15 88 49
Tif: +33 (0) 6 80 76 07 10
sales@kalstein.eu

<https://kalstein.net/>
<https://www.kalstein.fr/>

Official Websites for Latin America.

