



Kalstein®

A different accompaniment, at your service

All rights reserved © KALSTEIN France S. A. S.,

VETERINARY SECTOR

VETERINARY X-RAY



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Veterinary X-Ray Veterinary High-frequency Mobile C-arm System Model YR06165



Product Description:

Orthopedics, General Surgery, Pain management, Traumatology and so on for Veterinary..

Features:

- ✓ Integration digital acquisition and processing system, more convenient for operation
- ✓ High-quality combined HF-HV X-ray generator, greatly reducing the X-ray irradiation
- ✓ Electric rotatable adjustable collimator
- ✓ Unique hand-held remote controller, more convenient operation
- ✓ Graphical colorfully LCD touch screen interface, convenient operation, fashionable design
- ✓ Toshiba image intensifier is employed, stable and reliable in quality, high-definition image and display system
- ✓ With the perspective KV, MA automatic tracking feature, the image brightness and definition are automatically set at optimal function.



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 2 Rue Jean Lantier • 75001 Paris •
 +33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
 KALSTEIN FRANCE, S. A. S



KALSTEIN® Veterinary X-Ray
Veterinary High-frequency Mobile C-arm System
Model YR06165



Technical Description:

Model		YR06165
Electrical Performance	High Frequency Inverter Output Power	Output Power: 5kW Inverter Frequency:40kHz
	Automatic Fluoroscopy	Tube Voltage:40kV~120kV automatic adjustment Tube Current:0.3mA~4mA automatic adjustment
	Manual Fluoroscopy	Tube Voltage:40kV~120kV continuous adjustment Tube Current:0.3mA~4mA continuous adjustment
	Pulse Fluoroscopy	Tube Voltage:40kV~120kV continuous adjustment Tube Current:0.3mA~30mA continuous adjustment
	Photography Tube Voltage, mA, mAs	40 kV~120 kV 25-100mA 1-180mAs
	Film Holder Size	200mm×250mm(8"×10") or 250mm×300mm (10"×12")
X-Ray Tube	X-Ray Tube Special for High Frequency	Fixed anode Dual-focus, small focus: 0.6, large focus: 1.8mm Thermal capacity: 35.5KJ (47KHU)
Video System	Image Intensifier	Toshiba 9" Image Intensifier 3-fireds
	CCD Video camera	Medical Mega-pixels digital CCD camera
	24" LCD monitor	1920×1080, 60Hz
	Digital workstation	Storage without damage, multi-images display, image W/L real time adjustment, grayscale switch, region of interest balance, GAMMA correction, flip, noise reduction, intensifying, smoothing, sharpening, compression, zoom in, measurement, remark, picture-text typeset print, expert template, DICOM image transmit, DICOM image print, cine loop, image imprinting, Worklist registration etc.
Structure	Directive wheel	Direction wheel 360° rotation, main wheel ±90° rotation
	C-arm	Electric Ascending & Descending of pillar movement: 400mm Forward and Backward Movement: 200mm Revolution around Horizontal Axis: ±180° Revolution around Vertical Axis: ±15° Slip on orbit: 120°(+90°~ -30°) Monitor Rotation: 300° SID: 960mm Opening: 760mm Arc depth: 640mm



#Letsgivingmore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 2 Rue Jean Lantier • 75001 Paris •
 +33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



Fundación
Jacinto Convit

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale

<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed
The 3 articles of the week!
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW
How does the climate affect your microscope?
microscopes in a laboratory are one of the most sensitive pieces of equipment that require to maintain in a controlled environment (refrigerated space), without fluctuations in humidity and temperature.
Kalstein Blog
Article of the week:
How does the climate affect your microscope?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW
What is the importance of the laboratory refrigerator?
Refrigerator in a laboratory is one of the most sensitive pieces of equipment that require to maintain in a controlled environment (refrigerated space), without fluctuations in humidity and temperature.
Kalstein Blog
Article of the week:
What is the importance of the laboratory refrigerator?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Alliances
23 views • 3 years ago



Kalstein's Microscopes
445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER
2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89
Cel: +33 (0) 6 80 76 07 10
sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

