



A different accompaniment, at your service



 **Kalstein**®

All Right Reserved © KALSTEIN France S. A. S

LABORATORY SECTOR
ANALYZERS

www.kalstein.net



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





Benefits & Assistance

At Kalstein France we care about the full satisfaction of our customers, that's why we offer value-added services of the highest level based on our expertise.

Quick Response

Our team is always available to answer all your questions or queries, in order to support you in any situation...

Kalstein+ Platform

All of this would be presented in an intuitive and easy-to-use interface, providing a personalized and satisfying shopping experience for each customer.

Technical Support

Enjoy personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein manuals and articles, special catalogs and video tutorials.

Kalstein worldwide

With more than 25 years growing with our customers, Kalstein's multi-format and modern content is now present in more than 10 countries and growing.

#Letsgivemore

Thanks to your purchase, a donation is made to a non-profit foundation that fights cancer and helps vulnerable communities.

Delivery Logistics

We take care of all the logistics necessary for the shipment of your goods, whether by sea, land or air.

AI Model Recommendation

Model recommendation through AI is used to analyze user behavior, preferences and specific characteristics in order to recommend the most suitable product models for each customer.



KALSTEIN® Analyzers

5-Part Hematology Analyzer

Model YR06094



Product Description:

- ✓ The brand-new automated 5-part differentials hematology analyzer.
- ✓ Tri-angular laser scatter + Flow cytometry for 5-part WBC counting.
- ✓ Independent BASO counting channel.
- ✓ Whole Blood, Peripheral Blood, Pre-diluted blood.
- ✓ Throughput with 60 samples / hour.
- ✓ 29 parameters (include 4 research parameters).
- ✓ 3 reagents with original QC and Calibrators.

Principles:

- ✓ Electrical impedance, Colorimetry; Laser scattering method

Parameters:

- ✓ WBC, BAS#, BAS%, NEU#, NEU%, EOS#, EOS%, LYM#, LYM%, MON#, MON%, RBC, HGB, HCT, MCV, MCH, MCHC, RDW-CV, RDW-SD, PL T, MPV, PDW, PCT, P-LCR, P-LCC; AL Y#, AL Y%, LIC#, LIC%;
- ✓ WBC Histogram, RBC Histogram, PL T Histogram; 2D DIFF Scatter gram x 3; 3D DIFF Scatter gram x 1;

Reagent:


- ✓ RD-50 Diluent, RL-50D lyse, RL-50H lyse

Throughput :

- ✓ 60 S/H

System:

- ✓ Integrated windows system displayed by 10.4 inches color touch screen.
- ✓ Bi-directional LIS communication



Testing Mode:

- ✓ CBC, CBC+DIFF)

Sample Consumption:

- ✓ Whole Blood: 15ul (CBC+DIFF/ CBC)
- ✓ Pre-diluted 20ul(CBC+DIFF/ CBC)

Optional external devices:

- ✓ Bar-code reader, Mouse, Keyboard, External printer

Sample Types:

- ✓ Whole Blood, Peripheral Whole Blood, Pre-diluted blood



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® Analyzers

5-Part Hematology Analyzer

Model YR06094



Link Ports.

- ✓ LAN; USB x 4

Working Environment:

- ✓ 10°C ~30°(; humidity :S 85%

Data Storage:

- ✓ 100000 results

Power:


- ✓ 110V-240V; 50Hz/60Hz; 400VA

Printer:

- ✓ Internal thermal printer

Size (mm):

- ✓ 365x 470x 525;Weight: 31Kg



Test	Linear range	Carry-over	CV(whole blood)
WBC	0-500×10 ⁹ /L	≤0.5%	≤2.0%
RBC	0-8×10 ¹² /L	≤0.5%	≤ 1.5%
HGB	0-250g/L	≤0.5%	≤ 1.0%
PLT	0-5000×10 ⁹ /L	≤1.0%	≤ 4.0%



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





FOR ALL KALSTEIN CUSTOMERS

JOIN KALSTEIN+ NEW PLATFORM!

TECHNOLOGICAL SOLUTIONS!



Artificial Intelligence

State-of-the-art tool (integrated into all web pages) necessary for the rapid preparation of tenders or ordinary requests.



Technological Device

A cutting-edge device designed to seamlessly connect your lab equipment with our mobile app on Kalstein phones..



K+ Platform

The distributor client can build their business and make their quotes automatically with shipping costs to the final destination always with the identity of their company.

BUSINESS SOLUTIONS!



Marketing Service

Creation of content for the Blog and informative guides for SEO positioning



Dropshipping

For the reception and preparation of the order on behalf of your company and sent to the end customer connected to Kalstein tracking with product customization.



Special Discounts

Up to 41% plus 4% on shipping

INTRODUCING THE NEW ERA OF KALSTEIN INTELLIGENT EQUIPMENT

We are pleased to announce the development of a revolutionary addition to our product range: a cutting-edge device designed to seamlessly connect your laboratory equipment with our mobile app on Kalstein phones.

CONTACT US

(+33) 01 70 39 26 50 / (+33) 7 60 83 92 83

Mail: sales.department@kalstein.eu

www.kalstein.net

📍 5 rue de Castiglione 75001 Paris – France.

CONNECT WITH US ON KALSTEIN+

DO YOU NEED ASSISTANCE?



HOW TO REQUEST A VIDEO CALL?

At Kalstein+, we are committed to providing you with the best possible service, adapting to your needs in a personalized and efficient manner. That's why we cordially invite you to join us on a video call.



+(33) 7 77 66 24 64



HOW CAN YOU DOWNLOAD THE APPLICATION?

Download the Kalstein+ app now and take the Kalstein experience with you wherever you go! Don't miss the opportunity to stay up to date with the latest news and special offers.



Download it on
the Play Store



Download it on
the App Store

CONTACT US

+(33) 01 70 39 26 50 / +(33) 7 60 83 92 83
Mail: sales.department@kalstein.eu
www.kalstein.net

📍 5 rue de Castiglione 75001 Paris – France.