



A different kind of support at your service



Kalstein®

All Rights Reserved © KALSTEIN France S. A. S.,

**LABORATORY SECTOR
REFRIGERATORS
AND FREEZERS**



All rights reserved KALSTEIN France S.A.S.
5 Rue Castiglione 75001 Paris - France
+33 (0) 1 70 39 26 50 / +33 (0) 6 80 76 07
<https://kalstein.net>
KALSTEIN FRANCE, S.A.S

www.kalstein.net/en/



Benefits & Assistance

At Kalstein France we care about the full satisfaction of our customers, that's why we offer value-added services of the highest level based on our expertise.



Kalstein+ Platform

All of this would be presented in an intuitive and easy-to-use interface, providing a personalized and satisfying shopping experience for each customer.



Quick Response

Our team is always available to answer all your questions or queries, in order to support you in any situation...



#Letsgivemore ♡

Thanks to your purchase, a donation is made to a non-profit foundation that fights cancer and helps vulnerable communities.



Technical Support

Enjoy personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein manuals and articles, special catalogs and video tutorials.



Delivery Logistics

We take care of all the logistics necessary for the shipment of your goods, whether by sea, land or air.



Kalstein worldwide

With more than 25 years growing with our customers, Kalstein's multi-format and modern content is now present in more than 10 countries and growing.



AI Model Recommendation

Model recommendation through AI is used to analyze user behavior, preferences and specific characteristics in order to recommend the most suitable product models for each customer.

KALSTEIN® Refrigerators and Freezers -30°C Forced Air Cooling Biomedical Freezer

Model YR05308-1



Product Description:

Applied for the storage of critical and temperature sensitive biological samples, laboratory products and medical products. Suitable for blood banks, hospitals, disease control centers, research laboratories.

Features:

- ✓ Microprocessor-controlled system
- ✓ Multiple alarm functions
- ✓ Low noise design
- ✓ Alarm record query function
- ✓ Forced air and auto defrost design
- ✓ Adjustable temperature range -10°C~-30°C



#Letsgivemore 

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
5 Rue Castiglione 75001 Paris – France
+33 (0) 1 70 39 26 50 / +33 (0) 6 80 76 07 10
<https://kalstein.net>
KALSTEIN FRANCE, S.A.S



KALSTEIN® Refrigerators and Freezers

-30°C Forced Air Cooling Biomedical Freezer

Model YR05308-1



Technical Specifications:

Model	YR05308-1	
Technical Data	Cabinet Type	Upright
	Climate Class	SN N
	Cooling Type	Forced air cooling
	Defrost Mode	Automatic defrost
	Refrigerant	HC
	Sound Level((dB(A))	44
Performance	Cooling Performance (°C)	-30
	Temperature Range (°C)	-10-30
Control	Controller	Microprocessor
	Display	LED
Electrical Data	Power Supply(V/Hz)	220-240/50
	Power(W)	650
	Electrical Current(A)	3:1
Dimensions	Capacity(L/Cu.Ft)	520/18.4
	Net/Gross Weight (approx) (kg)	210/252
	Interior Dimension (W×D×H) (mm)	610*755*1006
	Exterior Dimension (W×D×H) (mm)	848*1011*1780
Alarms	High/Low Temperature	Y
	Remote Alarm	Y
	Power Failure	Y
	Sensor Error	Y
	Low Battery	Y
	High Ambient Temperature	Y
	Door Ajar	Y
Accessories	Castor	Y
	Foot	Y
	Porthole	Y
	Drawers/Inner Doors	/
	USB Interface	Optional
	RS485	Y
	NFC	Optional
	Temperature Recorder	Optional



#Letsgivemore ♥

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
5 Rue Castiglione 75001 Paris – France
+33 (0) 1 70 39 26 50 / +33 (0) 6 80 76 07 10
<https://kalstein.net>
KALSTEIN FRANCE, S.A.S



#LETSGIVEMORE 


**WITH THE
ACQUISITION
OF KALSTEIN
EQUIPMENT**

TO WHICH YOU CONTRIBUTE:



**We present the
prototype of the
Kalstein phone**

Stay ahead of
trends with
our advanced
mobile
technology.



www.kalstein.net

www.kalstein.net/en/



PARIS - FRANCE

HEADQUARTERS

5 Rue Castiglione
75001 Paris - France

Tlf: +33 (0) 1 70 39 26 50

Cel: +33 (0) 6 80 76 07 10

Cel: +33 (0) 6 63 81 00 23

sales@kalstein.eu

<https://www.kalstein.eu/>

<https://www.kalstein.net/>

Follow us on our Social Networks!



[@kalsteinlatinoamerica](https://www.instagram.com/kalsteinlatinoamerica)



[@kalsteinlab](https://www.twitter.com/kalsteinlab)



[@kalsteinlabfr](https://www.instagram.com/kalsteinlabfr)

Here you will find our Tutorials!



[@kalsteinenglish4897](https://www.youtube.com/channel/UCkALsteinenglish4897)

