




LABORATORY SECTOR

REFRIGERATORS AND
FREEZER



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Refrigerators and Freezers**22 °C Blood Platelet Incubator****Model YR06028 // YR06029****Product Description:**

Kalstein blood platelet incubator is integrated with agitator. The agitation system has been designed for continuous quiet operation. It is widely used to store blood platelet (20°C ~ 24°C), biological reagent (2°C ~ 40°C) and chemical product (2°C ~ 40°C). The applications could be found in hospital, blood bank, research institute, laboratory and universities. The platelet incubator resembles a simple storage cabinet with multiple levels of horizontal flatbed storing units, with a heating element and refrigerating element incorporated in the storage space.

Features:

- ✓ Constant Temperature: 2°C of temperature fluctuations
- ✓ Thermal Printer: Real-time automatic printing of temp, can be optional USB port
- ✓ Universal Castor Wheels: Castors with brake
- ✓ High-Quality Steel: 304 Stainless steel chamber and 304 stainless steel oscillator
- ✓ Microprocessor control system: Continuously alarm for 72h after power failure. Microprocessor electronic temperature controller, LED temperature display, the display accuracy is 0.1°C.
- ✓ 5 Types Of Alarms: High and low temperature alarm, power failure alarm, Sensor malfunction alarm, Door Ajar, Power failure backup system alarm.

**#Letsgivemore** ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® Refrigerators and Freezers

22 °C Blood Platelet Incubator

Model YR06028 // YR06029



Technical Description:

Model		YR06028	YR06029
Specificaitons	Shelf/drawer QTY	Drawer / 5	Drawer / 10
	Cooling method	Forced air cooling	
	Defrost	No frost	
	Refrigerant	R134a	
	Power consumption(kWh/24h)	2.69	4.4
	Noise (db)	49	50
	Ambient temperature (°C)	10~32 °C	
	Temperature range (°C)	20~24 °C	
Cooling system	Compressor brand/QTY	EMBRACO/1	
	Sensor	NTC	
	Temperature controller	Microprocessor	
	Display	Digital display	
Power	Voltage/frequency (V/Hz)	220/50Hz	
	Power (W)	197W	280W
	Current (A)	0.9A	1.27A
Material	Inside material	SUS	
	Outside material	Colored Spray Steel	
	Insulation	PU	
Size	Capacity (L/cu.ft)	5	10
	Bloodbags(450ml)/Layers	10 /5 layers	20 /10 layers
	NT./GT. (kg)	65/87	85/108.5
	Exterior size(W*D*H)(mm)	520*640*1030	520*640*1290
	Interior size(W*D*H)(mm)	427*360*467	432*462*728
	20GP/40GP/40HQ	24/48/96	32/68/68
Alarm System	High/low temperature	Y	
	Power failure	Y	
	Sensor malfunction	Y	
	Low batery	/	
	Door ajar	Y	
	Power failure backup system(alarm)	72h	
	Caster	Y	
Accessories	Outer door / Type	1 (Foamed glass door)	
	Test hole	Optional	
	Light	UV Light	
	Temperature recorder	Optional	
Other	Certification	ISO, CE	



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale

<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed
The 3 articles of the week!
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW
How does the climate affect your microscope?
microscope is an instrument used for observing objects that are too small to be seen by the naked eye. This instrument is composed of two lenses, which form an image of the object being observed.
Kalstein Blog
Article of the week:
How does the climate affect your microscope?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW
What is the importance of the laboratory refrigerator?
Refrigerator is a device used for the storage of the most sensitive equipment that requires to maintain, in a controlled environment (refrigerated space), various fluids and substances.
Kalstein Blog
Article of the week:
What is the importance of the laboratory refrigerator?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Alliances

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER

2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89

Cel: +33 (0) 6 80 76 07 10

sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

