




Kalstein®

A different accompaniment, at your service

All rights reserved © KALSTEIN France S. A. S.,

VETERINARY SECTOR

VETERINARY ELECTROSURGICAL UNIT



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Veterinary Electrosurgical Unit
26-inch Mobile High-Definition Endoscopic Imaging System
Model YR06156



Product Description:

- ✓ Monitor/camera/light source and image recorder in a machine
- ✓ Smaller size/Lighter weight/Easy to move
- ✓ Portable power
- ✓ Utilization in variable places and fields
- ✓ Variable laying style/Adjust monitor position any time as required by the operator
- ✓ Compatible with Optics/Electronic/Fiber endoscopy

Features:

- ✓ Camera/Display/Illumination
- ✓ Image recorder/Storage/Replay/Snapshot
- ✓ Simultaneous second transmission
- ✓ Network video transmission
- ✓ Network telecast live/Live browse on network terminal



Technical Description:

Model		YR06156
Medical LED Light Source	Power	≥30W
	Color Temperature	6500K
	Illumination	≥700,000Lux
	Color Render Index	93
	Adjustment Mode	Stepless Adjustment
Integrated Control System	Control Method	Foot Switch/Remote Control
	Wireless Remote Control	Infrared code full function remote control
	Foot Switch	Record/Snapshot
Camera	Sensor	1/3inch Full HD CMOS
	Resolution	1920*1080
	Scanning Method	Progressive Scanning
	Scanning Frame Rate	60fps
	Minimum Illumination	1Lux
	Electronic Shutter	1/60 ~ 1/10000s
	White Balance	Manual/Auto
Back Light Compensation	Auto	



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 2 Rue Jean Lantier • 75001 Paris •
 +33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale

<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed

The 3 articles of the week!

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

1
Read this ARTICLE NOW



How does the climate affect your microscope?

microscope is an instrument that allows observing and measuring small objects. This is achieved through a series of lenses that focus light on the object being observed.

Kalstein Blog

Article of the week:
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

2
Read this ARTICLE NOW



What is the importance of the laboratory refrigerator?

Refrigerator is a device that is one of the most common equipment that is used to maintain, in a controlled environment (refrigerated space), various fluids and substances.

Kalstein Blog

Article of the week:
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...
23 views • 3 years ago



Kalstein's Microscopes
445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER
2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89
Cel: +33 (0) 6 80 76 07 10
sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

