




LABORATORY SECTOR

REFRIGERATORS AND FREEZERS



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Refrigerators and Freezers**316L 2-8°C Upright Medical Grade Refrigerator Pharmaceutical****Model YR05288****Product Description:**

Kalstein Medical-grade Pharmacy Refrigerators are designed with the unique needs of the pharmacy environment in mind and safely store vaccines, pharmaceuticals, chemotherapy and other products requiring safe cold storage. The digital display control panel ensure the temperature display accurately in 0.1°C.

Features:

- ✓ Electrically Heated Glass Door: Fog or condensation often occurs when opening and closing the door. However, this can be easily prevented which the electrically heated glass door, which quickly evaporate any moisture.
- ✓ High Quality Materials: The Glass-door series products are made of high-quality coated cold-rolled steel both inside and outside, which is antibacterial, anti-corrosive, and easy to clean.
- ✓ Green Earth: Sustainable U.S. EPA, SNAP, and EU F-Gas compliant natural refrigerants and foaming agents support sustainability initiatives. The cooling system uses R290 which are environmentally friendly, having no impact on ozone depletion and a very low Global Warming Potential (GWP) grade.
- ✓ Ultimate Uniformity: The cabinet maintain temperature within +/-3 degree throughout the entire compartment, ensuring confidence that items are stored at the correct temperature regardless of their placement in the chamber.

**#Letsgivemore** ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® Refrigerators and Freezers**316L 2-8°C Upright Medical Grade Refrigerator Pharmaceutical****Model YR05288**

Model	YR05288	YR05289	YR05289-1	YR05289-2
Shelf/drawer QTY	Wired shelf/5			Wired shelf/4
Cooling system	Forced air			
Defrost	Frost-free			
Refrigerant	R290			R600a
Noise(db)	47	50	46	
Ambient temperature (C)	10~32°C			
Temperature range (C)	2~8°C			
Compressor QTY	1			
Sensor	NTC			
Temperature controller	Microprocessor			
Display	LED digital display			
Voltage(V)	187~242		187~242C	187~242
Frequency (Hz)	50Hz/60Hz			
Power Consumption (kWh/24h)	2.2	3.2	2.4	1.6
Power (W)	480W	280W	260W	140W
Internal material	Coated steel			HIPS
External material	Coated steel			
Insulation	CFC-free PURF			

**#Letsgivemore**

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® Refrigerators and Freezers**316L 2-8°C Upright Medical Grade Refrigerator Pharmaceutical****Model YR05288**

Capacity (L/cu.ft)	316/11.15	416/14.69		226/7.98
NT./GT.(kg)	92.5/100	120/133	118/130	73/80
Exterior size(W*D*H)(mm)	660*555*1915	782*552*1927		625*594*1707
Interior size(W*D*H)(mm)	540*425*1380	685*429*1380		534*406*1150
20GP/40GP/40HQ	27/54/54	18/39/39		24/51/51
High/low temperature	Yes			
Power failure	Yes			
Sensor error	Yes			
Low battery	N/A			
Condenser fault	N/A			
Ambient temperature	N/A			
Door ajar	Yes	N/A	Yes	
Power failure backup	8 hours			
Castor	Yes			
Leveling feet	Yes			
Glycerine box/Qty	Yes/1			N/A
Illumination	LED			
Access port Qty/diameter	1/25mm			
USB Interface	Optional			
Remote alarm port	Yes			
Door type/Qty	Heating glass/1		Self-closing Solid/1	Heating glass/1
Heated glass	Yes		N/A	Yes
Chart recorder	N/A	Optional	N/A	Optional
Data logger	Optional			
AISI 304 interior	Optional			N/A

**#Letsgivemore**

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale


<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed

The 3 articles of the week!

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

1
Read this ARTICLE NOW



How does the climate affect your microscope?

Kalstein Blog

Article of the week:
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

2
Read this ARTICLE NOW



What is the importance of the laboratory refrigerator?

Kalstein Blog

Article of the week:
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...
23 views • 3 years ago



Kalstein's Microscopes
445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER
2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89
Cel: +33 (0) 6 80 76 07 10
sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

