




LABORATORY LINES

BALANCES



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.




Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



Analytical balance

Model YR05779 // YR05781



Product Description

- ✓ Aluminum shell/Large space windproof cover
- ✓ LED (green)/LCD
- ✓ Tare/Count function
- ✓ Overload alarm/level indicator
- ✓ Option: Under weighing/Printer

Balances

The scale is a piece of equipment that measures the mass of a specific body, providing exact data on its weight. Among the many existing models, the analytical balance stands out, which has a high power of rigor, since with it a great precision can be obtained.

The precision is the number of decimal places used to express what is measured, while the accuracy indicates the degree of approximation between the measured value and the real value that is intended to be measured. On the other hand, if we refer to a series of measurements, the precision is the degree of grouping or repeatability of the measured values or dispersion of the results, while the accuracy represents the degree of approximation between the mean of the measured values and the real value.

[See products](#)



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél: +33 467158849 / +33 680760710/ +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S



Analytical balance
Model YR05779 // YR05781



Technical Description

Model	YR05779	YR05780	YR05781
Ability	110g	210g	310g
Readability	0.001g		
minimum weighing	0.004g		
repeatability	±0.002g		
linearity	±0.003g		
stable time	≤3S		
pan size	Ø90mm		
appearance size	340x215x350mm		
Draft Shield Size	240x190x265mm		
Communication	RS232 (standard)		
Calibration	external calibration		
Cal.Weight	100 grams	200g	
Option	Printer/Under weighing		
packing size	Inner box:415x285x430mm Outer box:475x315x485mm		
Gross weight	7.5kg		



#Letsgivingmore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 Optimum Business Center 450 Rue Baden Powell,
 • 34000 Montpellier, France.
 Tit: +33 467158849 / +33 680760710/ +33 663810023
<https://kalstein.eu>
 KALSTEIN FRANCE, S.A. S



Analytical balance

Model YR05779 // YR05781

Types of laboratory balances that a laboratory may need



Precision scale

- metallic base.
- Round stainless pan.
- LED (green).
- Tare/Count function.
- Overload alarm.
- Option: Under weighing/Interface/Printer/Windproof cover.

[see category](#)

analytical balances

- electromagnetic sensor.
- Aluminum casing Large space windproof cover.
- LCD (white backlight with black lettering).
- Tare/Counting/Unit conversion (g/mg/ct/oz) function.
- Sensitivity/Speed settings/Percentage.
- Overload alarm/Fault alarm/Level indicator.
- Option: Under weighing/Printer.

[see category](#)



#Letsgivingmore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél: +33 467158849 / +33 680760710 / +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S





#LETSGIVEMORE

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**


YOU MAKE A CONTRIBUTION TO:



<https://onetreepanted.org>



Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed

The 3 articles of the week!

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW

How does the climate affect your microscope?

microscope is an instrument used for observing small objects that are not visible to the naked eye. This is achieved through a series of lenses that magnify the object.

Kalstein Blog

Article of the week:
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW

What is the importance of the laboratory refrigerator?

Refrigerator is a device that is one of the most important equipment that is used to maintain, in a controlled environment (refrigerated space), various fluids and substances.

Kalstein Blog

Article of the week:
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



MONTPELLIER - FRANCE
HEADQUARTERS

Optimum Business Center,
450 Rue Baden Powell,
34000 Montpellier - France.

Fax: +33 (0) 4 67 15 88 49

Tlf: +33 (0) 6 80 76 07 10

sales@kalstein.eu

<https://kalstein.net/>
<https://www.kalstein.fr/>

Official Websites for Latin America.

