



A different accompaniment, at your service



 **Kalstein**®

All Right Reserved © KALSTEIN France S. A. S

LABORATORY SECTOR
ANALYZERS

www.kalstein.net/en/



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Analyzers

Automated Chemistry Analyzer

Model YR06073



Product Description:

- ✓ Automatic, Random Access
- ✓ Up to 100 tests per hour
- ✓ Reagent pre-heating, liquid level detection
- ✓ Labor saving, simple programming and real walk-away operation
- ✓ Micro-volume for sample and reagent
- ✓ User-friendly software, simplicity and flexibility for operation
- ✓ Multi-language software supported.



Technical Description:

Model:	YR06073
Performance:	up to 100 tests per hour
Test to perform:	chemical clinic, immunoturbidimetric
Analytical modes:	End of point, two points, kinetics
Programming:	Open system with user-defined programs
Sample position:	9 positions for samples
Reagent position:	Reagent position:
Refrigeration temperature:	4 – 15°C
Reaction cuvette positions:	90 cuvettes, 10 cuvettes / segment
System:	Reagent preheating Collision protection, liquid level detection and inventory control. Automatic internal and external system wash after each cycle, carry-over <0.1%
Printer:	Built-in thermal printer, optional external printer.
Optical and Measuring System:	Halogen lamp Absorbance Range: 0 – 3.500 Abs + 1% Wavelengths: 8 wavelengths from 340 – 700 nm Resolution: 0.0001 ABS



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



Fundación
Jacinto Convit

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale

<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed
The 3 articles of the week!
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW
How does the climate affect your microscope?
microscopes in a laboratory are one of the most sensitive pieces of equipment that require to maintain in a controlled environment (refrigerated space), without fluctuations in humidity and temperature.
Kalstein Blog
Article of the week:
How does the climate affect your microscope?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW
What is the importance of the laboratory refrigerator?
Refrigerator in a laboratory is one of the most sensitive pieces of equipment that require to maintain in a controlled environment (refrigerated space), without fluctuations in humidity and temperature.
Kalstein Blog
Article of the week:
What is the importance of the laboratory refrigerator?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...
23 views • 3 years ago



Kalstein's Microscopes
445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER
2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89
Cel: +33 (0) 6 80 76 07 10
sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

