



A different kind of support at your service



 **Kalstein**®

All Rights Reserved © KALSTEIN France S. A. S.,

LABORATORY SECTOR
pH METER



All rights reserved KALSTEIN France S.A.S.
5 Rue Castiglione 75001 Paris - France
+33 (0) 1 70 39 26 50 / +33 (0) 6 80 76 07
<https://kalstein.net>
KALSTEIN FRANCE, S.A.S

www.kalstein.net/en/



Benefits & Assistance

At Kalstein France we care about the full satisfaction of our customers, that's why we offer value-added services of the highest level based on our expertise.

Quick Response

Our team is always available to answer all your questions or queries, in order to support you in any situation...

Technical Support

Enjoy personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein manuals and articles, special catalogs and video tutorials.

#Letsgivemore

Thanks to your purchase, a donation is made to a non-profit foundation that fights cancer and helps vulnerable communities.

AI Model Recommendation

Model recommendation through AI is used to analyze user behavior, preferences and specific characteristics in order to recommend the most suitable product models for each customer.

Kalstein+ Platform

All of this would be presented in an intuitive and easy-to-use interface, providing a personalized and satisfying shopping experience for each customer.

Kalstein worldwide

With more than 25 years growing with our customers, Kalstein's multi-format and modern content is now present in more than 10 countries and growing.

Delivery Logistics

We take care of all the logistics necessary for the shipment of your goods, whether by sea, land or air.



KALSTEIN® Ion Meter Benchtop pH / Ion Meter

Model YR01826



Product Description:

- ✓ This sophisticated measurement device boasts a 5-inch LCD screen with a 150° viewing angle, ensuring clear readability from up to 5 meters away. It features a built-in microprocessor chip that provides intelligent functions such as automatic calibration, Automatic Temperature Compensation (ATC), data storage, clock display, USB output, function settings, automatic power-off, and a low voltage alarm.
- ✓ Equipped with two USB interfaces, the device supports PC connections, battery charging, and data transfer. The 3.5mm electrode interface simplifies electrode connection. The device can automatically recognize different electrodes, and those with built-in chips can store and retrieve calibration data upon connection.
- ✓ It supports the automatic recognition of 15 buffers with options for Europe & USA, NIST, and China, and allows for 1, 2, or 3-point calibration. The device can store up to 200 sets of data, including readings, GLP data, date, and time. Its rechargeable battery provides up to 30 hours of use, and the device will automatically power off after 10 minutes of inactivity.



#Let'sgivemore ♥

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
5 Rue Castiglione 75001 Paris – France
+33 (0) 1 70 39 26 50 / +33 (0) 6 80 76 07 10
<https://kalstein.net>
KALSTEIN FRANCE, S.A.S



KALSTEIN® Ion Meter Benchtop pH / Ion Meter

Model YR01826



Technical Specifications:

Model	YR01826	
ph	Measuring range	(-2~19.999)pH
	Resolution	0.1/0.01pH
	Accuracy	±0.02 pH
	Input current	≤2×10 ⁻¹² A
	Input impedance	21×10 ⁻¹² Ω
	Stability	±0.002 pH/3h
	Temp. compensation range	(0~100) °C, auto/manual
mV	Measuring range	-1999.9mV~0~1999.9mV
	Resolution	0.1mV
	Accuracy	±0.03% FS
Temperature	Measuring range	(-10~110) C
	Resolution	0.1°C
	Accuracy	5~ 60°C: range: ±0.4°C; other range: ±0.8°C
Calibration		1, 2, 3 point
Other Technical Parameters	Display	5-Inch white backlit LCD screen
	Data storage	200 sets
	Power	DC 5V/1A
	Communication connector	Micro USB & Standard USB
	Dimension/N.W.	180 x 130 x 14mm/280g
Working Condition	Environment Temp.	535 °C
	Environmental	≤85%
	IP Grade	IP54
Meter Kits	1) pH meter	
	2) Electrode holder	
	3) Two-in-one electrode that consists of pH combination electrode and temperature probe.	
	4) pH buffer solution	
	5) Power adapter	



#Letsgivemore ♥

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
5 Rue Castiglione 75001 Paris – France
+33 (0) 1 70 39 26 50 / +33 (0) 6 80 76 07 10
<https://kalstein.net>
KALSTEIN FRANCE, S.A.S





FOR ALL KALSTEIN CUSTOMERS

JOIN KALSTEIN+ NEW PLATFORM!

TECHNOLOGICAL SOLUTIONS!



Artificial Intelligence

State-of-the-art tool (integrated into all web pages) necessary for the rapid preparation of tenders or ordinary requests.



Technological Device

A cutting-edge device designed to seamlessly connect your lab equipment with our mobile app on Kalstein phones..



K+ Platform

The distributor client can build their business and make their quotes automatically with shipping costs to the final destination always with the identity of their company.

BUSINESS SOLUTIONS!



Marketing Service

Creation of content for the Blog and informative guides for SEO positioning



Dropshipping

For the reception and preparation of the order on behalf of your company and sent to the end customer connected to Kalstein tracking with product customization.



Special Discounts

Up to 41% plus 4% on shipping

INTRODUCING THE NEW ERA OF KALSTEIN INTELLIGENT EQUIPMENT

We are pleased to announce the development of a revolutionary addition to our product range: a cutting-edge device designed to seamlessly connect your laboratory equipment with our mobile app on Kalstein phones.

CONTACT US

(+33) 01 70 39 26 50 / (+33) 7 60 83 92 83

Mail: sales.department@kalstein.eu

www.kalstein.net

📍 5 rue de Castiglione 75001 Paris – France.

CONNECT WITH US ON KALSTEIN+

DO YOU NEED ASSISTANCE?



HOW TO REQUEST A VIDEO CALL?

At Kalstein+, we are committed to providing you with the best possible service, adapting to your needs in a personalized and efficient manner. That's why we cordially invite you to join us on a video call.



+ (33) 7 77 66 24 64



HOW CAN YOU DOWNLOAD THE APPLICATION?

Download the Kalstein+ app now and take the Kalstein experience with you wherever you go! Don't miss the opportunity to stay up to date with the latest news and special offers.



Download it on
the Play Store



Download it on
the App Store

CONTACT US

(+33) 01 70 39 26 50 / (+33) 7 60 83 92 83

Mail: sales.department@kalstein.eu

www.kalstein.net

📍 5 rue de Castiglione 75001 Paris – France.