



A different kind of support at your service



Kalstein®

All Rights Reserved © KALSTEIN France S. A. S.,

LABORATORY SECTOR
COD / BOD ANALYZER



All rights reserved KALSTEIN France S.A.S.
5 Rue Castiglione 75001 Paris - France
+33 (0) 1 70 39 26 50 / +33 (0) 6 80 76 07
<https://kalstein.net>
KALSTEIN FRANCE, S.A.S

www.kalstein.net/en/



Benefits & Assistance

At Kalstein France we care about the full satisfaction of our customers, that's why we offer value-added services of the highest level based on our expertise.

Quick Response

Our team is always available to answer all your questions or queries, in order to support you in any situation...

Technical Support

Enjoy personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein manuals and articles, special catalogs and video tutorials.

#Letsgivemore

Thanks to your purchase, a donation is made to a non-profit foundation that fights cancer and helps vulnerable communities.

AI Model Recommendation

Model recommendation through AI is used to analyze user behavior, preferences and specific characteristics in order to recommend the most suitable product models for each customer.

Kalstein+ Platform

All of this would be presented in an intuitive and easy-to-use interface, providing a personalized and satisfying shopping experience for each customer.

Kalstein worldwide

With more than 25 years growing with our customers, Kalstein's multi-format and modern content is now present in more than 10 countries and growing.

Delivery Logistics

We take care of all the logistics necessary for the shipment of your goods, whether by sea, land or air.



KALSTEIN® Analyzers

COD Analyzer Portable

Model YR06296



Product Description:

- ✓ This device features a 3-inch color screen with a user-friendly graphical interface. It provides direct concentration readings and has an automatic power-off function.
- ✓ It includes a built-in standard curve, allowing users to create custom curves for various measurement requirements. The device supports up to 7-point calibration and includes a calibration restoring function to revert to factory calibration data.
- ✓ With a large-capacity storage function, it can store up to 5000 sets of test data. Additionally, the rechargeable battery lasts for 10 hours, and the auto power-off time can be customized.

Technical Specifications:

Model	YR06296
Measurement Method	Fast digestion spectrophotometry
Measuring Range	10-10000mg/L
Accuracy	≤5%
Resolution	0.01
Light Source	Imported Silicon Photodiode
Display	3-inch color screen
Output	Micro USB
Power	Power adapter/Rechargeable lithium battery
Waterproof Level	IP54
Ambient Temperature	(5-40) °C
Ambient Humidity	Relative humidity ≤ 85% (No condensation)



#Letsgivemore

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
5 Rue Castiglione 75001 Paris – France
+33 (0) 1 70 39 26 50 / +33 (0) 6 80 76 07 10
<https://kalstein.net>
KALSTEIN FRANCE, S.A.S





FOR ALL KALSTEIN CUSTOMERS

JOIN KALSTEIN+ NEW PLATFORM!

TECHNOLOGICAL SOLUTIONS!



Artificial Intelligence

State-of-the-art tool (integrated into all web pages) necessary for the rapid preparation of tenders or ordinary requests.



Technological Device

A cutting-edge device designed to seamlessly connect your lab equipment with our mobile app on Kalstein phones..



K+ Platform

The distributor client can build their business and make their quotes automatically with shipping costs to the final destination always with the identity of their company.

BUSINESS SOLUTIONS!

✓ Marketing Service

Creation of content for the Blog and informative guides for SEO positioning

✓ Dropshipping

For the reception and preparation of the order on behalf of your company and sent to the end customer connected to Kalstein tracking with product customization.

✓ Special Discounts

Up to 41% plus 4% on shipping

INTRODUCING THE NEW ERA OF KALSTEIN INTELLIGENT EQUIPMENT

We are pleased to announce the development of a revolutionary addition to our product range: a cutting-edge device designed to seamlessly connect your laboratory equipment with our mobile app on Kalstein phones.

CONTACT US

(+33) 01 70 39 26 50 / (+33) 7 60 83 92 83

Mail: sales.department@kalstein.eu

www.kalstein.net

📍 5 rue de Castiglione 75001 Paris – France.

CONNECT WITH US ON KALSTEIN+

DO YOU NEED ASSISTANCE?



HOW TO REQUEST A VIDEO CALL?

At Kalstein+, we are committed to providing you with the best possible service, adapting to your needs in a personalized and efficient manner. That's why we cordially invite you to join us on a video call.



+ (33) 7 77 66 24 64



HOW CAN YOU DOWNLOAD THE APPLICATION?

Download the Kalstein+ app now and take the Kalstein experience with you wherever you go! Don't miss the opportunity to stay up to date with the latest news and special offers.



Download it on
the Play Store



Download it on
the App Store

CONTACT US

| (+33) 01 70 39 26 50 / (+33) 7 60 83 92 83

| Mail: sales.department@kalstein.eu

| www.kalstein.net

📍 5 rue de Castiglione 75001 Paris – France.