



A different accompaniment, at your service



 **Kalstein**®

All Right Reserved © KALSTEIN France S. A. S

**LABORATORY SECTOR**

**ANALYZERS**

[www.kalstein.net/en/](http://www.kalstein.net/en/)



All rights reserved © KALSTEIN France S. A. S.,  
2 Rue Jean Lantier • 75001 Paris •  
+33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
KALSTEIN FRANCE, S. A. S





# Kalstein®

A different accompaniment, at your service

OUR SERVICES

## Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



### Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



### Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



### #Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



### Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



### Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



### Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.





**KALSTEIN® Analyzers****Coagulation Analyzer****Model YR06081****Product Description:**

- ✓ Advanced scattered light principle and percentage analysis technology
- ✓ 4 durable LED detectors
- ✓ Large LCD and touch screen
- ✓ Low consumption of reagents and consumables
- ✓ QC and calibration program incorporated
- ✓ Multiform result output
- ✓ Auto Self-check when power on
- ✓ Electrically-linked pipette (optional)

**Technical Description:**

<b>Model</b>	YR06081
<b>Principle</b>	Scattered light with percentage analysis
<b>Test Channel</b>	4 channels
<b>Light Source</b>	Durable LED detector
<b>Test speed</b>	Normally 20-60 secs, Maximum lime up to 600 secs
<b>Sample position</b>	24 positions
<b>Reagent position</b>	6 positions (1 magnetic mixer included)
<b>Minimum reagent cost†</b>	25ul for TT, FIB, PT, APTT
<b>Memory‡</b>	10000 test results
<b>Display</b>	6" LCD (320 x240, 256 gray scale)
<b>Input †</b>	Touch panel or external mouse
<b>Output‡</b>	External printer(optional)
<b>Not weight</b>	6.5KG
<b>Dimensions L × W × H (mm)</b>	410 × 310 × 160
<b>Power supply</b>	AC 110V/220V±10%, 50/60Hz

**#Let'sgivemore** ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
2 Rue Jean Lantier • 75001 Paris •  
+33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
KALSTEIN FRANCE, S. A. S





#LETSGIVEMORE

**WITH THE  
ACQUISITION  
OF A KALSTEIN  
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



<https://onetreeplanted.org>





**Articles Kalstein**  
(Social Networks and Official Websites)



**Kalstein Blog**  
Not just be informed, be fully informed

**The 3 articles of the week!**

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>

**1**  
Read this ARTICLE NOW



**How does the climate affect your microscope?**

microscope is an instrument that allows observing and measuring small objects. This is achieved through a series of lenses that focus light on the object being observed.

**Kalstein Blog**

**Article of the week:**  
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>

**2**  
Read this ARTICLE NOW



**What is the importance of the laboratory refrigerator?**

Refrigerator is a device that is one of the most common equipment that is used to maintain, in a controlled environment (refrigerated space), various fluids and substances.

**Kalstein Blog**

**Article of the week:**  
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>

**Kalstein Training Videos**  
(Youtube Channel and Official Websites)



**Kalstein France - Jacinto Convit Foundation / Allianc...**  
23 views • 3 years ago



**Kalstein's Microscopes**  
445 views • 1 week ago



PARIS - FRANCE  
**HEADQUARTER**  
2 Rue Jean Lantier  
75001 Paris - France

**Tlf:** +33 (0) 01 78 95 87 89  
**Cel:** +33 (0) 6 80 76 07 10  
[sales@kalstein.eu](mailto:sales@kalstein.eu)

<https://kalstein.eu/>  
<https://www.kalstein.fr/>

