




A different kind of support at your service



Kalstein®

All Rights Reserved © KALSTEIN France S. A. S.,

MEDICAL SECTOR OPERATION CEILING PENDANT



All rights reserved KALSTEIN France S.A.S.
5 Rue Castiglione 75001 Paris - France
+33 (0) 1 70 39 26 50 / +33 (0) 6 80 76 07
<https://kalstein.net>
KALSTEIN FRANCE, S.A.S

www.kalstein.net/en/



Benefits & Assistance

At Kalstein France we care about the full satisfaction of our customers, that's why we offer value-added services of the highest level based on our expertise.

Quick Response

Our team is always available to answer all your questions or queries, in order to support you in any situation...

Technical Support

Enjoy personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein manuals and articles, special catalogs and video tutorials.

#Letsgivemore

Thanks to your purchase, a donation is made to a non-profit foundation that fights cancer and helps vulnerable communities.

AI Model Recommendation

Model recommendation through AI is used to analyze user behavior, preferences and specific characteristics in order to recommend the most suitable product models for each customer.

Kalstein+ Platform

All of this would be presented in an intuitive and easy-to-use interface, providing a personalized and satisfying shopping experience for each customer.

Kalstein worldwide

With more than 25 years growing with our customers, Kalstein's multi-format and modern content is now present in more than 10 countries and growing.

Delivery Logistics

We take care of all the logistics necessary for the shipment of your goods, whether by sea, land or air.



KALSTEIN® Operation Ceiling Pendant Double Arm Surgical Pendant (Electrical)



Model YR02086



Product Description:

- ✓ The installation uses a ceiling suspension method, making the equipment easily accessible for medical staff.
- ✓ It is suitable for operating rooms, emergency rooms, intensive care units, and recovery rooms.
- ✓ The pendant cantilever is made from imported aluminum alloy, treated to prevent oxidation and rust.
- ✓ Optional friction or pneumatic brakes and German-standard gas terminals are available.
- ✓ Customers can choose from various gas outlet adaptors, compatible with German, British, American, and Australian standards.
- ✓ The system supports suspensible infusion bags, bottles, or pumps, offering multiple configurations to meet your needs. Different adaptors for friction and pneumatic brakes are also available, following German standards.



#Letsgivingmore ♥

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
5 Rue Castiglione 75001 Paris – France
+33 (0) 1 70 39 26 50 / +33 (0) 6 80 76 07 10
<https://kalstein.net>
KALSTEIN FRANCE, S.A.S



KALSTEIN® Operation Ceiling Pendant

Double Arm Surgical Pendant (Electrical)



Model YR02086

Technical Specifications:

Model	YR02086
Technical Parameters	HFP-DS240/380 double arm surgical pendant
Arm rotating radius	1000mm
Arm rotating angle	340°
The length of the instrument column	1 m
Braking system	Electromagnetic brake and damping brake
Maximum load capacity	220kg
Rails for attachments along the entire length of the column	Available
Valve for connecting surgical pneumatic instruments:	1
Standard configuration DIN standard	2 oxygen(O2)
	2 pcs Air (3-5 bar)
	2 CO2
	2 vacuum suction(VAC)
	1 AGSS
	12 Electrical sockets
	4 RJ 45 network socket
	6 grounding terminal
Annex	2 Shelf
	1 drawer
	1 Universal monitor holder
	3 infusion stand
Accessories	1 Oxygen flow meter with humidifier (autoclavable)
	1 waste collector complete with secretion collection tank (1 litre, autoclavable)



#Let'sgive more

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
5 Rue Castiglione 75001 Paris – France
+33 (0) 1 70 39 26 50 / +33 (0) 6 80 76 07 10
<https://kalstein.net>
KALSTEIN FRANCE, S.A.S





FOR ALL KALSTEIN CUSTOMERS

JOIN KALSTEIN+ NEW PLATFORM!

TECHNOLOGICAL SOLUTIONS!



Artificial Intelligence

State-of-the-art tool (integrated into all web pages) necessary for the rapid preparation of tenders or ordinary requests.



Technological Device

A cutting-edge device designed to seamlessly connect your lab equipment with our mobile app on Kalstein phones..



K+ Platform

The distributor client can build their business and make their quotes automatically with shipping costs to the final destination always with the identity of their company.

BUSINESS SOLUTIONS!

✓ Marketing Service

Creation of content for the Blog and informative guides for SEO positioning

✓ Dropshipping

For the reception and preparation of the order on behalf of your company and sent to the end customer connected to Kalstein tracking with product customization.

✓ Special Discounts

Up to 41% plus 4% on shipping

INTRODUCING THE NEW ERA OF KALSTEIN INTELLIGENT EQUIPMENT

We are pleased to announce the development of a revolutionary addition to our product range: a cutting-edge device designed to seamlessly connect your laboratory equipment with our mobile app on Kalstein phones.

CONTACT US

(+33) 01 70 39 26 50 / (+33) 7 60 83 92 83

Mail: sales.department@kalstein.eu

www.kalstein.net

📍 5 rue de Castiglione 75001 Paris – France.

CONNECT WITH US ON KALSTEIN+

DO YOU NEED ASSISTANCE?



HOW TO REQUEST A VIDEO CALL?

At Kalstein+, we are committed to providing you with the best possible service, adapting to your needs in a personalized and efficient manner. That's why we cordially invite you to join us on a video call.



+ (33) 7 77 66 24 64



HOW CAN YOU DOWNLOAD THE APPLICATION?

Download the Kalstein+ app now and take the Kalstein experience with you wherever you go! Don't miss the opportunity to stay up to date with the latest news and special offers.



Download it on
the Play Store



Download it on
the App Store

CONTACT US

(+33) 01 70 39 26 50 / (+33) 7 60 83 92 83

Mail: sales.department@kalstein.eu

www.kalstein.net

📍 5 rue de Castiglione 75001 Paris – France.