



A different kind of support at your service



Kalstein®

All Rights Reserved © KALSTEIN France S. A. S.,

LABORATORY SECTOR

**ELECTROLYTE
ANALYZER**



All rights reserved KALSTEIN France S.A.S.
5 Rue Castiglione 75001 Paris - France
+33 (0) 1 70 39 26 50 / +33 (0) 6 80 76 07
<https://kalstein.net>
KALSTEIN FRANCE, S.A.S

www.kalstein.net/en/



Benefits & Assistance

At Kalstein France we care about the full satisfaction of our customers, that's why we offer value-added services of the highest level based on our expertise.

Quick Response

Our team is always available to answer all your questions or queries, in order to support you in any situation...

Technical Support

Enjoy personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein manuals and articles, special catalogs and video tutorials.

#Letsgivemore

Thanks to your purchase, a donation is made to a non-profit foundation that fights cancer and helps vulnerable communities.

AI Model Recommendation

Model recommendation through AI is used to analyze user behavior, preferences and specific characteristics in order to recommend the most suitable product models for each customer.

Kalstein+ Platform

All of this would be presented in an intuitive and easy-to-use interface, providing a personalized and satisfying shopping experience for each customer.

Kalstein worldwide

With more than 25 years growing with our customers, Kalstein's multi-format and modern content is now present in more than 10 countries and growing.

Delivery Logistics

We take care of all the logistics necessary for the shipment of your goods, whether by sea, land or air.



KALSTEIN® Electrolyte Analyzer

Electrolyte Analyzer

Model YR06391 // YR06395



Product Description:

- ✓ This device features a 7-inch color touch screen that facilitates interaction between the user and the machine, ensuring simple and intuitive operation. In addition, it offers low running costs, allowing a quick return on investment, as consumables are considerably cheaper than imported materials.
- ✓ With the ability to analyze up to six elements and nine parameters in a single sample, it provides a complete analysis in a single step. Automation of reagent use, sample handling and tube status simplify the testing process. In addition, it features power outage protection and can automatically store up to 10,000 test results.
- ✓ Using advanced ion-selective electrode (ISE) technology, it ensures accurate results, with automatic bubble filtration to avoid clogging and automatic correction of potential tracking to maintain instrument stability.



#Let'sgivemore 

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
5 Rue Castiglione 75001 Paris – France
+33 (0) 1 70 39 26 50 / +33 (0) 6 80 76 07 10
<https://kalstein.net>
KALSTEIN FRANCE, S.A.S



KALSTEIN® Electrolyte Analyzer

Electrolyte Analyzer

Model YR06391 // YR06395



Technical Specifications:

Model	YR06391	YR06392	YR06393	YR06394	YR06395
Analysis method	Ion selective electrode(ISE)				
Sample type	Serum, plasma, whole blood and dilute urine				
Measurement item	K、 Na、 Cl、 Ca、 pH、 Li、 Mg				
Measuring speed	≤40s				
Test Parameters	Na, K, Cl	Na, K, Cl, Ca, pH	Na, K, Li	Na, K, Cl, Li	Na, K, Cl, Ca, pH, Li



#Letsgivemore

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 5 Rue Castiglione 75001 Paris – France
 +33 (0) 1 70 39 26 50 / +33 (0) 6 80 76 07 10
<https://kalstein.net>
 KALSTEIN FRANCE, S.A.S



KALSTEIN® Electrolyte Analyzer

Electrolyte Analyzer

Model YR06396 // YR06406



Accessories:

Model	Description
YRA264	K Electrode
YRA265	K Electrode
YRA266	Cl Electrode
YRA267	Ca Electrode
YRA268	pH Electrode
YRA269	Li Electrode
YRA270	Mg Electrode
YRA271	Ref. Electrode
YRA272	Reagent Pack Calibration A 430ml; Calibration B 110ml; (For Model H900 only)
YRA273	Reference Electrode Internal Solution – 10ml
YRA274	Electrode Internal Solution – 3ml
YRA275	Activating Solution – 110ml
YRA276	QC Solution Level 2 (Water-based) – 110ml
YRA277	Cleaning Solution – 110ml
YRA278	Electrode Wash Solution (Protein Enzyme) -125mg (25mg/bottle x5)
YRA278-1	Electrode Wash Solution (Diluent) – 5ml (1 bottle)



#Letsgivemore

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
5 Rue Castiglione 75001 Paris – France
+33 (0) 1 70 39 26 50 / +33 (0) 6 80 76 07 10
<https://kalstein.net>
KALSTEIN FRANCE, S.A.S





FOR ALL KALSTEIN CUSTOMERS

JOIN KALSTEIN+ NEW PLATFORM!

TECHNOLOGICAL SOLUTIONS!



Artificial Intelligence

State-of-the-art tool (integrated into all web pages) necessary for the rapid preparation of tenders or ordinary requests.



Technological Device

A cutting-edge device designed to seamlessly connect your lab equipment with our mobile app on Kalstein phones..



K+ Platform

The distributor client can build their business and make their quotes automatically with shipping costs to the final destination always with the identity of their company.

BUSINESS SOLUTIONS!

✓ Marketing Service

Creation of content for the Blog and informative guides for SEO positioning

✓ Dropshipping

For the reception and preparation of the order on behalf of your company and sent to the end customer connected to Kalstein tracking with product customization.

✓ Special Discounts

Up to 41% plus 4% on shipping

INTRODUCING THE NEW ERA OF KALSTEIN INTELLIGENT EQUIPMENT

We are pleased to announce the development of a revolutionary addition to our product range: a cutting-edge device designed to seamlessly connect your laboratory equipment with our mobile app on Kalstein phones.

CONTACT US

(+33) 01 70 39 26 50 / (+33) 7 60 83 92 83

Mail: sales.department@kalstein.eu

www.kalstein.net

📍 5 rue de Castiglione 75001 Paris – France.

CONNECT WITH US ON KALSTEIN+

DO YOU NEED ASSISTANCE?



HOW TO REQUEST A VIDEO CALL?

At Kalstein+, we are committed to providing you with the best possible service, adapting to your needs in a personalized and efficient manner. That's why we cordially invite you to join us on a video call.



+ (33) 7 77 66 24 64



HOW CAN YOU DOWNLOAD THE APPLICATION?

Download the Kalstein+ app now and take the Kalstein experience with you wherever you go! Don't miss the opportunity to stay up to date with the latest news and special offers.



Download it on
the Play Store



Download it on
the App Store

CONTACT US

(+33) 01 70 39 26 50 / (+33) 7 60 83 92 83

Mail: sales.department@kalstein.eu

www.kalstein.net

📍 5 rue de Castiglione 75001 Paris – France.