



A different accompaniment at your service



 **Kalstein**®

All Right Reserved © KALSTEIN France S. A. S

LABORATORY SECTOR

CENTRIFUGES

www.kalstein.net/en/



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Centrifuges**Floor-standing High-Speed cold Centrifuge****Model YR0111-5****Product Description:**

It has micro-computer control, electronic lock and AC variable-conversion motor which delivers high torque.

Technical performance:

- ✓ The max speed is 12000rpm, and the largest capacity of fixed rotor is 4X2000ml.
- ✓ It has micro-computer control, electronic lock and AC variable-conversion motor which delivers high torque.
- ✓ Front LCD touch screen make operation conveniently in any installed position.
- ✓ Imported compressor unit with Non-CFC refrigerant R404a.
- ✓ The operating data and failure record can be stored to read and print via USB interface.
- ✓ Automatic rotor recognition make safety for multiple rotors choice.
- ✓ Up to 100 user programs storage.
- ✓ Multistep-stages program supplied meet more demands.
- ✓ It has self-locking, over-speed detection, over-temperature detection, imbalance detection that is automatic alarm device.
- ✓ It has 10 levels acceleration and 11 levels deceleration to be setting according to the different requirement. Slow down freely is available.
- ✓ 9 levels acceleration and 9 levels deceleration.
- ✓ Imported PC transparent bottles is suitable for high temperature sterilization.
- ✓ Silastic sealed ring comfort to GMP, US FDA certificated.



#Let'sgive more ♥

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél: +33 467 158849 / +33 690760710 / +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S



KALSTEIN® Centrifuges
Floor-standing High-Speed cold Centrifuge
Model YR0111-5



Technical Specifications :

Model	YR0111-5
Max speed	12000r/min
Max RCF	19830xg
Max capacity	4x2000ml
Speed accuracy	± 10r/min
Time setting	1s~99h 59min59s
Temperature setting	-20 °C ~ +40 °C
Temperature accuracy	±1 °C
Compressor unit	French Tecumseh Compressor Unit R404a
Noise	<65db(A)
Power supply	AC220V, 50Hz 30A
Consumption	5.5kW



#Letsgivemore ♥

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 Optimum Business Center 450 Rue Baden Powell,
 • 34000 Montpellier, France.
 Tel: +33 467 158849 / +33 690760710 / +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S



KALSTEIN® Centrifuges
Floor-standing High-Speed cold Centrifuge
Model YR0111-5



 <p>No.1</p>	<p>Capacity: 8x100ml/85ml Max Speed: 12000r/min Max RCF: 19830xg Adaptor: 10ml、20ml</p>	 <p>No.2</p>	<p>Capacity: 6x1000ml Max Speed: 7000r/min Max RCF: 11600xg Adaptor: 100ml、250ml、300ml、500ml</p>
 <p>No.3</p>	<p>Capacity: 4x2000ml Max Speed: 6000r/min Max RCF: 8410xg Adaptor: 100ml、500ml、1000ml、500ml</p>	 <p>No.4</p>	<p>Capacity: 6x500ml/400ml Max Speed: 8000r/min Max RCF: 11620xg Adaptor: 50ml、100ml、250ml</p>



#Let'sgive more ♡
 Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 Optimum Business Center 450 Rue Baden Powell,
 • 34000 Montpellier, France.
 Tel: +33 467 158849 / +33 690760710 / +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A.S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale

<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed
The 3 articles of the week!
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW
How does the climate affect your microscope?
microscope is an instrument that uses a series of lenses to magnify images of small objects. This is a very important tool in many fields, such as biology, medicine, and materials science.
Kalstein Blog
Article of the week:
How does the climate affect your microscope?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW
What is the importance of the laboratory refrigerator?
Refrigeration is a very important part of many laboratory processes. It is used to store samples, reagents, and equipment that require low temperatures to maintain their stability and performance.
Kalstein Blog
Article of the week:
What is the importance of the laboratory refrigerator?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Alliances
23 views • 3 years ago



Kalstein's Microscopes
445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER
2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89
Cel: +33 (0) 6 80 76 07 10
sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

