



A different accompaniment, at your service



 **Kalstein**®

All Right Reserved © KALSTEIN France S. A. S

**LABORATORY SECTOR**  
**ANALYZERS**

[www.kalstein.net/en/](http://www.kalstein.net/en/)



All rights reserved © KALSTEIN France S. A. S.,  
2 Rue Jean Lantier • 75001 Paris •  
+33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
KALSTEIN FRANCE, S. A. S





# Kalstein®

A different accompaniment, at your service

OUR SERVICES

## Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



### Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



### Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



### #Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



### Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



### Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



### Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.





# KALSTEIN® Analyzers

## Fully Automatic Chemistry Analyzer

### Model YR06074



#### Product Description:

- ✓ Automatic, Random Access.
- ✓ Automatic 8 steps washing system, Low carry over.
- ✓ Labor saving, simple operation and real walk-away system.
- ✓ 24 hours non-stop reagent cooling compartment.
- ✓ Pre-diluent function for sample.
- ✓ Low consumption of the reagent.
- ✓ User-friendly software.
- ✓ Multi-language supported.



#### Technical Description:

Model	YR06074
<b>System Function</b>	Throughput: Up to 160 tests per hour
	Contents of the testing: Clinical chemistry, Immunoturbidimetric
	Analytical modes: End-point, Two-point, Kinetic
	Programming: Open system with user defined programs
	Automatic Random Access
	STAT sample priority
<b>Sample Tray</b>	Sample position: 40 positions for samples
	Sample volume: 3- 100µL, 0.5µL adjustable
	Any position for standard and quality control setting
<b>Reagent Tray</b>	Reagent position: 40 positions in refrigerated compartment
	Refrigerated temperature: ≤16C
	Reagent volume
	R1:180~450µL, 1µL adjustable
	R2: 30~250µL, 1µL adjustable
<b>Working Condition</b>	Power supply: a.c. 110V/220V, 50Hz
	Temperature: 0~40°C
	Humidity: ≤80%
	Water consumption: Max 3 W/hour
	Dimensions: 810mm (L)*680mm (W)*1150mm (H)
	Gross weight: 145KG



#### #Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
2 Rue Jean Lantier • 75001 Paris •  
+33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
KALSTEIN FRANCE, S. A. S





#LETSGIVEMORE 

**WITH THE  
ACQUISITION  
OF A KALSTEIN  
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



Fundación  
Jacinto Convit

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



**Humatem**

Solidarité et coopération  
biomédicale

<https://www.humatem.org/>



**MANIAPURE  
FOUNDATION USA**

<https://www.maniapurefoundation.org/>



**Articles Kalstein**  
(Social Networks and Official Websites)



**Kalstein Blog**  
Not just be informed, be fully informed

**The 3 articles of the week!**

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>

**1**  
Read this ARTICLE NOW



**How does the climate affect your microscope?**

microscope is an instrument that allows observing and measuring small objects. This is achieved through a series of lenses that focus light on the object being observed.

**Kalstein Blog**

**Article of the week:**  
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>

**2**  
Read this ARTICLE NOW



**What is the importance of the laboratory refrigerator?**

Refrigerator is a device that is one of the most essential equipment that allows to maintain, in a controlled environment (refrigerated space), various fluids and substances.

**Kalstein Blog**

**Article of the week:**  
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>

**Kalstein Training Videos**  
(Youtube Channel and Official Websites)



**Kalstein France - Jacinto Convit Foundation / Allianc...**  
23 views • 3 years ago



**Kalstein's Microscopes**  
445 views • 1 week ago



PARIS - FRANCE  
**HEADQUARTER**  
2 Rue Jean Lantier  
75001 Paris - France

**Tlf:** +33 (0) 01 78 95 87 89  
**Cel:** +33 (0) 6 80 76 07 10  
[sales@kalstein.eu](mailto:sales@kalstein.eu)

<https://kalstein.eu/>  
<https://www.kalstein.fr/>

