



A different kind of support at your service



**Kalstein**®

All Rights Reserved © KALSTEIN France S. A. S.,

**LABORATORY SECTOR**

**Hba1C  
ANALYZER**



All rights reserved KALSTEIN France S.A.S.  
5 Rue Castiglione 75001 Paris - France  
+33 (0) 1 70 39 26 50 / +33 (0) 6 80 76 07  
<https://kalstein.net>  
KALSTEIN FRANCE, S.A.S

[www.kalstein.net/en/](http://www.kalstein.net/en/)





## Benefits & Assistance

At Kalstein France we care about the full satisfaction of our customers, that's why we offer value-added services of the highest level based on our expertise.

### Quick Response

Our team is always available to answer all your questions or queries, in order to support you in any situation...

### Technical Support

Enjoy personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein manuals and articles, special catalogs and video tutorials.

### #Letsgivemore

Thanks to your purchase, a donation is made to a non-profit foundation that fights cancer and helps vulnerable communities.

### AI Model Recommendation

Model recommendation through AI is used to analyze user behavior, preferences and specific characteristics in order to recommend the most suitable product models for each customer.

### Kalstein+ Platform

All of this would be presented in an intuitive and easy-to-use interface, providing a personalized and satisfying shopping experience for each customer.

### Kalstein worldwide

With more than 25 years growing with our customers, Kalstein's multi-format and modern content is now present in more than 10 countries and growing.

### Delivery Logistics

We take care of all the logistics necessary for the shipment of your goods, whether by sea, land or air.



# KALSTEIN® Hba1c Analyzer

## HbA1c Analyzer

Model YR05129



### Product Description:

- ✓ HbA1c Analyzer is an advanced and compact device designed to measure HbA1c levels in whole blood samples. It provides fast and accurate results in just a few minutes, facilitating diabetes monitoring. This device is easy to use, with a user-friendly interface and a portable design, making it ideal for various clinical settings.
- ✓ Additionally, it features data storage capabilities and connectivity for hospital information systems, enhancing efficiency in patient result management. An optional accessory for the HbA1c Analyzer, is available to further extend its functionality.



### Technical Specifications:

<b>Model</b>	YR05129
<b>Test parameters</b>	HbA1c, eAG
<b>Sample type</b>	whole blood
<b>Sample volume</b>	10µL
<b>Test time</b>	3-4min.
<b>Measurement range</b>	HbA1c 4.0%-15.0%
<b>CV</b>	<5%
<b>Storage</b>	10,000 sample; 12 months QC value, power failure protection
<b>COM</b>	RS-232, HIS/LIS
<b>Power requirement</b>	AC220V±10%, 50Hz, 50VA, with Li battery option
<b>Dimension</b>	205mm*130mm*70mm, 1KG
<b>Working temperature</b>	5-40°C
<b>Consumables</b>	HbA1c cartridge, kits, IC card



#Let'sgive more 

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
5 Rue Castiglione 75001 Paris – France  
+33 (0) 1 70 39 26 50 / +33 (0) 6 80 76 07 10  
<https://kalstein.net>  
KALSTEIN FRANCE, S.A.S







# FOR ALL KALSTEIN CUSTOMERS

## JOIN KALSTEIN+ NEW PLATFORM!

### TECHNOLOGICAL SOLUTIONS!



#### Artificial Intelligence

State-of-the-art tool (integrated into all web pages) necessary for the rapid preparation of tenders or ordinary requests.



#### Technological Device

A cutting-edge device designed to seamlessly connect your lab equipment with our mobile app on Kalstein phones..



#### K+ Platform

The distributor client can build their business and make their quotes automatically with shipping costs to the final destination always with the identity of their company.

### BUSINESS SOLUTIONS!



#### Marketing Service

Creation of content for the Blog and informative guides for SEO positioning



#### Dropshipping

For the reception and preparation of the order on behalf of your company and sent to the end customer connected to Kalstein tracking with product customization.



#### Special Discounts

Up to 41% plus 4% on shipping

### INTRODUCING THE NEW ERA OF KALSTEIN INTELLIGENT EQUIPMENT

We are pleased to announce the development of a revolutionary addition to our product range: a cutting-edge device designed to seamlessly connect your laboratory equipment with our mobile app on Kalstein phones.

### CONTACT US

(+33) 01 70 39 26 50 / (+33) 7 60 83 92 83

Mail: [sales.department@kalstein.eu](mailto:sales.department@kalstein.eu)

[www.kalstein.net](http://www.kalstein.net)

📍 5 rue de Castiglione 75001 Paris – France.



# CONNECT WITH US ON KALSTEIN+

## DO YOU NEED ASSISTANCE?



### HOW TO REQUEST A VIDEO CALL?

At Kalstein+, we are committed to providing you with the best possible service, adapting to your needs in a personalized and efficient manner. That's why we cordially invite you to join us on a video call.



+ (33) 7 77 66 24 64



### HOW CAN YOU DOWNLOAD THE APPLICATION?

Download the Kalstein+ app now and take the Kalstein experience with you wherever you go! Don't miss the opportunity to stay up to date with the latest news and special offers.



Download it on  
the Play Store



Download it on  
the App Store

## CONTACT US

(+33) 01 70 39 26 50 / (+33) 7 60 83 92 83

Mail: [sales.department@kalstein.eu](mailto:sales.department@kalstein.eu)

[www.kalstein.net](http://www.kalstein.net)

📍 5 rue de Castiglione 75001 Paris – France.