



A different accompaniment at your service



Kalstein®

All Right Reserved © KALSTEIN France S. A. S

MEDICAL SECTOR

INFANT

PHOTOTHERAPY INCUBATOR

www.kalstein.net/en/



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Infant Phototherapy Incubator

Infant Incubator

Model YR02188



Product Description:

- ✓ 7-inch LCD touch screen control panel
- ✓ Air temperature and skin temperature servo-controlled by computer
- ✓ Various and self-check alarms
- ✓ Free step mattress tilting adjustment;
- ✓ Removable humidity trough, easy to clean
- ✓ >37°C temperature setting
- ✓ With the transfusion shelf and tray
- ✓ Added with guardrail
- ✓ RS-232 connector
- ✓ With drawer and cabinet
- ✓ Has humidity display and control function
- ✓ Phototherapy is optional
- ✓ Electronic lift is optional

Model	YR02188
Electrical requirement	~220V 50Hz
Heater Power	550VA
Air mode temperature range	25.0°C-37.0°C
Skin mode temperature range	34.0°C-37.0°C
Air mode temperature alarm	<38.0°C
Skin mode temperature alarm	< 40.0°C(Air temperature)
Air mode high temperature range	37.1°C~39.0°C
Skin mode high temperature range	37.1°C~38.0°C
Hight temperature mode alarm	< 40.0°C (Air temperature)
Air mode temperature display range	20.0°C-45.0°C
Skin mode temperature display range	25.0°C~45.0°C
Air mode alarm	±3.0°C
Skin mode alarm	±1.0°C
Skin temperature sensor accuracy	±0.3°C
Temperature control accuracy	≤0.5°C
Mattress temperature uniformity	≤0.8°C
Humidity control range	30-90%RH
Humidity display range	0-99.9%RH



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



Fundación
Jacinto Convit

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale

<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed

The 3 articles of the week!

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

1
Read this ARTICLE NOW



How does the climate affect your microscope?

microscope is an instrument that allows observing and measuring small objects. This is achieved through a series of lenses that focus light on the object to be observed.

Kalstein Blog

Article of the week:
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

2
Read this ARTICLE NOW



What is the importance of the laboratory refrigerator?

Refrigerator is a device that is one of the most common equipment that is used to maintain, in a controlled environment (refrigerated space), various fluids and substances.

Kalstein Blog

Article of the week:
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...
23 views • 3 years ago



Kalstein's Microscopes
445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER
2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89
Cel: +33 (0) 6 80 76 07 10
sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

