



A different accompaniment at your service



Kalstein®

All Right Reserved © KALSTEIN France S. A. S

MEDICAL SECTOR

**INFANT RADIANT
WARMER**

www.kalstein.net/en/



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Infant Radiant Warmer
LCD Touch Screen Infant Radiant Warmer
Model YR06216



Product Description:

- ✓ Skin temperature and manual control by micro-computer
- ✓ Separate setting temperature, real temperature and timer display windows
- ✓ Available adjustment of radiator and mattress tilting
- ✓ Four self-check alarms
- ✓ With one illumination light
- ✓ RS-232 connector
- ✓ With X-ray cassette
- ✓ LED phototherapy is optional
- ✓ Electronic lift is optional

Optional Accesories:

- ✓ Phototherapy Light
- ✓ LED Phototherapy Light
- ✓ Suction Unit
- ✓ X-ray tray
- ✓ Oxygen Bottle
- ✓ Oxygen Flow Meter

Technical Description:

Model	YR06216
Electrical requirement	~ 220V 50Hz or ~ 230V 60Hz
Power consumption	≤700VA
Air mode temperature range	25.0°C-37.0°C
Over temperature alarm	<38.0°C
High mode temperature range	37.1°C ~ 38.0°C
High mode temperature alarm	< 39.0°C
Temperature deviation alarm	±1.0°C
Skin sensor accuracy	≤±0.3°C
Mattress temperature Uniformity	≤2.0°C
Heating rate	≤30min
APGAR timing	1min. 5min. 10min



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 2 Rue Jean Lantier • 75001 Paris •
 +33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale

<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed
The 3 articles of the week!
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW
How does the climate affect your microscope?
microscope is an instrument that allows observing and measuring small objects. This instrument is composed of several parts, which form a system that allows observing and measuring small objects.
Kalstein Blog
Article of the week:
How does the climate affect your microscope?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW
What is the importance of the laboratory refrigerator?
Refrigerator is a device that is one of the most important equipment that allows to maintain, in a controlled environment (refrigerated space), various fluids and substances.
Kalstein Blog
Article of the week:
What is the importance of the laboratory refrigerator?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Alliances

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER

2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89

Cel: +33 (0) 6 80 76 07 10

sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

