




LABORATORY SECTOR

REFRIGERATORS AND FREEZERS



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Refrigerators and Freezers**Laboratory Medical Combined Refrigerator and Freezer For Vaccine****Model YR05321-1****Product Description:**

Medical Combined Refrigerator and Freezer is widely used in hospitals, blood stations and laboratories. The product has double chamber, double display, double control system, two sets of individual compressor and refrigeration system. Temperature in both chambers can be independently adjusted.

Features:

- ✓ **Dual Cooling System:** This unit is equipped with two independent cooling systems, one for the upper compartment and one for the lower compartment. Each system has its own independent control system, allowing for individual temperature settings.
- ✓ **Ultimate Uniformity:** The chambers maintain temperatures within $\pm 3^{\circ}\text{C}$ throughout the entire compartment, ensuring confidence that items are stored at the correct temperature regardless of their placement in the chamber.
- ✓ **Fast Recovery:** Fast temperature recovery even after prolonged door openings. The high efficient cooling system maintains consistent temperatures, ensuring that preservations are stored at the right temperature even when the units is opened frequently.
- ✓ **Optimized Fluctuation:** Maintain superior temperature stability, avoiding rapid and significant changes in temperature. This ensures that samples, medications, vaccines, and Blood therapies are consistently stored in their optimal environment, without compromise.
- ✓ **Natural Refrigerants:** Sustainable U.S. EPA, SNAP, and EU F-Gas compliant natural refrigerants and foaming agents support sustainability initiatives. The cooling system uses R290 or R600a Which are environmentally friendly, having no impact on ozone depletion and a very low Global warming potential(GWP) grade.

**#Letsgivemore**

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® Refrigerators and Freezers**Laboratory Medical Combined Refrigerator and Freezer For Vaccine****Model YR05321-1****Technical Specifications:**

Model	YR05321-1
Shelf And Drawer	Shelves*2 and Drawers*3
Capacity (L/cu.ft)	Upper: 198/6.99
	Lower: 102/3.60
Temperature	Upper: 2°C~8°C
	Lower: -10°C~-25°C
Chamer/Out Shell	HIPS/SPCC sprayed steel
Doot Type	Upper: Solid Door
	Lower: Solid Door
Data Logger	Optional
Dismension(W*D*H)(mm)	650*620*1980

**#Letsgivemore**

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:

 **Fundación
Jacinto Convit**
<https://www.jacintoconvit.org/>

 **ONE
TREE
PLANTED**

<https://onetreepanted.org>

 **Humatem**
Solidarité et coopération
biomédicale
<https://www.humatem.org/>


**MANIAPURE
FOUNDATION USA**
<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed

The 3 articles of the week!

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW

How does the climate affect your microscope?

Kalstein Blog

Article of the week:
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW

What is the importance of the laboratory refrigerator?

Kalstein Blog

Article of the week:
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER
2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89
Cel: +33 (0) 6 80 76 07 10
sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

