




# LABORATORY SECTOR

REFRIGERATORS AND  
FREEZERS



**Kalstein**®

A different accompaniment, at your service





# Kalstein®

A different accompaniment, at your service

OUR SERVICES

## Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



### Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



### Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



### #Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



### Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



### Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



### Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



**KALSTEIN® Refrigerators and Freezers****Large Capacity Biomedical Blood Bag Transport Cool Box with Adjustable Temperature****Model YR06027****Product Description:**

- ✓ Built-in protective basket to store the samples
- ✓ PC lock for easy operation and good appearance ,Two power ports for AC/DC connection.
- ✓ Equipped with protective fuse to avoid the destruction caused by overcurrent. Able to be stacked, saving space.

**Features:**

- ✓ Temperature range: +30°C, +4°C
- ✓ High efficiency air cooling system to ensure good temperature uniformity.
- ✓ Auto switch cooling and heating system to ensure the inside temperature in all climatic conditions.
- ✓ Wide voltage design, able to be used when the voltage is between 187V~242V AC and 12V/24V DC.
- ✓ Accurate temperature control, the temperature accuracy to 0.1°C in digital display.
- ✓ Integrated and strong PE body to resist impact and corrosion and Smooth and flat surface for easy cleaning

<b>Model</b>	YR06027
<b>Capacity (L/cu.ft)</b>	700/24.72
<b>Cooling Method</b>	Forced Air Cooling
<b>Defrost</b>	No Frost / Frost Free
<b>Power(W)</b>	160W
<b>Voltage/frequency (V/Hz)</b>	DC12V DC24V AC220V
<b>Refrigerant</b>	R134a
<b>Material</b>	PE
<b>NT./GT. (kg)</b>	182/209

**#Letsgivemore** ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
2 Rue Jean Lantier • 75001 Paris •  
+33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
**KALSTEIN FRANCE, S. A. S**







 **Kalstein**<sup>®</sup>  
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE  
ACQUISITION  
OF A KALSTEIN  
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación  
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



**Humatem**

Solidarité et coopération  
biomédicale

<https://www.humatem.org/>



**MANIAPURE  
FOUNDATION USA**

<https://www.maniapurefoundation.org/>



**Articles Kalstein**  
(Social Networks and Official Websites)



**Kalstein Blog**  
Not just be informed, be fully informed  
**The 3 articles of the week!**  
Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**1**  
Read this ARTICLE NOW  
**How does the climate affect your microscope?**  
microscope is an instrument that allows observing and measuring small objects. This instrument requires a controlled environment to work properly.

**Kalstein Blog**  
**Article of the week:**  
How does the climate affect your microscope?  
Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**2**  
Read this ARTICLE NOW  
**What is the importance of the laboratory refrigerator?**  
Refrigerator is a piece of equipment that is one of the most important in a laboratory. It is used to maintain, in a controlled environment (refrigerated space), various fluids and substances.

**Kalstein Blog**  
**Article of the week:**  
What is the importance of the laboratory refrigerator?  
Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>

**Kalstein Training Videos**  
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Alliances

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



PARIS - FRANCE  
**HEADQUARTER**

2 Rue Jean Lantier  
75001 Paris - France

**Tlf:** +33 (0) 01 78 95 87 89

**Cel:** +33 (0) 6 80 76 07 10

[sales@kalstein.eu](mailto:sales@kalstein.eu)

<https://kalstein.eu/>  
<https://www.kalstein.fr/>

