




# LABORATORY LINES

VORTEX MIXER



# Kalstein®

A different accompaniment, at your service





# Kalstein®

A different accompaniment, at your service

OUR SERVICES

## Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



### Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



### Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



### #Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



### Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



### Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



### Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



# KALSTEIN® Vortex Mixer

## Mini Vortex Mixer

Model YR06042 // YR06043



### Product Description:

- ✓ Smart and portable / Touch to throw / Brushless motor / Stepless speed regulation.
- ✓ It can mix whatever you want, such as nail polishes, eyelash adhesives, paints, glues, tattoo ink, pigments and other difficult liquids.
- ✓ Applicable to nail salons, beauty salons, artist studios, tattoo stores, schools, businesses and home use.
- ✓ Suitable for centrifuge tube, Eppendorf tube and so on..

### Technical Specifications:

Model	YR06043	YR06042
Control method	Speed regulation by push button	
Motor	Brushless motor	
Trembling orbit	4,5 mm (circle)	
Load capacity	50 ml	
Agitation speed	0 – 3000RPM	3000RPM
Enclosure protection level	IP42	
Power	12W	
Voltage	100 – 240V 50/60HZ	
Permissible ambient temperature	50 ~ 40°C, 80%RH	
Dimension	100*93*72mm	93*93*72mm
Packing size	190*126*113mm	
Weight	0.56/0.37kg	0.56/0.36kg



#Letsgivemore ♥

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
Optimum Business Center 450 Rue Baden Powell,  
• 34000 Montpellier, France.  
Tél: +33 467 158849 / +33 690760710/ +33 663810023  
<https://kalstein.eu>  
KALSTEIN FRANCE, S.A. S







 **Kalstein**<sup>®</sup>  
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE  
ACQUISITION  
OF A KALSTEIN  
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación  
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



**Humatem**

Solidarité et coopération  
biomédicale

<https://www.humatem.org/>



**MANIAPURE  
FOUNDATION USA**

<https://www.maniapurefoundation.org/>



**Articles Kalstein**  
(Social Networks and Official Websites)



**Kalstein Blog**  
Not just be informed, be fully informed  
**The 3 articles of the week!**  
Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**1**  
Read this ARTICLE NOW  
**How does the climate affect your microscope?**  
microscope is an instrument that uses a series of lenses to magnify images of small objects. This is a very important tool in many fields, such as biology, medicine, and materials science.  
**Kalstein Blog**  
**Article of the week:**  
How does the climate affect your microscope?  
Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**2**  
Read this ARTICLE NOW  
**What is the importance of the laboratory refrigerator?**  
Refrigeration is a process of removing heat from a space to maintain it at a controlled temperature. In a laboratory, this is a critical process to maintain the stability of samples and equipment.  
**Kalstein Blog**  
**Article of the week:**  
What is the importance of the laboratory refrigerator?  
Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>

**Kalstein Training Videos**  
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



PARIS - FRANCE  
**HEADQUARTER**  
2 Rue Jean Lantier  
75001 Paris - France

**Tlf:** +33 (0) 01 78 95 87 89  
**Cel:** +33 (0) 6 80 76 07 10  
[sales@kalstein.eu](mailto:sales@kalstein.eu)

<https://kalstein.eu/>  
<https://www.kalstein.fr/>

