



A different accompaniment, at your service



 **Kalstein**®

All Right Reserved © KALSTEIN France S. A. S

LABORATORY SECTOR

G E L D O C U M E N T A T I O N

www.kalstein.net/en/



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Gel Documentation Multifunctional Imager Model YR06210 // YR06211



Product Description:

The Multifunctional Imagers provide high performance with ease of use for visible light (RGB) and far red/near infrared (FR/NIR) fluorescence and chemiluminescence detection and all general gel documentation applications. They adopt electronic lifting platforms to adjust the best capture area for diverse sizes of the different samples. The removable PC, which is an optional accessory, runs quickly and feasibly, and avoids any after sales risks caused by embedded Windows Pad..

Features:

- ✓ Super sensitive Cooled digital Camera
- ✓ Wide range of applications
- ✓ Electronic lifting platform
- ✓ Optional removable PC
- ✓ Intuitive Software for unlimited users



#Letsgivemore ❤️

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® Gel Documentation Multifunctional Imager Model YR06210 // YR06211



Technical Description:

Model		YR06210	YR06211
Application	Chemiluminescence Imaging, ECL Western Blotting	YES	YES
	Fluorescence Imaging, Fluorescent Western Blotting	YES	YES
	Gel Doc Application	YES	YES
	Protein gel ; SDS Pages	YES	YES
Standard Parameters	Sensor	Back illuminated SCMOS Sensor	Front illuminated CCD Sensor
	Resolution	9.0Megapixels	6 Megapixels
	Pixel Size	3.76 × 3.76 µm	4.54 x 4.54 µm
	Temperature	-30 °C regulated	
	Dynamic Range	4.80D	
	Pixel density	16bit(65536 Grey scales)	
	Lens	F0.95 as standard (optional F0.8) Auto focus	
	Chemiluminescence Sample tray	Electronic Lifting Platform for diverse sizes of western blotting membranes	
	DNA/RNA detection	Electronic Lifting Blue LED Transmission for DNA/RNA detection	
	Protein gel ; SDS Pages	Electronic Lifting White LED Transmission for Protein Gels	
Optional Accessories	Software	SHST Capture and Analysis Software for unlimited number of users	
		Epi-UV, RGB and NIR fluorescent detection Removable industrial computer, 13.1 Inch Touch Screen, Intel processor	



#Letsgivemore ❤️

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





#LETSGIVEMORE

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



Fundación
Jacinto Convit

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale

<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed

The 3 articles of the week!

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

1
Read this ARTICLE NOW



How does the climate affect your microscope?

microscope is an instrument that allows observing and measuring small objects. This is achieved through a series of lenses that focus light on the object being observed.

Kalstein Blog

Article of the week:
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

2
Read this ARTICLE NOW



What is the importance of the laboratory refrigerator?

Refrigerator is a device that is one of the most essential equipment that is used to maintain, in a controlled environment (refrigerated space), various fluids and substances.

Kalstein Blog

Article of the week:
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...
23 views • 3 years ago



Kalstein's Microscopes
445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER
2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89
Cel: +33 (0) 6 80 76 07 10
sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

