




# LABORATORY SECTOR

REFRIGERATORS AND  
FREEZERS



**Kalstein**®

A different accompaniment, at your service





# Kalstein®

A different accompaniment, at your service

OUR SERVICES

## Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



### Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



### Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



### #Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



### Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



### Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



### Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



**KALSTEIN® Refrigerators and Freezers**  
**Multipurpose Medical Transporter**  
**Model YR06025 // YR06026**



**Product Description:**

YR06025 – YR06026 is designed to transport and transfer vaccines, blood, enzymes, hormones, stem cells, extracted RNA, gene libraries and important biological and chemical reagents. The temperature in the chamber can be set to 2~8°C or -25°C degrees. Supporting AC 110V & 220V and vehicle DC 12V/24V power supply, This transport can be widely used in hospitals, laboratories and research institutes.

**Features:**

- ✓ Fast Pull Down
- ✓ Multiple Temp. Options
- ✓ AC&DC Power Supply
- ✓ Multi-Layer Stacking
- ✓ Military-Grade Material
- ✓ Precise Temp. Control
- ✓ 4 Type Of Alarm
- ✓ Eco-Friendly



**#Letsgivemore** ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
 2 Rue Jean Lantier • 75001 Paris •  
 +33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
**KALSTEIN FRANCE, S. A. S**



# KALSTEIN® Refrigerators and Freezers

## Multipurpose Medical Transporter

### Model YR06025 // YR06026



#### Technical Description:

Model	YR06025	YR06026
<b>Performance</b>		
Cooling method	Forced Air Cooling	
Defrost	Manual	
Basket	Basket around Evaporator	
Temperature range	-10°C~-25°C/2~8°C	
Ambient temperature	10~32°C	
<b>Refrigeration system</b>		
Refrigerant	R134a	
Temp. controller	Microprocessor	
Sensor	NTC	
Display	Digital display	
External material	PE	
<b>Alarms</b>		
High low temperature	Yes	
Power failure	Yes	
Sensor error	Yes	
Controller failure	Yes	
Power failure backup	8 hours	
<b>Electrical Data</b>		
Voltage frequency( V/Hz)	AC 110V 220V/DC 12V 24V	
Power (W)	80W	100W
Current (A)	DC 6.5A	DC 8A
Power consumption	1.2kWh/24h	1.6kWh/24h
<b>Dimensions</b>		
Capacity (L/cu.ft)	30/1.06	80/2.83
Exterior size(W*D*H)(mm)	615*370*471.5	937*565*580.5
Interior size(W*D*H)(mm)	364*253*352.5	547*378*369.5
<b>Other</b>		
Data logger	Optional	
Insulation	PURF	
Noise (db)	45	46



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
2 Rue Jean Lantier • 75001 Paris •  
+33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
KALSTEIN FRANCE, S. A. S







 **Kalstein**<sup>®</sup>  
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE  
ACQUISITION  
OF A KALSTEIN  
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación  
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



**Humatem**

Solidarité et coopération  
biomédicale

<https://www.humatem.org/>



**MANIAPURE  
FOUNDATION USA**

<https://www.maniapurefoundation.org/>



**Articles Kalstein**  
(Social Networks and Official Websites)



**Kalstein Blog**  
Not just be informed, be fully informed  
**The 3 articles of the week!**  
Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**1**  
Read this ARTICLE NOW  
**How does the climate affect your microscope?**  
microscope is an instrument used for observing objects too small to be seen by the naked eye. This instrument is composed of two lenses, which form an image of the object being observed.  
**Kalstein Blog**  
**Article of the week:**  
How does the climate affect your microscope?  
Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**2**  
Read this ARTICLE NOW  
**What is the importance of the laboratory refrigerator?**  
Refrigerator is a device that is one of the most common equipment that is used to maintain, in a controlled environment (refrigerated space), various fluids and substances.  
**Kalstein Blog**  
**Article of the week:**  
What is the importance of the laboratory refrigerator?  
Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>

**Kalstein Training Videos**  
(Youtube Channel and Official Websites)



**Kalstein France - Jacinto Convit Foundation / Alliances**  
23 views • 3 years ago



**Kalstein's Microscopes**  
445 views • 1 week ago



PARIS - FRANCE  
**HEADQUARTER**  
2 Rue Jean Lantier  
75001 Paris - France

**Tlf:** +33 (0) 01 78 95 87 89  
**Cel:** +33 (0) 6 80 76 07 10  
[sales@kalstein.eu](mailto:sales@kalstein.eu)

<https://kalstein.eu/>  
<https://www.kalstein.fr/>

