




Kalstein®

A different accompaniment, at your service

LABORATORY LINES

VORTEX MIXER





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Vortex Mixer

Multi-tube Vortex Mixer

Model YR06048



Product Description:

- ✓ Digital display available; Stepless speed regulation; Easy touch/run switching.
- ✓ Microprocessor control, simplistic design appearance, LED display of speed and time.
- ✓ Humanized program design, built-in short mixing and timed operation modes. Stable and quiet operation.
- ✓ Large orbital vortex, efficient performance. Various kinds of foam tube frames and pads for option.
- ✓ Maximum processing 50 samples at the same time to make the experiment more convenient and effective.

Technical Specifications :

Model	YR06048	
RPM	500 ~ 2800rpm	
Amplitude	4mm horizontal swing	
Time Setting	1min ~ 99h59min	
Max Capacity	4.5kg100w	
Input Power	100w	
Dimensions	410 * 276 * 432 mm	
Sample Capacity	Test Tube Rack Specs	Applicable Specs
		(outer diameter mm)
	50*φ10mm	φ9~φ10mm
	50*φ12mm	φ11~φ13mm
		(1.5/2ml Centrifuge tube)
	50*φ15mm	φ14~φ17mm
		(5/15ml Centrifuge tube)
	50*φ20mm	φ18~φ21mm
15*φ25mm	φ24~φ27mm	
15*φ29mm	φ28~φ32mm	
		(50ml Centrifuge tube)



#Letsgivemore ♥

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél: +33 467 158849 / +33 690760710/ +33 663810023
https://kalstein.eu
KALSTEIN FRANCE, S.A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale

<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed

The 3 articles of the week!

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

1
Read this ARTICLE NOW



How does the climate affect your microscope?

Kalstein Blog

Article of the week:
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

2
Read this ARTICLE NOW



What is the importance of the laboratory refrigerator?

Kalstein Blog

Article of the week:
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER
2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89
Cel: +33 (0) 6 80 76 07 10
sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

