



A different kind of support at your service



Kalstein®

All Rights Reserved © KALSTEIN France S. A. S.,

MEDICAL SECTOR OPERATION LAMPS



All rights reserved KALSTEIN France S.A.S.
5 Rue Castiglione 75001 Paris - France
+33 (0) 1 70 39 26 50 / +33 (0) 6 80 76 07
<https://kalstein.net>
KALSTEIN FRANCE, S.A.S

www.kalstein.net/en/



Benefits & Assistance

At Kalstein France we care about the full satisfaction of our customers, that's why we offer value-added services of the highest level based on our expertise.

Quick Response

Our team is always available to answer all your questions or queries, in order to support you in any situation...

Technical Support

Enjoy personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein manuals and articles, special catalogs and video tutorials.

#Letsgivemore

Thanks to your purchase, a donation is made to a non-profit foundation that fights cancer and helps vulnerable communities.

AI Model Recommendation

Model recommendation through AI is used to analyze user behavior, preferences and specific characteristics in order to recommend the most suitable product models for each customer.

Kalstein+ Platform

All of this would be presented in an intuitive and easy-to-use interface, providing a personalized and satisfying shopping experience for each customer.

Kalstein worldwide

With more than 25 years growing with our customers, Kalstein's multi-format and modern content is now present in more than 10 countries and growing.

Delivery Logistics

We take care of all the logistics necessary for the shipment of your goods, whether by sea, land or air.



KALSTEIN® Operation Lamps

Operation Lamp

Model YR02125 // YR02140



Product Description

- ✓ Our ergonomic, sterile operating concept allows for the simple and flexible control of all functions by the surgical team.
- ✓ Almost unlimited LED life expectancy delivers tremendous product reliability and investment protection.
- ✓ Environmentally proactive: low power consumption and durable LEDs minimize the impact on our precious resources.

Features

- ✓ Endo mode.
- ✓ Deep illumination.
- ✓ PWM dimming technology, no flickering.
- ✓ LED cold light, head temperature rises less than 1°C.
- ✓ High color rendering index ensures clear views of the tissue and Organ.
- ✓ Removable sterilizer handle.
- ✓ Light field adjusting.



#Letsgivemore

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
5 Rue Castiglione 75001 Paris – France
+33 (0) 1 70 39 26 50 / +33 (0) 6 80 76 07 10
<https://kalstein.net>
KALSTEIN FRANCE, S.A.S



KALSTEIN® Operation Lamps

Operation Lamp

Model YR02125 // YR02140



Technical Description:

Model	YR02125	YR02126	YR02127	YR02128	YR02129	YR02130	YR02131	YR02132
Head Diameter (mm)	700	500				500/500	500	300
No. Of LED Blubs	80pcs	48pcs		80pcs		160pcs	80pcs	25 pcs
Illumination (Lux)	160000	140000		100000		≥100,000/100,000	≥60,000	≥70,000
CRI (Ra)	Ra>95			Ra>85				92
Color Temperature (K)	3500-5000			4800±300			4500±500	4800±200
Illumination Depth (mm)	1200			≥500				300
Spot Diameter (mm)	120-300			160-220				150
LED Life (h)	>50,000 Hours							
Model	YR02133	YR02134	YR02135	YR02136	YR02137	YR02138	YR02139	YR02140
Head Diameter (mm)	300			250		300		100
No. Of LED Blubs	25 pcs	7pc led bulb		12pc led bulb		19pc led bulb		Halogen Bulb
Illumination (Lux)	≥70,000	≥30000		60000		≥80000		≥50,000
CRI (Ra)	92			100≥Ra≥85		≥85	100≥Ra≥85	-
Color Temperature (K)	4800±200	4500±500		4500±500			3500,4000,4500	4000±500
Illumination Depth (mm)	300	≥700						300
Spot Diameter (mm)	150	100-280						150
LED Life (h)	>50,000 Hours	≥40,000 Hours						



#Letsgivemore

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
5 Rue Castiglione 75001 Paris – France
+33 (0) 1 70 39 26 50 / +33 (0) 6 80 76 07 10
<https://kalstein.net>
KALSTEIN FRANCE, S.A.S





FOR ALL KALSTEIN CUSTOMERS

JOIN KALSTEIN+ NEW PLATFORM!

TECHNOLOGICAL SOLUTIONS!



Artificial Intelligence

State-of-the-art tool (integrated into all web pages) necessary for the rapid preparation of tenders or ordinary requests.



Technological Device

A cutting-edge device designed to seamlessly connect your lab equipment with our mobile app on Kalstein phones..



K+ Platform

The distributor client can build their business and make their quotes automatically with shipping costs to the final destination always with the identity of their company.

BUSINESS SOLUTIONS!

✓ Marketing Service

Creation of content for the Blog and informative guides for SEO positioning

✓ Dropshipping

For the reception and preparation of the order on behalf of your company and sent to the end customer connected to Kalstein tracking with product customization.

✓ Special Discounts

Up to 41% plus 4% on shipping

INTRODUCING THE NEW ERA OF KALSTEIN INTELLIGENT EQUIPMENT

We are pleased to announce the development of a revolutionary addition to our product range: a cutting-edge device designed to seamlessly connect your laboratory equipment with our mobile app on Kalstein phones.

CONTACT US

(+33) 01 70 39 26 50 / (+33) 7 60 83 92 83

Mail: sales.department@kalstein.eu

www.kalstein.net

📍 5 rue de Castiglione 75001 Paris – France.

CONNECT WITH US ON KALSTEIN+

DO YOU NEED ASSISTANCE?



HOW TO REQUEST A VIDEO CALL?

At Kalstein+, we are committed to providing you with the best possible service, adapting to your needs in a personalized and efficient manner. That's why we cordially invite you to join us on a video call.



+ (33) 7 77 66 24 64



HOW CAN YOU DOWNLOAD THE APPLICATION?

Download the Kalstein+ app now and take the Kalstein experience with you wherever you go! Don't miss the opportunity to stay up to date with the latest news and special offers.



Download it on
the Play Store



Download it on
the App Store

CONTACT US

(+33) 01 70 39 26 50 / (+33) 7 60 83 92 83

Mail: sales.department@kalstein.eu

www.kalstein.net

📍 5 rue de Castiglione 75001 Paris – France.