




LABORATORY SECTOR

PCR THERMAL CYCLERS



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® PCR Thermal Cycler
PCR Thermal Cycler
Model YR096TG



Product Description:

The PCR thermal cycler YR096TG from Kalstein is easy to operate with a 7-inch color LCD touch screen. By taking advantage of its excellent thermal block, brings you homogeneous and accurate temperatures, rapidly and precisely controlled ramp rates, and fast and reproducible PCR cycles. The programming is also remarkably simple and intuitive. Thanks to all these advantages, is your ideal equipment far PCR.

Features:

Flexible to share experiment protocols

1000 experiment protocols can be saved. Protocols can also share easily between different via USB flash disk. Personal protocols can be saved in your USB flash disk to quickly set up your own experiment.

Power failure protection design

Unique power-off protection function can save all the set configurations after power-off, and allow the experiment continues when power-on.

intuitive software design

intuitive software design and friendly interface make it simple even for first time users. Specialized user training is not required..

7-inch touch screen operation

With 7-inch LCD touch screen, can function quickly through simple, one-touch commands. Compact design with small footprint, it is also easy to move.



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.






All rights reserved © KALSTEIN France S. A. S.,
 2 Rue Jean Lantier • 75001 Paris •
 +33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® PCR Thermal Cycler
PCR Thermal Cycler
Model YR096TG



Specifications:

Model	YR096TG		
Reaction Volume	0-100 uL		
Thermo Block	Aluminium		
Compatible Consumables			
	0.2ml single tube	0.2ml 8 strip tubes	0.2ml 96-well plate (skirted, semi-skirted, unskirted)
Temperature Control Range of The Block	4° (-99°)		
Temperature Control Mode	Tube mode & Bloc Mode		
Heating Technology of the Block	Peltier		
Gradient Block	12 row		
Gradient Temperature Span	1 °C -40°C		
Gradient Temperature Range	30° (99.9°)		
Lid Temperature Range	40-110°		
Temperature Uniformity	± 0.2°		
Block Temperature Accuracy	± 0.1°C		
Heating Rate	3 S°C/s		
Cooling Rate	2.S°C/s		
Interface	USB, Ethernet		
Dimensions (W*D*H)	260mm*400mm*260mm		
Weight	11 kg		
Power Supply	AC 110-240V, 50-60Hz		
Max. Power Consumption	600VA		
Running Noise	<55dB		



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 2 Rue Jean Lantier • 75001 Paris •
 +33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



Fundación
Jacinto Convit

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale

<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed

The 3 articles of the week!

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW

How does the climate affect your microscope?

Kalstein Blog

Article of the week:
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW

What is the importance of the laboratory refrigerator?

Kalstein Blog

Article of the week:
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Alliances
23 views • 3 years ago



Kalstein's Microscopes
445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER
2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89
Cel: +33 (0) 6 80 76 07 10
sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

