



A different accompaniment at your service



 **Kalstein**®

All Right Reserved © KALSTEIN France S. A. S

LABORATORY SECTOR

CENTRIFUGES

www.kalstein.net/en/



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Centrifuges**PRP Centrifuge****Model YR469****Product Description:**

PRP produced from a person's own blood. It is a concentration of one type of cell, known as platelets, which circulate through the blood and are critical for blood clotting.

Product Details:

PRP produced from a person's own blood. It is a concentration of one type of cell, known as platelets, which circulate through the blood and are critical for blood clotting. Platelets and the liquid plasma portion of the blood contain many factors that are essential for the cell recruitment, multiplication and specialization that are required for healing.

**Technical Performance:**

- ✓ It has micro-computer control, electronic lock and DC brushless motor which is low noise and high speed accuracy.
- ✓ Digital LCD screen make visualized operation conveniently.
- ✓ Special program customized is made for any kinds of PRP kit.
- ✓ Up to 9 user programs storage meet more demands.
- ✓ It has self-locking, over-speed detection, over-temperature detection, imbalance detection that is automatic alarm device.
- ✓ It has 9 levels acceleration and 10 levels deceleration to be setting according to the different requirement. Slow down freely is available.
- ✓ Silastic sealed ring comfort to GMP, US FDA certificated.

Technical Specifications:

Model	YR469 (Maintain Beauty)	YR469 (Fat Transfer)
Max speed	5000r/min	5000r/min
Max RCF	4100xg	4390xg
Max capacity	4/8x50ml (PRP Tube)	8x50ml (Syringe)
Speed accuracy	± 30r/min	
Time setting	1min~99min59s	
Noise	<65db(A)	
Power supply	AC110V/AC220V, 50Hz-60Hz	
Consumption	250w	
Dimension(L×W×H)	410mm×490mm×310mm	
Packing size(L×W×H)	500mm×580mm×400mm	
Net weight (No rotor)	25kg	

**#Let'sgive more** ♥

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél: +33 467 158849 / +33 690760710 / +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S



KALSTEIN® Centrifuges

PRP Centrifuge

Model YR469



CENTRIFUGES

 <p>No.1</p>	<p>PRP 4x50ml/8 x50ml Adaptor: 50ml、 10ml 、 2 x 10ml</p>  <p>50ml</p>  <p>10ml</p>  <p>2x10ml</p>	 <p>No.1</p>	<p>Syringe 4x50ml/8x50ml Adaptor: 4x50ml、 8x50ml、 4x20ml、 4x10ml</p>  <p>50ml</p>  <p>2x4x50ml</p>
---	---	---	---

Kits:

- ✓ Blood Collection
- ✓ Transfer Blood to PRP Kit
- ✓ PRP Kit Ready to Spin
- ✓ Red Blood Cell PRP Plasma
- ✓ PRP Extraction



Blood Collection
Transfer Blood to PRP Kit



PRP Kit Ready to Spin



Red Blood Cell
PRP



PRP Extraction



#Let'sgive more ♥

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél: +33 467 158849 / +33 690760710 / +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale


<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed

The 3 articles of the week!

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW

How does the climate affect your microscope?

Kalstein Blog

Article of the week:
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW

What is the importance of the laboratory refrigerator?

Kalstein Blog

Article of the week:
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Alliances

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER

2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89

Cel: +33 (0) 6 80 76 07 10

sales@kalstein.eu

<https://kalstein.eu/>

<https://www.kalstein.fr/>

