



A different accompaniment at your service



 **Kalstein**®

All Right Reserved © KALSTEIN France S. A. S

LABORATORY SECTOR

CENTRIFUGES

www.kalstein.net/en/



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Centrifuges**Portable Centrifuge****Model YR142-6****Product Description:**

- ✓ It is popular choice in any clinic and laboratory.
- ✓ Microcomputer control, brushless motor.
- ✓ It is programmable and digital display.

Technical Performance:

- ✓ It is popular choice in any clinic and laboratory.
- ✓ Microcomputer control, brushless motor.
- ✓ It is programmable and digital display.
- ✓ It has speed/RCF switch easily.
- ✓ There is lid safety device, speed safety device, and error detection.
- ✓ Silastic sealed ring comfort to GMP.US FDA certificated.



Model	YR142-6
Max speed	4000r/min
Max RCF	1980xg
Max capacity	8x20ml
Speed accuracy	± 30r/min
Time setting	1s~99min
Noise	<50db(A)
Power supply	AC110V/AC220V
Consumption	80w
Dimension(L×W×H)	310mm×260mm×220mm
Packing size(L×W×H)	370mm×320mm×280mm
Net weight (No rotor)	7kg

**#Letsgivingore** ♥

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél: +33 467 158849 / +33 680760710/ +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S




KALSTEIN® Centrifuges

Portable Centrifuge

Model YR142-6



<p>No.1</p> 	<p>Capacity: 12x7ml Max Speed:4000r/min Max RCF: 1980xg</p>	<p>No.2</p> 	<p>Capacity: 6x15ml Max Speed: 400r/min Max RCF: 1980xg</p>
<p>No.3</p> 	<p>Capacity: 8x15ml Max Speed: 4000r/min Max RCF: 1980xg</p>	<p>No.3</p> 	<p>Capacity:8x20ml Max Speed: 4000r/min Max RCF: 1980xg</p>
<p>No.3</p> 	<p>Capacity: 12x10ml Max Speed:4000r/min Max RCF: 1980xg</p>		



#Letsgivemore 

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved (© KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tit: +33 467 158849 / +33 690760710/ +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale


<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed
The 3 articles of the week!
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW
How does the climate affect your microscope?
microscope is an instrument used for observing small objects. This is a common biological instrument used for...
Kalstein Blog
Article of the week:
How does the climate affect your microscope?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW
What is the importance of the laboratory refrigerator?
Refrigerator is a device used for one of the most important equipment that is used to maintain, in a controlled environment (refrigerated space), various fluids and substances...
Kalstein Blog
Article of the week:
What is the importance of the laboratory refrigerator?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Alliances

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER

2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89

Cel: +33 (0) 6 80 76 07 10

sales@kalstein.eu

<https://kalstein.eu/>

<https://www.kalstein.fr/>

