



A different kind of support at your service



 **Kalstein**®

All Rights Reserved © KALSTEIN France S. A. S.,

LABORATORY SECTOR
pH METER



All rights reserved KALSTEIN France S.A.S.
5 Rue Castiglione 75001 Paris - France
+33 (0) 1 70 39 26 50 / +33 (0) 6 80 76 07
<https://kalstein.net>
KALSTEIN FRANCE, S.A.S

www.kalstein.net/en/



Benefits & Assistance

At Kalstein France we care about the full satisfaction of our customers, that's why we offer value-added services of the highest level based on our expertise.

Quick Response

Our team is always available to answer all your questions or queries, in order to support you in any situation...

Technical Support

Enjoy personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein manuals and articles, special catalogs and video tutorials.

#Letsgivemore

Thanks to your purchase, a donation is made to a non-profit foundation that fights cancer and helps vulnerable communities.

AI Model Recommendation

Model recommendation through AI is used to analyze user behavior, preferences and specific characteristics in order to recommend the most suitable product models for each customer.

Kalstein+ Platform

All of this would be presented in an intuitive and easy-to-use interface, providing a personalized and satisfying shopping experience for each customer.

Kalstein worldwide

With more than 25 years growing with our customers, Kalstein's multi-format and modern content is now present in more than 10 countries and growing.

Delivery Logistics

We take care of all the logistics necessary for the shipment of your goods, whether by sea, land or air.



KALSTEIN® Ion Meter

Portable pH / Ion Meter

Model YR01817 // YR01817-15



Product Description:

- ✓ The Conductivity Meter YR01817 This advanced multi-functional meter features a large LCD screen with a controllable white backlight. It includes a built-in microprocessor chip with intelligent functions such as auto calibration, ATC, data storage, function settings, auto diagnosis, auto power off, and undervoltage prompt.
- ✓ pH Mode: Supports auto recognition of 25 buffers and up to 3-point calibration.
- ✓ Conductivity Mode: Measures with high accuracy from 0.5 $\mu\text{S}/\text{cm}$ to 200 mS/cm and supports TDS, salinity, and resistivity.
- ✓ DO Mode: Offers ATC, auto salinity compensation, and manual pressure compensation with a quick polarization time of 3-5 minutes.
- ✓ ISE Mode: Measures electric potential, pH, PX, ORP, ion concentration, and temperature, supporting 6 concentration units and up to 5-point calibration.
- ✓ The YR01817-1 versatile measurement instrument features a 6.5-inch LCD screen with a user-friendly interface. It includes a built-in microprocessor chip offering intelligent functions like automatic calibration, Automatic Temperature Compensation (ATC), data storage, clock display, USB output, function settings, wireless printing, and self-diagnosis.
- ✓ The device tests electric potential, pH, ORP, ion concentration, and temperature, and supports unit switching between ppm, g/L, and mol/L. Optional Bluetooth printing and real-time data transfer to mobile phones add convenience. It has a smart electrode status display, durable Omron keys rated for over 100,000 uses, and can store up to 1,000 sets of test data, transferable to a USB stick for easy access with Excel.
- ✓ The YR01817-2 model is an advanced measurement instrument equipped with a 7-inch colored capacitive touch screen, offering high resolution (1024*600) and sensitivity for an enhanced user experience. It features a built-in microprocessor chip with intelligent functions like automatic calibration, Automatic Temperature Compensation (ATC), data storage, clock display, USB output, function settings, wireless printing, and self-diagnosis.
- ✓ This device is capable of testing electric potential, pH, DX, ORP, ion concentration, and temperature, with support for PC connections. Optional wireless Bluetooth printing and real-time data transfer to mobile phones add convenience.
- ✓ A smart electrode status display ensures accuracy, and the instrument can store up to 1,000 sets of test data, transferable to a USB memory stick for easy access with Excel.



#Letsgivemore 

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
5 Rue Castiglione 75001 Paris – France
+33 (0) 1 70 39 26 50 / +33 (0) 6 80 76 07 10
<https://kalstein.net>
KALSTEIN FRANCE, S.A.S



KALSTEIN® Ion Meter

Portable pH / Ion Meter

Model YR01817 // YR01817-15



Model		YR01817	YR01817-1	YR01817-2
Measurement Range	pH/pX	(-2.00~19.99)pH	(-2.00~20.00) pH, (0.00~14.00) PX	
	mV	±1999.9mV	(-1999.9~1999.9) mV	
	Concentration	10-19990), optional units: mol/L, ppm, g/L		
	Temperature	(0~100)°C (Auto/manual)		(-10~110.0) °C
Resolution	PH/pX	0.01pH/pX		
	mV	1mV		
	Concentration	4 significant digits expressed in scientific notation)		
	Temperature	0.1°C		
Accuracy	PH/pX	±0.02pH/pX		
	mV	±0.1%FS		
	Concentration	±0.1%	±0.5%	
	Temperature	0.4°C	±0.1°C	
Other	Data storage	200 sets	1000 sets	
	Power	1.5V x2	12V/1A	
	Output	USB		
	Ambient Temp.	5-35°C		
	Ambient Humidity	<85%		
	IP grade	IP57	IP54	



#Letsgivemore

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
5 Rue Castiglione 75001 Paris – France
+33 (0) 1 70 39 26 50 / +33 (0) 6 80 76 07 10
<https://kalstein.net>
KALSTEIN FRANCE, S.A.S





FOR ALL KALSTEIN CUSTOMERS

JOIN KALSTEIN+ NEW PLATFORM!

TECHNOLOGICAL SOLUTIONS!



Artificial Intelligence

State-of-the-art tool (integrated into all web pages) necessary for the rapid preparation of tenders or ordinary requests.



Technological Device

A cutting-edge device designed to seamlessly connect your lab equipment with our mobile app on Kalstein phones..



K+ Platform

The distributor client can build their business and make their quotes automatically with shipping costs to the final destination always with the identity of their company.

BUSINESS SOLUTIONS!

✓ Marketing Service

Creation of content for the Blog and informative guides for SEO positioning

✓ Dropshipping

For the reception and preparation of the order on behalf of your company and sent to the end customer connected to Kalstein tracking with product customization.

✓ Special Discounts

Up to 41% plus 4% on shipping

INTRODUCING THE NEW ERA OF KALSTEIN INTELLIGENT EQUIPMENT

We are pleased to announce the development of a revolutionary addition to our product range: a cutting-edge device designed to seamlessly connect your laboratory equipment with our mobile app on Kalstein phones.

CONTACT US

(+33) 01 70 39 26 50 / (+33) 7 60 83 92 83

Mail: sales.department@kalstein.eu

www.kalstein.net

📍 5 rue de Castiglione 75001 Paris – France.

CONNECT WITH US ON KALSTEIN+

DO YOU NEED ASSISTANCE?



HOW TO REQUEST A VIDEO CALL?

At Kalstein+, we are committed to providing you with the best possible service, adapting to your needs in a personalized and efficient manner. That's why we cordially invite you to join us on a video call.



+ (33) 7 77 66 24 64



HOW CAN YOU DOWNLOAD THE APPLICATION?

Download the Kalstein+ app now and take the Kalstein experience with you wherever you go! Don't miss the opportunity to stay up to date with the latest news and special offers.



Download it on
the Play Store



Download it on
the App Store

CONTACT US

(+33) 01 70 39 26 50 / (+33) 7 60 83 92 83

Mail: sales.department@kalstein.eu

www.kalstein.net

📍 5 rue de Castiglione 75001 Paris – France.