




# LABORATORY SECTOR

CENTRIFUGE



## Kalstein®

A different accompaniment, at your service





# Kalstein®

A different accompaniment, at your service

OUR SERVICES

## Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



### Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



### Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



### #Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



### Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



### Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



### Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



**KALSTEIN® Centrifuges**  
**Refrigerated High-RPM Centrifuge**  
**Model YR0137**



**Product Description:**

- ✓ Microprocessor control, touch panel for easy visibility and access, LED or LCD display (user's choice), AC frequency conversion motor without carbon dust contamination and CFC-free cooling system (France Tecumseh compressor unit).
- ✓ Automatic rotor identification system to prevent overspeed.
- ✓ Electronic door interlock for safe operation.
- ✓ It is programmable and can store ten programs in memory.
- ✓ Fast access and Dec: the fastest access time is 20 s, and the shortest Dec time is 18 s. Acc and Dec time can be adjusted.
- ✓ The RCF value can be calculated automatically.
- ✓ High speed centrifugation; Low speed centrifugation and high capacity centrifugation all in one.

<b>Model</b>	YR0137
<b>Maximum Speed</b>	20 500 r/min
<b>Maximum RCF</b>	29 200 x g
<b>Maximum Capacity</b>	4 x 750ml Spin Rotor
<b>Accuracy Speed Control</b>	± 50 r/min
<b>Temperature Range</b>	-20 ~ +40°C
<b>Temperature Accuracy</b>	±1.5°C
<b>Timer</b>	9h 59min
<b>Noise</b>	<65db (A)
<b>Engine</b>	AC Frequency Conversion Motor
<b>Power Supply</b>	110V / 60Hz, 220V/50 Hz
<b>Dimension (L x W x H)</b>	650 x 710 x 1100 (mm)
<b>Net Weight</b>	170 Kg



#Letsgivingore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
 2 Rue Jean Lantier • 75001 Paris •  
 +33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
**KALSTEIN FRANCE, S. A. S**







 **Kalstein**<sup>®</sup>  
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE  
ACQUISITION  
OF A KALSTEIN  
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación  
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



**Humatem**

Solidarité et coopération  
biomédicale

<https://www.humatem.org/>



**MANIAPURE  
FOUNDATION USA**

<https://www.maniapurefoundation.org/>



**Articles Kalstein**  
(Social Networks and Official Websites)



**Kalstein Blog**  
Not just be informed, be fully informed  
**The 3 articles of the week!**  
Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**1**  
Read this ARTICLE NOW  
**How does the climate affect your microscope?**  
Microscopes are an essential tool for scientific research, allowing researchers to observe and analyze small objects. This is why it is important to maintain your microscope in a controlled environment, such as a laboratory refrigerator.

**Kalstein Blog**  
**Article of the week:**  
How does the climate affect your microscope?  
Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**2**  
Read this ARTICLE NOW  
**What is the importance of the laboratory refrigerator?**  
Refrigerators in a laboratory are one of the most important pieces of equipment that you need to maintain in a controlled environment (refrigerated space), without fluctuations in temperature.

**Kalstein Blog**  
**Article of the week:**  
What is the importance of the laboratory refrigerator?  
Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>

**Kalstein Training Videos**  
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Alliances

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



PARIS - FRANCE  
**HEADQUARTER**

2 Rue Jean Lantier  
75001 Paris - France

**Tlf:** +33 (0) 01 78 95 87 89

**Cel:** +33 (0) 6 80 76 07 10

[sales@kalstein.eu](mailto:sales@kalstein.eu)

<https://kalstein.eu/>  
<https://www.kalstein.fr/>

