



A different accompaniment at your service



 **Kalstein**®

All Right Reserved © KALSTEIN France S. A. S

LABORATORY SECTOR

CENTRIFUGES

www.kalstein.net/en/



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Centrifuges

Refrigerated Low-speed Centrifuge

Model YR122



Product Description:

It has self-locking device, over-speed safety device, over-temperature safety device, imbalance safety device and automatic alarm device. It has 10 files of acceleration and deceleration level for your choice. Silastic sealed ring comfort to GMP. US FDA certificated.

Technical performance:

- ✓ The max capacity of swing rotor is 4X1000ml.
- ✓ It has micro-computer control, electronic lock and AC variable-conversion motor which delivers high torque.
- ✓ Digital LCD screen make visualized operation conveniently.
- ✓ Imported compressor unit with Non-CFC refrigerant R404a.
- ✓ It has self-locking, over-speed detection, over-temperature detection, imbalance detection that is automatic alarm device.
- ✓ Essential configuration included rotor exchange, speed/RCF, time, temperature, 9 levels acceleration and deceleration is the most economical choice.
- ✓ Silastic sealed ring comfort to GMP, US FDA certificated..



#Let'sgivemore ♥

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél: +33 467 158849 / +33 690760710 / +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A.S



KALSTEIN® Centrifuges
Refrigerated Low-speed Centrifuge
Model YR122



Technical Specifications :

Model	YR122	YR0140
Max speed	5000r/min	5000r/min
Max RCF	4100xg 4640xg	
Max capacity	4x1000ml	
Speed accuracy	± 10r/min	
Time setting	1s~99h 59min59s	
Temperature setting	-20 °C ~ +40 °C	No
Temperature accuracy	±1 °C	No
Compressor unit	Yes	No
Noise	<65db(A)	
Power supply	AC220V, 50Hz-60Hz	
Consumption	2Kw	1.5Kw



#Letsgivingore ♥

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 Optimum Business Center 450 Rue Baden Powell,
 • 34000 Montpellier, France.
 Tit: +33 467 158849 / +33 690760710/ +33 663810023
<https://kalstein.eu>
 KALSTEIN FRANCE, S.A. S




KALSTEIN® Centrifuges
Refrigerated Low-speed Centrifuge
Model YR122



 <p>No.1</p>	<p>Capacity: 6x300ml Max Speed: 5000r/min Max RCF:3850xg Adaptor: 50ml、100ml</p>	<p>YR0140</p>	<p>Capacity: 6x500ml Max Speed: 5000r/min Max RCF: 4640xg Adaptor: 50ml、100ml、250ml</p>
 <p>No.2</p>	<p>Capacity: 4 x500ml Max Speed: 4000r/min Max RCF:3780 xg</p>	<p>Adaptor:4 x20 x1.5ml、4 x20 x5ml/7ml(Blood Vacuum Tube)、4 x 15 x 10ml (Internal Thread)、4 x8 x15ml、4 x9 x20ml、4 x9 x20ml、4 x4 x50ml、4 x2 x100ml、4 x250ml</p>	
 <p>No.2</p>	<p>Capacity: 4 x750ml Max Speed: 4000r/min Max RCF:3780</p>	<p>Adaptor:4 x21 x1.5ml、4 x24 x5ml/7ml(Blood Vacuum Tube)、4 x 37 x 5ml/7ml (G RIA Tube)、4 x 20 x 10ml (Internal Thread)、4 x14 x15ml、4 x12 x20ml、4 x5 x50ml、4 x3 x100ml、4 x250ml、4 x500ml</p>	
 <p>No.3</p>	<p>Capacity: 4 x1000ml Max Speed: 4000r/min Max RCF:4100 xg</p>	<p>Adaptor:4 x21 x1.5ml、4 x24 x5ml/7ml(Blood Vacuum Tube)、4 x 37 x 5ml/7ml (G RIA Tube)、4 x 20 x 10ml (Internal Thread)、4 x14 x15ml、4 x12 x20ml、4 x5 x50ml、4 x3 x100ml、4 x250ml、4 x500ml</p>	
 <p>No.5</p>	<p>Capacity: Square Bucket Max Speed: 4000r/min Max RCF:3040 xg</p>	<p>Adaptor:1.5mlx25x4、5mlx28x4(Blood Vacuum Tube)、10m x 20 x 4l(Internal Thread)、4 x10 x15ml、4 x4 x50ml、4 x250ml</p>	
 <p>No.5</p>	<p>Capacity: 4 x2 x96 Max Speed: 4000r/min Max RCF:3040 xg (Microplate Bucket)</p>		



#Letsgivemore 
 Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 Optimum Business Center 450 Rue Baden Powell,
 • 34000 Montpellier, France.
 Tit: +33 467 158849 / +33 690760710/ +33 663810023
<https://kalstein.eu>
 KALSTEIN FRANCE, S.A.S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale

<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed

The 3 articles of the week!

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW

How does the climate affect your microscope?

Kalstein Blog

Article of the week:
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW

What is the importance of the laboratory refrigerator?

Kalstein Blog

Article of the week:
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER

2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89
Cel: +33 (0) 6 80 76 07 10
sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

