



A different accompaniment at your service



 **Kalstein**®

All Right Reserved © KALSTEIN France S. A. S

# LABORATORY SECTOR

---

# CENTRIFUGES

[www.kalstein.net/en/](http://www.kalstein.net/en/)



All rights reserved © KALSTEIN France S. A. S.,  
2 Rue Jean Lantier • 75001 Paris •  
+33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
KALSTEIN FRANCE, S. A. S





# Kalstein®

A different accompaniment, at your service

OUR SERVICES

## Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



### Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



### Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



### #Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



### Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



### Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



### Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.





## KALSTEIN® Centrifuges

### Refrigerated Low-speed Centrifuge

#### Model YR124-1



#### Product Description:

It has micro-computer control, electronic lock and AC variable conversion motor which delivers high torque.

#### Technical performance:

- ✓ The max capacity of swing rotor is 4X750ml.
- ✓ It has micro-computer control, electronic lock and AC variable conversion motor which delivers high torque.
- ✓ Digital LCD screen make visualized operation conveniently.
- ✓ Imported compressor unit with Non-CFC refrigerant R404a.
- ✓ It has self-locking, over-speed detection, over-temperature detection, imbalance detection that is automatic alarm device.
- ✓ Essential configuration included rotor exchange, speed/RCF, time, temperature, 9 levels acceleration and deceleration is the most economical choice.
- ✓ Silastic sealed ring comfort to GMP, US FDA certificated.



#Let'sgive more ♥

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
Optimum Business Center 450 Rue Baden Powell,  
• 34000 Montpellier, France.  
Tél: +33 467 158849 / +33 680760710 / +33 663810023  
<https://kalstein.eu>  
KALSTEIN FRANCE, S.A.S



**KALSTEIN® Centrifuges**  
**Refrigerated Low-speed Centrifuge**  
**Model YR124-1**



**Technical Specifications :**

<b>Model</b>	YR124	YR124-1
<b>Max speed</b>	5000r/min	5000r/min
<b>Max RCF</b>	3780 xg	
<b>Max capacity</b>	4x750ml	
<b>Speed accuracy</b>	± 10r/min	
<b>Time setting</b>	1s~99h 59min59s	
<b>Temperature setting</b>	-20 °C ~ +40 °C	No
<b>Temperature accuracy</b>	±1 °C	No
<b>Compressor unit</b>	Yes	No
<b>Noise</b>	<65db(A)	
<b>Power supply</b>	AC220V, 50Hz-60Hz	
<b>Consumption</b>	2Kw	1Kw



#Letsgivemore ♥

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
 Optimum Business Center 450 Rue Baden Powell,  
 • 34000 Montpellier, France.  
 Tel: +33 467 158849 / +33 690760710/ +33 663810023  
<https://kalstein.eu>  
 KALSTEIN FRANCE, S.A.S



**KALSTEIN® Centrifuges**  
**Refrigerated Low-speed Centrifuge**  
**Model YR124-1**



<p>No.2</p> 	<p>Capacity: 6x50ml          Max Speed: 5000r/min          Max RCF:3780xg          Adaptor: 15ml、20ml</p>		
<p>No.3</p> 	<p>Capacity: 4 x500ml          Max Speed: 4000r/min          Max RCF:3780 xg</p>  <p>500 ml Bucket</p>  <p>500ml PE Bottle</p>	<p>Adaptor:4 x20 x1.5ml、 4 x20 x5ml/7ml(Blood Vacuum Tube)、 4 x 15 x 10ml (Internal Thread)、 4 x8 x15ml、 4 x9 x20ml、 4 x9 x20ml、 4 x4 x50ml、 4 x2 x100ml、 4 x250ml</p>     <p>4x20x1.5ml    4x20x5ml/7ml Blood Vacuum Tube    Internal Thread    4x8x15ml</p>     <p>4x9x20ml    4x4x50ml    4x2x100ml    4x250ml</p>	



#Letsgivemore ♡  
 Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
 Optimum Business Center 450 Rue Baden Powell,  
 • 34000 Montpellier, France.  
 Tel: +33 467 158849 / +33 690760710/ +33 663810023  
<https://kalstein.eu>  
**KALSTEIN FRANCE, S.A. S**



**KALSTEIN® Centrifuges**  
**Refrigerated Low-speed Centrifuge**  
**Model YR124-1**



 <p data-bbox="414 950 502 1000"><b>No.3</b></p>	<p data-bbox="864 672 1079 819">Capacity: 4 x750ml  Max Speed:  4000r/min  Max RCF:3780</p>	<p data-bbox="1088 687 1738 840">Adaptor: 4 x21 x1.5ml、 4 x24 x5ml/7ml (Blood Vacuum Tube)、 4 x 37 x 5ml/7ml (G RIA Tube)、 4 x 20 x 10ml (Internal Thread)、 4 x14 x15ml、 4 x12 x20ml、 4 x5 x50ml、 4 x3 x100ml、 4 x250ml、 4 x500ml</p>  <p data-bbox="1142 978 1652 1179">4x21x1.5ml   4x24x5ml/7ml   4x37x5ml/7ml   4x20x10ml   4x14x15ml  Blood Vacuum Tube   G RIA Tube   Internal Thread</p> <p data-bbox="1142 1146 1652 1179">4x12x20ml   4x5x50ml   4x3x100ml   4x250ml   4x500ml</p>
 <p data-bbox="451 1703 538 1753"><b>No.5</b></p>	<p data-bbox="864 1452 1079 1648">Capacity: Square Bucket  Max Speed:  4000r/min  Max RCF:3040 xg</p>	<p data-bbox="1088 1496 1738 1605">Adaptor: 4 x25 x1.5ml、 4 x28 x5ml (Blood Vacuum Tube)、 4 x 20 x 10ml (Internal Thread)、 4 x 10x 15ml、 4 x4 x50ml、 4 x250ml</p>  <p data-bbox="1252 1714 1597 1747">1.5mlx25x4   5mlx28x4   10mlx20x4  Blood Vacuum Tube   Internal Thread</p> <p data-bbox="1252 1834 1597 1867">4x10x15ml   4x4x50ml   4x250ml</p>



**#Letsgivemore** ♥  
Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
Optimum Business Center 450 Rue Baden Powell,  
• 34000 Montpellier, France.  
Tit: +33 467 158849 / +33 690760710/ +33 663810023  
https://kalstein.eu  
**KALSTEIN FRANCE, S.A.S**



**KALSTEIN® Centrifuges**  
**Refrigerated Low-speed Centrifuge**  
**Model YR124-1**



 <p>No.2</p>	<p>Capacity: 6x50ml          Max Speed: 5000r/min          Max RCF:3780xg          Adaptor: 15ml、20ml</p>	 <p>2x2x96          (Microplate Bucket)</p>	
 <p>No.3</p>	<p>Capacity: 4 x500ml          Max Speed: 4000r/min          Max RCF:3780 xg</p>	<p>Adaptor:4 x20 x1.5ml、 4 x20 x5ml/7ml(Blood Vacuum Tube)、 4 x 15 x 10ml (Internal Thread)、 4 x8 x15ml、 4 x9 x20ml、 4 x9 x20ml、 4 x4 x50ml、 4 x2 x100ml、 4 x250ml</p>	
 <p>No.3</p>	<p>Capacity: 4 x750ml          Max Speed: 4000r/min          Max RCF:3780</p>	<p>Adaptor:4 x21 x1.5ml、 4 x24 x5ml/7ml(Blood Vacuum Tube)、 4 x 37 x 5ml/7ml (G RIA Tube)、 4 x 20 x 10ml (Internal Thread)、 4 x14 x15ml、 4 x12 x20ml、 4 x5 x50ml、 4 x3 x100ml、 4 x250ml、 4 x500ml</p>	
 <p>No.5</p>	<p>Capacity: Square Bucket          Max Speed: 4000r/min          Max RCF:3040 xg</p>	<p>Adaptor:4 x25 x1.5ml、 4 x28 x5ml (Blood Vacuum Tube)、 4 x 20 x 10ml (Internal Thread)、 4 x 10x 15ml、 4 x4 x50ml、 4 x250ml</p>	
 <p>No.4</p>	<p>Capacity: 2x2 x96          Max Speed: 4000r/min          Max RCF:2200 xg          (Microplate Bucket)</p>		
 <p>No.5</p>	<p>Capacity: 4x2 x96          Max Speed: 4000r/min          Max RCF:3040 xg          (Microplate Bucket)</p>		



**#Letsgivemore** ♥  
 Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
 Optimum Business Center 450 Rue Baden Powell,  
 • 34000 Montpellier, France.  
 Tel: +33 467 158849 / +33 690760710 / +33 663810023  
<https://kalstein.eu>  
**KALSTEIN FRANCE, S.A. S**







 **Kalstein**<sup>®</sup>  
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE  
ACQUISITION  
OF A KALSTEIN  
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación  
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



**Humatem**

Solidarité et coopération  
biomédicale

<https://www.humatem.org/>



**MANIAPURE  
FOUNDATION USA**

<https://www.maniapurefoundation.org/>



**Articles Kalstein**  
(Social Networks and Official Websites)



**Kalstein Blog**  
Not just be informed, be fully informed

**The 3 articles of the week!**

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**1**  
Read this ARTICLE NOW

**How does the climate affect your microscope?**

**Kalstein Blog**

**Article of the week:**  
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**2**  
Read this ARTICLE NOW

**What is the importance of the laboratory refrigerator?**

**Kalstein Blog**

**Article of the week:**  
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>

**Kalstein Training Videos**  
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



PARIS - FRANCE  
**HEADQUARTER**  
2 Rue Jean Lantier  
75001 Paris - France

**Tlf:** +33 (0) 01 78 95 87 89  
**Cel:** +33 (0) 6 80 76 07 10  
[sales@kalstein.eu](mailto:sales@kalstein.eu)

<https://kalstein.eu/>  
<https://www.kalstein.fr/>

