



A different accompaniment, at your service



 **Kalstein**®

All Right Reserved © KALSTEIN France S. A. S

LABORATORY SECTOR

ANALYZERS

www.kalstein.net/en/



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Analyzer

Semi-auto Chemistry Analyzer

Model YR05114



Product Description:

- ✓ Easy Windows operation system with mouse, large color LCD display
- ✓ Reagent open system supports flow cell and cuvette mode. close system on requesting
- ✓ Analytical modes include kinetic, fixed time, biochromatic end point with or without reagent blank or sample blank, linear or non-linear calibration
- ✓ Real-time curve monitoring enhances the reliability of result
- ✓ Large memory, up to 200 test protocols can be programmed
- ✓ Powerful QC function: Westgard Multi-rule, Levey Jennings plot, 2 controls per testing
- ✓ Auto power on self-check
- ✓ Additional lamp seminal
- ✓ Multiform result output including patient comprehensive report
- ✓ Software optional for date management in PC
- ✓ Easy maintenance, no service engineers required to change lamp, tube etc.



Technical Description:

Model	YR05114
Photometric Range	-0.500-3.500Abs
Resolution	0.001Abs(displayed), 0.0001Abs(calculated)
Light Source	Halogen lamp, 6V/10W
Wavelengths	340,405,500,546,578,620,670nm, 1 more filter optional
Wavelength Accuracy	±2nm
Band Width	≤ 10nm
Aspirate Volume	100-3000ul
Cross Contamination Rate	<1%
Memory	200 test protocols, up to 3000 test results
Interface	RS-232,PSTN (optional)
Cuvette Mode	12.5mmx 12.5mm cuvette optional
Temperature Control	25 °C,30 °C;37 °C 10.5°C
Display	7" color LCD
Input	External mouse and optional external keyboard
Output	Built-in printer or external printer (optional)
Net Weight	11,5KG
Dimensions L × W × H (mm)	460×330×190
Power Supply	AC 110V/220V±10%, 50Hz/60Hz



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:

 **Fundación
Jacinto Convit**
<https://www.jacintoconvit.org/>



<https://onetreepanted.org>

 **Humatem**
Solidarité et coopération
biomédicale
<https://www.humatem.org/>


**MANIAPURE
FOUNDATION USA**
<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed

The 3 articles of the week!

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

1
Read this ARTICLE NOW



How does the climate affect your microscope?

Kalstein Blog

Article of the week:
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

2
Read this ARTICLE NOW



What is the importance of the laboratory refrigerator?

Kalstein Blog

Article of the week:
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...
23 views • 3 years ago



Kalstein's Microscopes
445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER
2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89
Cel: +33 (0) 6 80 76 07 10
sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

