



A different accompaniment, at your service



 **Kalstein**®

All Right Reserved © KALSTEIN France S. A. S

**LABORATORY SECTOR**

**ANALYZERS**

[www.kalstein.net/en/](http://www.kalstein.net/en/)



All rights reserved © KALSTEIN France S. A. S.,  
2 Rue Jean Lantier • 75001 Paris •  
+33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
KALSTEIN FRANCE, S. A. S





# Kalstein®

A different accompaniment, at your service

OUR SERVICES

## Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



### Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



### Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



### #Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



### Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



### Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



### Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.





## KALSTEIN® Analyzers

### Semi-auto Chemistry Analyzer

#### Model YR05115



#### Product Description:

- ✓ Easy keypad operation with LCD display
- ✓ Wavelengths range is 330-800nm, 5 standard filters, 3 more filters optional
- ✓ Reagent open system supports flow cell and cuvette mode, |
- ✓ close system on request
- ✓ Analytical modes include kinetic, fixed time, bichromatic end point with or
- ✓ without reagent blank or sample blank, linear or non-linear calibration curves
- ✓ Up to 60 programs and 2200 test results can be stored
- ✓ Excellent QC function, 2 controls per test
- ✓ Internal thermal-sensitive printer
- ✓ Compact design enables USER-DIY maintenance
- ✓ Multi-language software available on request



<b>Model</b>	YR05115
<b>Absorbance Range</b>	-0,500-3,500 Abs
<b>Resolution</b>	0,001Abs (displayed), 0,0001Abs (calculated)
<b>Light Source</b>	Halogen lamp
<b>Wavelengths</b>	340,405,500,546,620nm,3 more filters optional
<b>Wavelength Accuracy</b>	± 2nm
<b>Band Width</b>	≤ 10nm
<b>Cross Contamination Rate</b>	€ 1%
<b>interface</b>	RS-232 serial cable
<b>Flow cell</b>	Metal – quartz flow cell
<b>Temperature Control</b>	25 C, 30°C,37 C: £0.5 C and ambient temperature
<b>Display</b>	LCD display
<b>Output</b>	internal printer
<b>Power Supply</b>	AC 110V/220V ± 10%, 50Hz/60H2
<b>Net Weight</b>	11,5KG
<b>Dimensions L x W x H (mm)</b>	360 × 318 x 160



#### #Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
2 Rue Jean Lantier • 75001 Paris •  
+33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
KALSTEIN FRANCE, S. A. S





 **Kalstein**<sup>®</sup>  
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE  
ACQUISITION  
OF A KALSTEIN  
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



Fundación  
Jacinto Convit

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



**Humatem**

Solidarité et coopération  
biomédicale

<https://www.humatem.org/>



**MANIAPURE  
FOUNDATION USA**

<https://www.maniapurefoundation.org/>



**Articles Kalstein**  
(Social Networks and Official Websites)



**Kalstein Blog**  
Not just be informed, be fully informed

**The 3 articles of the week!**

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>

**1**  
Read this ARTICLE NOW



**How does the climate affect your microscope?**

microscope is an instrument that allows observing and measuring small objects. This is achieved through a series of lenses that focus light on the object being observed.

**Kalstein Blog**

**Article of the week:**  
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>

**2**  
Read this ARTICLE NOW



**What is the importance of the laboratory refrigerator?**

Refrigerator is a device that is one of the most essential equipment that allows to maintain, in a controlled environment (refrigerated space), various fluids and substances.

**Kalstein Blog**

**Article of the week:**  
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>

**Kalstein Training Videos**  
(Youtube Channel and Official Websites)



**Kalstein France - Jacinto Convit Foundation / Allianc...**  
23 views • 3 years ago



**Kalstein's Microscopes**  
445 views • 1 week ago



PARIS - FRANCE  
**HEADQUARTER**  
2 Rue Jean Lantier  
75001 Paris - France

**Tlf:** +33 (0) 01 78 95 87 89  
**Cel:** +33 (0) 6 80 76 07 10  
[sales@kalstein.eu](mailto:sales@kalstein.eu)

<https://kalstein.eu/>  
<https://www.kalstein.fr/>

