


# Tabletop Class N Autoclave

YR03391/YR03392/YR03393/YR03394



### Default Program

Program	Temperature	Pressure	Holding time (min)
Solid 134°C	134°C	210 Kpa	4
Solid 121°C	121°C	110 Kpa	20
Wrapped 134°C	134°C	210 Kpa	4
Wrapped 121°C	121°C	110 Kpa	20
Textile 134°C	134°C	210 Kpa	8
Textile 121°C	121°C	110 Kpa	30
Prion 134°C	134°C	210 Kpa	18
Liquid 121°C	121°C	110 Kpa	30
Drying(optional)	—	—	—



Model	Chamber Size	Capacity	Type	Weight	Overall Size (WxHxD)	Packing Size (WxHxD)	Number of Trays
YR03391	Ø170x320mm	8L	Class N	35kg	420x370x525mm	503x438x640mm	2pcs, Max 3 pcs(optional)
YR03392	Ø200x360mm	12L	Class N	38kg	420x370x595mm	503x438x640mm	2pcs, Max 3 pcs(optional)
YR03393	Ø247x350mm	18L	Class N	45kg	490x455x600mm	570x520x690mm	3pcs, Max 5 pcs(optional)
YR03394	Ø247x450mm	23L	Class N	49kg	490x455x690mm	575x520x780mm	3pcs, Max 5 pcs (optional)

Note: Internal printer plus 90USD