



A different accompaniment at your service



 **Kalstein**®

All Right Reserved © KALSTEIN France S. A. S

LABORATORY SECTOR

CENTRIFUGES

www.kalstein.net/en/



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Centrifuges**Tabletop High Capacity Refrigerated Centrifuge****Model YR120-1****Product Description:**

It has micro-computer control, electronic lock and AC variable-conversion motor which delivers high torque.

Technical performance:

- ✓ It has micro-computer control, electronic lock and AC variable-conversion motor which delivers high torque.
- ✓ Digital display make visualized operation conveniently.
- ✓ Imported compressor unit with Non-CFC refrigerant R134a.
- ✓ It has self-locking, over-speed detection, overtemperature detection, imbalance detection that is automatic alarm device.
- ✓ It has 9 levels acceleration and deceleration to be setting according to the different requirement.
- ✓ Silastic sealed ring comfort to GMP, US FDA certificated.



#Letsgivingore ♥

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél: +33 467 158849 / +33 690760710 / +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A.S



KALSTEIN® Centrifuges**Tabletop High Capacity Refrigerated Centrifuge****Model YR120-1****Technical Specifications :**

Model	YR120-1	YR124-2
Max speed	8000r/min	5000r/min
Max RCF	9100xg	4390xg
Max capacity	4x250ml	72x7ml
Speed accuracy	± 30r/min	
Time setting	1min~99min	
Temperature setting	-5 °C ~ +35 °C	
Temperature accuracy	±1 °C	
Compressor unit	French Tecumseh Compressor Unit R134a	
Noise	<65db(A)	
Power supply	AC110V/AC220V, 50Hz-60Hz	
Consumption	1.2kw	
Dimension(L×W×H)	460mm×710mm×420mm	
Packing size(L×W×H)	580mm×730mm×520mm	
Net weight (No rotor)	80kg	

**#Letsgivingore** ♥

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél: +33 467 158849 / +33 680760710 / +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S



KALSTEIN® Centrifuges

Tabletop High Capacity Refrigerated Centrifuge

Model YR120-1



Rotor for YR120-1

 <p>No.1</p>	<p>Capacity: 4 x250ml Max Speed: 4300r/min Max RCF:3582xg</p>	<p>Adaptor:4 x9 x1.5ml、 4 x7 x5ml、 4 x9 x5ml /7ml (Blood Vacuum Tube)、 4 x 7 x 10ml (Internal Thread)、 4 x5 x15ml、 4 x4 x20ml、 4 x50ml、 4 x100ml</p>  <p>250 ml Bucket 250ml PE Bottle 4x9x1.5ml 4x7x5ml 4x9x5/7ml Blood Vacuum Tube 4x7x10ml Internal Thread 4x5x15ml 4x4x20ml 4x50ml 4x100ml</p>	
 <p>No.3</p>	<p>Capacity: 8x100ml Max Speed: 8000r/min Max RCF:9100xg Adaptor: 10ml、 20ml</p>		



#Letsgivemore 

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 Optimum Business Center 450 Rue Baden Powell,
 • 34000 Montpellier, France.
 Tel: +33 467 158849 / +33 690760710/ +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A.S



KALSTEIN® Centrifuges

Tabletop High Capacity Refrigerated Centrifuge

Model YR120-1



Rotor for YR120-1 and YR124-2

 <p>No.1</p>	<p>Capacity: 4 x100ml Max Speed: 5000r/min Max RCF:4390 xg</p>	<p>Adaptor: 4x4x1.5ml、 4 x4 x5ml/7ml、 4 x2x 5ml、 4 x3 x10ml、 4 x 15ml、 4 x20ml、 4 x50ml</p>	
 <p>No.2</p>	<p>Max Speed: 4000r/min Max RCF:2810 xg</p>	<p>Adaptor: 4x2x1.5ml、 4x4x1.5ml、 4x6x1.5ml、 8x4x1.5ml、 4x2x5ml/7ml、 4 x4x5ml/7ml、 4 x6x5ml/7ml、 8x4x5ml/7ml、 4 x 2 x 10ml、 4 x 4 x 10ml、 4 x 6 x 10ml、 8x4x10ml、 4 x 2 x15ml、 4 x2 x20ml、 4 x4 x15ml、 4 x6 x15ml、 8x4x15ml、 4 x2 x50ml</p>	
 <p>No.2</p>	<p>Capacity: 4 x18x7ml Max Speed: 4000r/min Max RCF:2810xg 75mm length tube only</p>	 <p>No.2</p>	<p>Capacity: 4 x12x7ml Max Speed: 4000r/min Max RCF:2810xg 75mm &100mm length tube</p>
 <p>No.4</p>	<p>Capacity: 2x2 x96 Max Speed: 4000r/min Max RCF:2200 xg (Microplate Bucket)</p>		



#Letsgivemore 

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél: +33 467 158849 / +33 690760710/ +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale

<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed

The 3 articles of the week!

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW

How does the climate affect your microscope?

Kalstein Blog

Article of the week:
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW

What is the importance of the laboratory refrigerator?

Kalstein Blog

Article of the week:
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Alliances

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER

2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89

Cel: +33 (0) 6 80 76 07 10

sales@kalstein.eu

<https://kalstein.eu/>

<https://www.kalstein.fr/>

