



A different accompaniment at your service



 **Kalstein**®

All Right Reserved © KALSTEIN France S. A. S

LABORATORY SECTOR

CENTRIFUGES

www.kalstein.net/en/



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Centrifuges**Tabletop High Capacity Refrigerated Centrifuge****Model YR124-2****Product Description:**

Microcomputer control, DC brushless motor, stable operation, low noise and high speed accuracy. Touching panel, programmable operation, The operation data can be set according to demands and stored automatically..

Large capacity refrigerated centrifuge technical performance:

- ✓ Microcomputer control, DC brushless motor, stable operation, low noise and high speed accuracy.
- ✓ Touching panel, programmable operation, The operation data can be set according to demands and stored automatically.
- ✓ LED display, user-friendly interface, simple and convenient operations.
- ✓ The real-time conversation of Speed and RCF is convenient for operation.
- ✓ It has self-locking device, over-speed safety device, over-temperature safety device, imbalance safety device and automatic alarm device.
- ✓ It has 9 files of acceleration and deceleration level for your choice.
- ✓ Silastic sealed ring comfort to GMP. US FDA certificated

**#Letsgivingore** ♥

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tel: +33 467 158849 / +33 690760710 / +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S



KALSTEIN® Centrifuges**Tabletop High Capacity Refrigerated Centrifuge****Model YR124-2****Technical Specifications :**

Model	YR124-2
Max speed	5000r/min
Max RCF	4390×g
Max capacity	15ml×32
Speed accuracy	± 30r/min
Time setting	1min~99min
Temperature Setting	-5 ~ +35°C
Temperature Accuracy	±1°C
Compressor Unit	Imported refrigeration Units. Non-CFC R134a
Noise	<65db(A)
Power supply	AC110V/60Hz AC220V/50Hz
Dimension(L×W×H)	460mm×710mm×420mm
Packing size(L×W×H)	580mm×730mm×520mm
Net weight	80kg

**#Letsgivemore** ♥

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél: +33 467 158849 / +33 690760710 / +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S



KALSTEIN® Centrifuges

Tabletop High Capacity Refrigerated Centrifuge

Model YR124-2



<p>No.1</p>  <p>Rotor Capacity: 250mlx4 Max Speed: 4300r/min Max RCF: 3582xg</p>	 250ml Bucket  250ml Bottle	 1.5mlx9x4  5mlx7x4	 5/7mlx9x4 vacuum tube  15mlx5x4  20mlx4x4	 50mlx4  100mlx4	<p>No.3</p>  <p>Rotor Capacity: 100mlx8 Max Speed: 8000r/min Max RCF: 9100xg Adapters: 10ml, 20ml</p>
<p>No.1</p>  <p>Rotor Capacity: 100mlx4 Max Speed: 5000r/min Max RCF: 4390xg</p>	<p>No.2</p>  <p>Max Speed: 4000r/min Max RCF: 2810xg</p>	 50mlx2x4 Bucket  15mlx4x4 Bucket	<p>Adapters: 1.5mlx2 5ml/7ml vacuum tube 5ml, 10ml 15ml, 20ml</p> <p>Adapters: 1.5ml, 5ml 5ml/7ml vacuum tube 10ml</p>	 15mlx6x4 Bucket  15mlx8x4 Bucket	<p>Adapters: 1.5ml, 5ml 5ml/7ml vacuum tube 10ml</p> <p>Adapters: 1.5ml, 5ml 5ml/7ml vacuum tube 10ml</p>
<p>No.2</p>  <p>Rotor Capacity: 7mlx72 Max Speed: 4000r/min Max RCF: 2810xg Apply to tubes with 81mm length</p>	<p>No.2</p>  <p>Rotor Capacity: 7mlx48 Max Speed: 4000r/min Max RCF: 2810xg</p>	<p>No.4</p>  <p>Rotor Capacity: 8x96 well Max Speed: 4000r/min Max RCF: 2200xg</p>			



#Letsgivemore 

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél: +33 467 158849 / +33 680760710 / +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale


<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed

The 3 articles of the week!

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW

How does the climate affect your microscope?

Kalstein Blog

Article of the week:
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW

What is the importance of the laboratory refrigerator?

Kalstein Blog

Article of the week:
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Alliances

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER

2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89

Cel: +33 (0) 6 80 76 07 10

sales@kalstein.eu

<https://kalstein.eu/>

<https://www.kalstein.fr/>

