



A different accompaniment at your service



 **Kalstein**®

All Right Reserved © KALSTEIN France S. A. S

LABORATORY SECTOR

CENTRIFUGES

www.kalstein.net/en/



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Centrifuges
Tabletop High Speed Centrifuge
Model YR129



Product Description:

The main working principle of high-speed microcentrifuge: high-speed microcentrifuge is the powerful centrifugal force generated by high-speed rotation of the high-speed microcentrifuge frozen microcentrifuge rotor to accelerate the sedimentation speed of particles in the liquid, and separate the sedimentation coefficient and buoyant density of the sample.

Technical Performance:

- ✓ This device is small, light and easy to carry that is perfect for Eppendorf tubes.
- ✓ It has micro-computer control, electronic lock and DC brushless motor which is low noise and high speed accuracy.
- ✓ Digital display make visualized operation conveniently.
- ✓ It has self-locking, over-speed detection, over-temperature detection, imbalance detection that is automatic alarm device.
- ✓ Silastic sealed ring comfort to GMP, US FDA certificated..



#Letsgivemore ♥

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 Optimum Business Center 450 Rue Baden Powell,
 • 34000 Montpellier, France.
 Tel: +33 467 158849 / +33 690760710 / +33 663810023
<https://kalstein.eu>
 KALSTEIN FRANCE, S.A. S



KALSTEIN® Centrifuges
Tabletop High Speed Centrifuge
Model YR129



Technical Specifications :

Model	YR129
Max speed	16000r/min
Max RCF	17800xg
Max capacity	10x5ml
Speed accuracy	± 30r/min
Time setting	1s~99min59s
Noise	<65db(A)
Power supply	AC110V/AC220V, 50Hz-60Hz
Consumption	250w
Dimension(L×W×H)	230mm×310mm×240mm
Packing size(L×W×H)	290mm×400mm×330mm
Net weight (No rotor)	9kg



#Letsgivingore ♥

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 Optimum Business Center 450 Rue Baden Powell,
 • 34000 Montpellier, France.
 Tel: +33 467 158849 / +33 690760710/ +33 663810023
<https://kalstein.eu>
 KALSTEIN FRANCE, S.A. S.




KALSTEIN® Centrifuges
Tabletop High Speed Centrifuge
Model YR129



 <p>No.1</p>	<p>Capacity: 12x1.5ml/2.2ml Max Speed:16000r/min Max RCF: 17800xg Adaptor: 0.2ml、 0.5ml</p>	 <p>No.2</p>	<p>Capacity: 10x5ml Max Speed: 13000r/min Max RCF: 11400xg Adaptor: 0.5ml、 1.5ml</p>
 <p>No.3</p>	<p>Capacity: 18x1.5ml/2.2ml Max Speed: 13200r/min Max RCF: 16110xg Adaptor: 0.2ml、 0.5ml</p>	 <p>No.4</p>	<p>Capacity: 24x1.5ml/2.2ml Max Speed: 10000r/min Max RCF: 10500xg Adaptor: 0.2ml、 0.5ml</p>



#Let'sgive more 
 Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 Optimum Business Center 450 Rue Baden Powell,
 • 34000 Montpellier, France.
 Tel: +33 467 158849 / +33 690760710/ +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale

<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed

The 3 articles of the week!

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW

How does the climate affect your microscope?

Kalstein Blog

Article of the week:
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW

What is the importance of the laboratory refrigerator?

Kalstein Blog

Article of the week:
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...
23 views • 3 years ago



Kalstein's Microscopes
445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER

2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89

Cel: +33 (0) 6 80 76 07 10

sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

