



A different accompaniment at your service



 **Kalstein**®

All Right Reserved © KALSTEIN France S. A. S

LABORATORY SECTOR

CENTRIFUGES

www.kalstein.net/en/



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Centrifuges
Tabletop High Speed Centrifuge
Model YR137-3



Product Description:

We are a professional enterprise engaged in the manufacture and sale of YR137-3 Tabletop High Capacity Centrifuge. The unique design and stable quality win praises from customers in the field of colleges and universities, research units, pharmaceuticals, medical treatment and petrochemical engineering.

We are a professional enterprise engaged in the manufacture and sale of YR137-3 Tabletop High Capacity Centrifuge. The unique design and stable quality win praises from customers in the field of colleges and universities, research units, pharmaceuticals, medical treatment and petrochemical engineering.

Technical Performance:

- ✓ Microcomputer control, DC brushless motor, stable operation, low noise and high speed accuracy.
- ✓ Touch panel, programmable operation, operation data can be set according to demand and stored automatically.
- ✓ LED or LCD display. user-friendly interface, simple and convenient operations.
- ✓ The real-time conversation Of Speed and RCF is convenient for operation.
- ✓ It has self-locking device, lid safety device, over-speed safety device, and automatic alarm device, etc.
- ✓ It has 10 files of acceleration and deceleration lever for you choice.
- ✓ Silastic sealed ring comfort to GMP.US FDA certificated.



#Let'sgive more ♥

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 Optimum Business Center 450 Rue Baden Powell,
 • 34000 Montpellier, France.
 Tel: +33 467 158849 / +33 690760710 / +33 663810023
<https://kalstein.eu>
 KALSTEIN FRANCE, S.A. S



KALSTEIN® Centrifuges
Tabletop High Speed Centrifuge
Model YR137-3



Technical Specifications :

Model	YR137-3
Max speed	5500r/min
Max RCF	5310×g
Max capacity	500ml×4
Speed accuracy	± 30r/min
Time setting	1min~99min/1min~99h59min
Noise	<65dB(A)
Power supply	110V/60Hz 220V/50Hz
Dimension(LxWxH)	460mm×540mm×340mm
Packing size(LxWxH)	550mm×630mm×430mm
Net weight	31kg



#Letsgivemore ♥

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.
















All rights reserved © KALSTEIN France S. A. S.,
 Optimum Business Center 450 Rue Baden Powell,
 • 34000 Montpellier, France.
 Tel: +33 467 158849 / +33 690760710/ +33 663810023
<https://kalstein.eu>
 KALSTEIN FRANCE, S.A. S



KALSTEIN® Centrifuges Tabletop High Speed Centrifuge Model YR137-3



<p>No.4</p>  <p>Rotor capacity: 100mlx4 Max Speed: 5000r/min Max RCF: 4390xg</p>	 <p>100mlx4 Bucket</p>	<p>Adapters: 1.5mlx4x4</p> <hr/> <p>Adapters: 5mlx2x4</p> <hr/> <p>Adapters: 15mlx4</p> <hr/> <p>Adapters: 50mlx4</p>	<p>Adapters: 5ml/7mlx4x4 Blood Vacuum Tube</p> <hr/> <p>Adapters: 10mlx3x4 G RIA Tube</p> <hr/> <p>Adapters: 20mlx4</p>
<p>No.5</p>  <p>Rotor capacity: 4200r/min Max RCF: 3100xg</p>	 <p>50mlx2x4 Bucket</p>  <p>15/10mlx4x4 Bucket</p>	<p>Adapters: 1.5mlx2 5ml/7ml Blood Vacuum Tube 5ml, 10ml 15ml, 20ml</p> <hr/> <p>Adapters: 1.5ml 5ml/7m Blood Vacuum Tube 5ml, 10ml</p>	 <p>15/10mlx6x4 Bucket</p> <hr/>  <p>15/10mlx8x4 Bucket</p> <hr/> <p>Adapters: 1.5ml 5ml/7m Blood Vacuum Tube 5ml, 10ml</p>
<p>No.6</p>  <p>Rotor capacity: 4000r/min Max Speed: 3310xg Max RCF: 500mlx4</p>	 <p>500ml Conical Tube 500mlx4</p>	 <p>500mlx4 V-bottom</p> <hr/>  <p>50mlx4x4 V-bottom</p>	 <p>250mlx4 V-bottom</p> <hr/>  <p>15mlx7x4 V-bottom</p>



#Let'sgive more ♥
Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.




















All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél: +33 467 158849 / +33 680760710 / +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S



KALSTEIN® Centrifuges Tabletop High Speed Centrifuge Model YR137-3



<p>No.2</p>  <p>Rotor capacity: 500mlx4 Max Speed: 4000r/min Max RCF: 3220xg</p>	 <p>500ml Bucket</p>	 <p>1.5mlx20x4</p>	 <p>5ml/7mlx20x4 Blood Vacuum Tube</p>	 <p>10mlx15x4 GRIA Tube</p>
	 <p>500ml Bottle</p>	 <p>15mlx8x4</p>	 <p>20mlx9x4</p>	 <p>50mlx4x4</p>
	 <p>Bucket</p>	 <p>100mlx2x4</p>	 <p>250mlx4</p>	
	 <p>4x2x96 Microplate</p>	 <p>1.5mlx25x4</p>	 <p>5mlx28x4 Blood Vacuum Tube</p>	 <p>10mlx20x6</p>
	 <p>4x24x7ml Blood Vacuum Tube</p>	 <p>15mlx10x4</p>	 <p>50mlx4x4</p>	 <p>250mlx4</p>



#Letsgivemore ♥
Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.
















All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél: +33 467 158849 / +33 680760710 / +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A.S.



KALSTEIN® Centrifuges Tabletop High Speed Centrifuge Model YR137-3



<p>No.4</p>  <p>Rotor capacity: 100mlx4 Max Speed: 5000r/min Max RCF: 4390xg</p>	 <p>100mlx4 Bucket</p>	<p>Adapters: 1.5mlx4x4</p> <p>Adapters: 5mlx2x4</p> <p>Adapters: 15mlx4</p> <p>Adapters: 50mlx4</p>	<p>Adapters: 5ml/7mlx4x4 Blood Vacuum Tube</p> <p>Adapters: 10mlx3x4 G RIA Tube</p> <p>Adapters: 20mlx4</p>
<p>No.5</p>  <p>Rotor capacity: 4200r/min Max RCF: 3100xg</p>	 <p>50mlx2x4 Bucket</p>  <p>15/10mlx4x4 Bucket</p>	<p>Adapters: 1.5mlx2 5ml/7ml Blood Vacuum Tube 5ml, 10ml 15ml, 20ml</p> <p>Adapters: 1.5ml 5ml/7m Blood Vacuum Tube 5ml, 10ml</p>	 <p>15/10mlx6x4 Bucket</p>  <p>15/10mlx8x4 Bucket</p> <p>Adapters: 1.5ml 5ml/7m Blood Vacuum Tube 5ml, 10ml</p> <p>Adapters: 1.5ml 5ml/7m Blood Vacuum Tube 5ml, 10ml</p>
<p>No.6</p>  <p>Rotor capacity: 4000r/min Max Speed: 3310xg Max RCF: 500mlx4</p>	 <p>500ml Conical Tube 500mlx4</p>	 <p>500mlx4 V-bottom</p>  <p>50mlx4x4 V-bottom</p>	 <p>250mlx4 V-bottom</p>  <p>15mlx7x4 V-bottom</p>



#Letsgivemore ♥
Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél: +33 467 158849 / +33 680760710 / +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale

<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed

The 3 articles of the week!

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW

How does the climate affect your microscope?

Kalstein Blog

Article of the week:
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW

What is the importance of the laboratory refrigerator?

Kalstein Blog

Article of the week:
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...
23 views • 3 years ago



Kalstein's Microscopes
445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER

2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89

Cel: +33 (0) 6 80 76 07 10

sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

